

SUPERase • In ™ RNase Inhibitor

Catalog Number AM2694, AM2696

Pub. No. 4393876 Rev. E

Contents	Quantity	Storage conditions
SUPERase•In™ RNase inhibitor, 20 U/μL	Cat. no. AM2694: 2,500 Units	Store at -20°C. Do not store in a frost-free
	Cat. no. AM2696: 10,000 Units	freezer.



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **www.lifetechnologies.com/support**.

Product Description

Nuclease-free SUPERase • In™ RNase Inhibitor is a protein-based ribonuclease inhibitor which noncovalently binds and inactivates a wide variety of RNases in a range of temperature (37–65°C) and pH (5.5–8.5) conditions. This product is distinct from ANTI-RNase Inhibitor (Cat. nos. AM2690, AM2692) and placental RNase Inhibitor protein (Cat. nos. AM2682, AM2684) in that it inactivates RNases I and T1 in addition to RNases A, B and C. Both ANTI-RNase Inhibitor and SUPERase • In™ RNase Inhibitor are distinct from RNase Inhibitor protein in that they have more robust interactions with RNAses, and do not release active RNases in the absence of DTT or other reducing agents. SUPERase • In™ RNase Inhibitor does not interfere with the activities of SP6, T7, and T3 RNA Polymerases, MMLV Reverse Transcriptase, or Taq DNA Polymerase.

Unit (U) definition: SUPERase • In™ RNase Inhibitor at 1 U/ μ L will block the degradation of 0.1 μ g/ μ L labeled RNA by 2.5 pg/ μ L of RNase A, 2.5 pg/ μ L of RNase I and 0.0075 U/ μ L of RNase T1, for 4 hours at 37°C, in 20 mM Tris-HCl, pH 7.5, 50 mM NaCl, 1 mM EDTA. Analysis is by denaturing PAGE.

Storage buffer (not included): 2 mM KH₂PO₄, 8 mM Na₂HPO₄, 2.7 mM KCl, 137 mM NaCl, pH 7.4 in 50% glycerol.

Using Superase In RNase Inhibitor

Use SUPERase \bullet In ** RNase Inhibitor at a final concentration of 1 U/µL to prevent RNA degradation in applications including cDNA synthesis, RT-PCR, *in vitro* transcription, *in vitro* translation, preparation of cell lysates, RNA isolation and storage, and in any application where ribonuclease inhibitor protein is used. In experiments requiring intact RNA, (e.g., transcription), avoid denaturation of the SUPERase \bullet In ** protein by heat, SDS or other detergents, urea, etc., which could release active RNases. The addition of DTT is not recommended for storage.

Note: Heating samples containing SUPERase • In™ RNase Inhibitor and DTT to 95°C may cause RNA to hang up in the wells of denaturing urea/polyacrylamide gels. Load directly onto gels without heating.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT

© 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For Research Use Only. Not for use in diagnostic procedures. The information in this guide is subject to change without notice.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

LIMITED TO THE USE THEREOF.

lifetechnologies.com

7 January 2013