

## 293 SFM II (1X)

### Description

293 SFM II is a serum-free, human, and animal-origin free, low-protein (<10 µg/mL) formulation optimized for high-density suspension culture of 293 (human embryonic kidney) cells. 293 SFM II supports HEK 293 cellular production of adenovirus and/or glycosylated recombinant protein, comparable to levels produced in serum-supplemented media. 293 SFM II is not recommended for adherent 293 cell cultures. 293 SFM II has been demonstrated to support the growth of Per.C6® cells and the high-density suspension culture of HeLa S3 cells.

Product	Catalog no.	Amount	Storage	Shelf life*
293 SFM II (1X), liquid	11686-029	1000 mL	2°C to 8°C; Protect from light	6 months

\* Shelf life duration is determined from Date of Manufacture.

### Product use

For Research Use Only. Not for use in diagnostic procedures.

### Important information

293 SFM II facilitates the adaptation of monolayer dependent 293 cells to growth in suspension culture.

### Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

### Prepare medium

- 293 SFM II Medium requires aseptic supplementation with 4 mM L-glutamine or GlutaMAX™-I prior to use.
- Antibiotics are not recommended; however, 5 mL/L of Penicillin-Streptomycin (100X) may be used when required.
- Addition of a surfactant such as Pluronic® F-68 is not required.

### Culture conditions

**Media:** Complete 293 SFM II Medium.

**Cell line:** HEK 293 cells.

**Culture type:** Suspension

**Culture vessels:** Shake flask, spinner bottle, or bioreactor.

**Temperature range:** 36°C to 38°C.

**Incubator atmosphere:** Humidified atmosphere of 8% CO<sub>2</sub> in air. Ensure proper gas exchange and minimize exposure of cultures to light.

### Adapt to suspension culture in 293 SFM II

**Note:** We offer 293-F and 293-H cells that have been pre-adapted to growth in 293 SFM II.

1. Aspirate media from cell monolayer and displace 293 cells from the flask surface by rapping the flask sharply against your hand or a protected surface several times.  
**Note:** Do not use trypsin or other proteolytic agents to dislodge cells.
2. Resuspend dislodged cells in 5 mL of 293 SFM II.  
**Note:** 293 cells cultured in 293 SFM II may grow as clusters of 2–10 cells.
3. Disperse clusters into a single-cell suspension by triturating with a small bore pipette or vortexing before passaging or counting. Optimal vortexing conditions must be determined based upon speed and duration versus viability.
4. Determine viable cell density using a Countess® Automated Cell Counter. Alternate methods (e.g., Coulter counter or hemocytometer) may also be used.

5. Dilute cells in prewarmed complete 293 SFM II to a viable cell density of  $1 \times 10^6$  cells/mL. Dispense cell suspension, up to a maximum of 30 mL/flask, into sterile 125-mL Erlenmeyer shake flasks.

6. Incubate the shake flask(s) on a rotary shaker (125–130 rpm) at 37°C in a humidified atmosphere of 8% CO<sub>2</sub> in air.

**Note:** Cells will not thrive in 293 SFM II at lower CO<sub>2</sub> levels (e.g., 5% CO<sub>2</sub>).

7. When the viable cell density reaches  $\sim 1.5 \times 10^6$  cells/mL, dilute the cells with prewarmed, complete 293 SFM II to  $2.5 \times 10^5$ – $3.0 \times 10^5$  cells/mL.

8. After the first passage, dilute cells to  $2.5 \times 10^5$ – $3.0 \times 10^5$  cells/mL whenever the viable cell density exceeds  $7.5 \times 10^6$ – $1 \times 10^6$  cells/mL. After several passages of consistent growth and viability in 293 SFM II, the culture is considered to be adapted.

**Note:** After adaptation to growth in serum-free suspension culture, it is possible to scale-up the cultures in spinner flasks or bioreactors. The appropriate spinner or impeller speed should be individually determined.

**Caution:** Some spinner apparatus emit significant heat and water-jacketed incubators usually cannot readily equilibrate to temperature variations. Temperatures >40°C are lethal to HEK 293 cells.

### Transfection protocol

**Important:** Complex formation of DNA with transfection reagents such as Lipofectamine® Plus and Lipofectamine® 2000 are inhibited by constituents of 293 SFM II. Therefore, these transfection reagents and media should not be used together.

### Prepare DNA-liposome complexes:

1. To each well of a six-well plate, add 1 mL Opti-MEM® Reduced Serum Medium.
2. For each transfection, add 3–4 µg of DNA to each well. Gently swirl the plate to mix.
3. For each transfection, add 10–20 µL of Lipofectamine® 2000 transfection reagent to each well. Gently swirl to mix the DNA and lipid.
4. Incubate at room temperature for 20–45 minutes, and repeatedly gently swirl the plate allowing DNA-liposome complexes to form. Although the solution may appear somewhat cloudy it will not impair transfection.

### Prepare 293 cells for transfection:

1. Determine the viable cell density of a suspension culture of 293 cells growing in 293 SFM II.  
**Note:** Cells should be in mid-log phase growth prior to transfection for optimal results.
2. Calculate the volume of cell suspension required to plate 293 cells into six-well plates at  $2 \times 10^6$  cells/well. Include two wells ( $4 \times 10^6$  cells) overage. Transfer this volume of cell suspension to a sterile 50-mL centrifuge tube. Bring the total volume to 50 mL with Opti-MEM® Reduced Serum Medium.
3. Centrifuge cell suspension at  $200 \times g$  for 5 minutes. Aspirate and discard the supernatant, being sure not to disturb the cell pellet.
4. Resuspend the cell pellet in 10 mL Opti-MEM® Reduced Serum Medium, pipetting up and down to achieve a single-cell suspension. Add Opti-MEM® Reduced Serum Medium to a total volume of 50 mL.
5. Centrifuge cell suspension at  $200 \times g$  for 5 minutes. Aspirate and discard the supernatant, being sure not to disturb the cell pellet.  
**Important:** This second wash step is imperative for optimal transfection efficiency.
6. Resuspend the cell pellet in 200  $\mu$ L of Opti-MEM® Reduced Serum Medium for every  $2 \times 10^6$  cells contained in the tube. This will yield a suspension of  $1 \times 10^7$  cells/mL.
7. Add 200  $\mu$ L of cell suspension to each well of the six-well plate containing DNA-liposome complexes. Pipet up and down to mix well.
8. Incubate at 37°C for 5 hours in a humidified 8% CO<sub>2</sub> incubator. There is no need to remove the transfection mixture, or to feed with growth medium.

**Note:** Transfection procedure can be scaled up or down by adjusting the amount of DNA, Lipofectamine® 2000 transfection reagent and cell concentration in proportion to the difference in surface area of the culture vessel.

### Transient expression

Harvest and assay cell extracts or stain cells *in situ* for reporter gene activity at 24 hours after the start of transfection.

### Stable expression

1. Passage the cells at 24 hours post-transfection using the same seeding density or split ratio that is normally used.
2. At 48 hours post-transfection, replace spent medium with medium containing the appropriate selective antibiotic (e.g., Geneticin®).

**Note:** Suspension cells will need to be collected by centrifugation ( $100 \times g$  for 4 minutes) before changing spent medium for medium containing selective antibiotic.

### Related products

Product	Catalog no.
L-Glutamine-200mM (100X), Liquid	25030
GlutaMAX™-I, 200mM (100X), Liquid	35050
293 F Cells, SFM Adapted	11625
293 H Cells, SFM Adapted	11631
Geneticin® (G-418 Sulfate)	11811
Opti-MEM® I Reduced Serum Medium (1X), Liquid	31985
Penicillin-Streptomycin 100X Solution	15070
293Fectin™ Transfection Reagent	12347
Lipofectamine® 2000 Transfection Reagent	11668
Trypan Blue Stain	15250
Countess® Automated Cell Counter	C10227

### Explanation of symbols and warnings

The symbols present on the product label are explained below:

				
Temperature Limitation	Manufacturer	Batch code	Use By:	Catalog number
				
Caution, consult accompanying documents	Consult instructions for use	Keep away from light	Sterilized using aseptic processing techniques	

### Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

### Important licensing information

This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support). For further assistance, email [techsupport@lifetech.com](mailto:techsupport@lifetech.com).

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Per.C6 is a trademark of Crucell Holland B.V.

**DISCLAIMER:** LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

[www.lifetechnologies.com](http://www.lifetechnologies.com)

 life technologies