



PRODUCT INFORMATION

Thermo Scientific Phire Reaction Buffer

Pub. No. MAN0013373

Rev. Date 18 March 2019 (Rev. C.00)

Store at -25 °C to -15 °C

Ordering information

Component	#F-524L
5X Phire Reaction Buffer	4 × 1.5 mL
50 mM MgCl ₂	2 × 1.5 mL
100 % DMSO	0.5 mL

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

Description

Thermo Scientific™ Phire™ Reaction Buffer is a 5X buffer solution optimized for DNA amplification using Phire Hot Start DNA polymerases. The buffer contains 7.5 mM MgCl₂ which provides 1.5 mM MgCl₂ in the final reaction conditions.

Important Notes

- Repeated freezing and thawing of the buffer can result in the precipitation or accumulation of MgCl₂ in insoluble form. For consistent results heat the buffer to 90 °C for 10 min and vortex prior to use if needed or store refrigerated.
- The freezing point of DMSO is 18-19 °C, so it does not melt on ice.

Important Licensing Information:

This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermofisher.com for Material Safety Data Sheet of the product.

© 2019 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.