

USER GUIDE

invitrogen
life technologies

ChargeSwitch[®] gDNA Rendered Meat Purification Kit

Purification of genomic DNA (gDNA) from cattle
feed, meal, and heparin products

Catalog Number CS400-100

Publication Number MAN0000574

Revision A.0

For testing of Food and Environmental samples only.

life
technologies

The information in this guide is subject to change without notice.

DISCLAIMER

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

LIMITED USE LABEL LICENSE No. 492: Environmental Testing, Quality Control/Quality Assurance Testing, Food and Agricultural Testing

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product (a) to perform internal research for the sole benefit of the purchaser; and (b) for environmental testing, quality control/quality assurance testing, food and agricultural testing, including reporting results of purchaser's activities in environmental testing, quality control/quality assurance testing, food and agricultural testing for a fee or other commercial consideration. No other right is hereby granted expressly, by implication, or by estoppel. This product is for environmental testing, quality control/ quality assurance testing, food and agricultural testing and research purposes only.

The purchase of this product does not grant the purchaser any additional rights, including (without limitation) the right to transfer or resell the product in any form or the right to use the product as a therapeutic agent or diagnostics test component. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Trademarks

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Tween is a registered trademarks of Croda Americas, LLC.

©2014 Thermo Fisher Scientific Inc. All rights reserved.

Contents

| | |
|--|-----------|
| About this guide | 4 |
| Revision history | 4 |
| ■ Product information | 5 |
| Product description | 5 |
| Kit contents and storage | 6 |
| Required materials not provided with the kit | 6 |
| ■ Methods | 8 |
| Workflow | 8 |
| Important procedural guidelines | 8 |
| Handling magnetic beads | 8 |
| Sample processing | 9 |
| Prepare the sample lysate | 9 |
| Bind the DNA to the beads | 10 |
| Wash the DNA on the beads | 11 |
| Elute the DNA from the beads | 11 |
| Quantify the DNA concentration | 12 |
| ■ APPENDIX A Troubleshooting | 13 |
| ■ APPENDIX B Supplemental information | 15 |
| System specifications | 15 |
| Documentation and support | 16 |
| Customer and technical support | 16 |
| Food Safety support | 16 |
| Limited product warranty | 16 |

About this guide

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Revision history

| Revision | Date | Description |
|----------|---------------|--|
| A.0 | December 2014 | <ul style="list-style-type: none">• Added heparin protocol.• Updated protocol organization to align with current style.• Updated user guide template with associated updates to covers, legal, document support, and safety sections.• New document control nomenclature. |
| 1.0 | 2006 | New document |



Product information

Product description

The ChargeSwitch[®] gDNA Rendered Meat Purification Kit allows rapid and efficient purification of PCR-ready genomic DNA (gDNA) from:

- Animal feed (500 mg)
- Animal meal (250 mg)
- Heparin (200–250 mg)

The kit uses magnetic bead-based technology for purification of gDNA without centrifugation, vacuum manifolds, or organic solvents. Genomic DNA can be prepared from sample lysates in less than 15 minutes, when processing 1–5 samples. The purified DNA has minimal RNA contamination, and it is suitable for analysis using real-time quantitative PCR (qPCR) or another method of choice, for identification of mammalian DNA in samples.

ChargeSwitch[®] technology is a magnetic bead-based technology that utilizes a switchable surface whose charge is dependent on the pH of the surrounding buffer (Figure 1).

- In low pH conditions, the ChargeSwitch[®] Magnetic Beads have a positive charge that binds the negatively charged nucleic acid backbone. Proteins and other contaminants remain unbound and are removed in an aqueous wash buffer.
- To elute the nucleic acid, the charge on the surface of the ChargeSwitch[®] Magnetic Beads is neutralized by raising the pH to >8.5 using a low salt elution buffer. Purified DNA elutes instantly into the elution buffer and is ready for use in downstream applications.



Figure 1 ChargeSwitch[®] technology



Kit contents and storage

| Component | Amount ^[1] | Storage ^[2] |
|---|-----------------------|------------------------|
| ChargeSwitch [®] Magnetic Beads (25 mg/ml in 10 mM MES, pH 5.0, 10 mM NaCl, 0.1% Tween 20) | 1.1 mL | Room temperature |
| ChargeSwitch [®] SDS | 10 mL | |
| ChargeSwitch [®] 10% Detergent | 10 mL | |
| ChargeSwitch [®] Lysis Buffer | 100 mL | |
| ChargeSwitch [®] Precipitation Buffer (N5) | 38.5 mL | |
| ChargeSwitch [®] Wash Buffer (W12) | 200 mL | |
| ChargeSwitch [®] Elution Buffer (E5; 10 mM Tris HCl, pH 8.5) | 15 mL | |

^[1] 50 purifications from feed or meal samples; 25 purifications from heparin samples.

^[2] All components are guaranteed stable for 6 months when stored properly.

Required materials not provided with the kit

Unless otherwise indicated, all materials are available from Life Technologies (www.lifetechnologies.com). MLS: Fisher Scientific (www.fisherscientific.com) or other major laboratory supplier.

| Product | Catalog number |
|---|----------------|
| MagnaRack [™] Magnetic Rack | CS15000 |
| For animal feed samples: No. 6 wire-mesh sieve (3.35 mm opening) | MLS |
| Sterile 2.0-mL microcentrifuge tubes, with locking lids Note: Locking lids are recommended to prevent the lid from popping open during the 95°C incubation. | MLS |
| Adjustable pipettes and aerosol barrier pipette tips | MLS |
| 1.5-mL microcentrifuge tubes | MLS |
| Microcentrifuge | MLS |
| Water bath or heat block set to 95°C | MLS |
| Smart Spatula | MLS |
| Materials for DNA quantification using one of these methods: | |
| UV absorbance: spectrophotometer and accessories | MLS |




| Product | Catalog number |
|---|--------------------------------|
| Fluorescence technology: Quant-iT™ DNA assay kits <ul style="list-style-type: none"><li data-bbox="505 317 1003 352">• Quant-iT™ DNA Assay Kit, High Sensitivity<li data-bbox="505 359 964 394">• Quant-iT™ DNA Assay Kit, Broadrange<li data-bbox="505 401 943 436">• Quant-iT™ PicoGreen® dsDNA Assay | — Q33120 Q33130 P7589 |



Methods

Workflow



Important procedural guidelines


Handling magnetic beads

For best results:

- Do not freeze ChargeSwitch[®] Magnetic Beads—freezing damages the magnetic property of the beads. Store the beads at room temperature.
- Always keep the ChargeSwitch[®] Magnetic Beads in solution. Do not allow the beads to dry, including during washing procedures, as this renders the beads non-functional.
- Vortex the ChargeSwitch[®] Magnetic Beads to resuspend thoroughly before pipetting.
- During mixing steps, to avoid forming bubbles:
 - Use an adjustable pipette set to a specific volume as directed in the protocol.
 - Pipet up and down gently with the pipette tip submerged in the solution.



- To aspirate the supernatant after bead washing, point the pipette tip away from the beads, and carefully remove the supernatant without disturbing the beads.



- Discard ChargeSwitch[®] Magnetic Beads after use.

Sample processing

To maximize DNA yield:

- Maintain a sterile environment when handling DNA to avoid any contamination from DNases.
- Ensure that no DNases are introduced into the solutions supplied with the kit.
- Make sure all equipment that contacts DNA is sterile, including pipette tips and tubes.
- Perform all steps at room temperature, unless otherwise noted.
- Make sure that:
 - The ChargeSwitch[®] Wash Buffer is removed completely before elution.
 - The ChargeSwitch[®] Magnetic Beads are fully resuspended during the elution step.

Prepare the sample lysate

1. Transfer the indicated amount of sample to sterile 2.0-mL microcentrifuge tubes.

| Sample type | Amount |
|-------------|---|
| Animal feed | 500 mg; transfer to two tubes. Note: Sift feed with a No. 6 wire-mesh sieve or equivalent before weighing sample, to remove large components of the feed (for example, whole corn kernels and whole grains) that can reduce the sensitivity of downstream assays. |
| Animal meal | 250 mg |
| Heparin | 200–250 mg |

2. At room temperature, add 1 mL of ChargeSwitch[®] Lysis Buffer (L15) and 100 µL of ChargeSwitch[®] SDS to each tube.



- Mix the sample with a ChargeSwitch® Lysis Buffer and ChargeSwitch® SDS either with Smart Spatula or by inversion

| For these samples... | To mix... |
|-------------------------------------|---|
| Heparin, some feed and meal samples | Use a Smart Spatula, then lock the tube lid. |
| All other samples | Lock the tube lid, invert five times, then gently tap the tube on a hard surface to remove sample that is stuck to the lid. |

- Incubate in a 95°C water bath or heat block for 5 minutes.



CAUTION! Be careful when opening the tubes after heating.

- Carefully open each tube, add 400 µL of ChargeSwitch® Precipitation Buffer (N5) to the lysate, then mix by inversion.
- Place each tube on ice for 5 minutes to precipitate the proteins.
- Centrifuge the tube at 17,000 × g for 5 minutes at room temperature to pellet the debris.

Note: The supernatant contains the lysate.

Bind the DNA to the beads

- Carefully transfer the indicated volume of lysate to a sterile 1.5-mL microcentrifuge tube.

| Sample type | Volume |
|-------------|--|
| Animal feed | 500 µL from each tube of lysate for a total of ~1000 µL |
| Animal meal | 500 µL |
| Heparin | 1200 µL Note: It may not be possible to recover 1200 µL of lysate, due to differences in heparin products. |

- Add 100 µL of ChargeSwitch® 10% Detergent (D1) to the tube of lysate.
- Thoroughly vortex the tube of ChargeSwitch® Magnetic Beads, and add the indicated volume of beads to the sample.

| Sample type | Volume |
|-------------|--------|
| Animal feed | 20 µL |
| Animal meal | 20 µL |
| Heparin | 40 µL |



- Mix gently, without forming bubbles, by pipetting up and down 5 times using a 1-mL pipette set to 900 μ L.

IMPORTANT! Avoid forming bubbles by ensuring that the pipette tip is fully submerged during mixing, and by pipetting up and down gently.

- Incubate the tube at room temperature for 1 minute.
- Place the tube on the MagnaRack™ Magnetic Rack until the beads have formed a tight pellet and the supernatant has cleared (~ 1 minute).
- Without removing the tube from the magnet, carefully aspirate and discard the supernatant without disturbing the bead pellet.

Note: The lysate is viscous; angle the pipette tip so that it is pointed away from the pellet, and aspirate slowly, to avoid disturbing the bead pellet.

Wash the DNA on the beads

- Remove the tube from the MagnaRack™ magnet, add 1 mL of ChargeSwitch® Wash Buffer (W12) to the tube, and pipet up and down 5 times using a 1-mL pipette set to 900 μ L.

Note: Pipet up and down to gently mix without forming bubbles.

- Place the tube on the magnet for ~1 minute until the beads have formed a tight pellet and the supernatant is clear.
- Without removing the tube from the magnet, carefully aspirate and discard the supernatant without disturbing the bead pellet.

Note: Heparin samples are very viscous; pipet slowly to ensure that beads are not removed from the magnet. Removal of beads will decrease the DNA yield.

- Remove the tube from the MagnaRack™ magnet, add 750 μ L of ChargeSwitch® Wash Buffer (W12) and mix up and down 5 times.
- Place the tube on the magnet for ~1 minute until the beads have formed a pellet and the supernatant is clear.
- Without removing the tube from the magnet, carefully aspirate and discard the supernatant without disturbing the bead pellet.
- Repeat step 4–step 6 one more time.

Remove the supernatant as completely as possible after the final wash.

Note: Keep the pelleted beads, which contain the bound DNA.

Elute the DNA from the beads

- Remove the tube containing the pelleted magnetic beads from the magnet.

Note: If any supernatant is visible, carefully remove it before proceeding.



2. Add the indicated amount of ChargeSwitch® Elution Buffer (E5) to the bead pellet.

| Sample type | Volume |
|-------------|--------|
| Animal feed | 100 µL |
| Animal meal | 100 µL |
| Heparin | 75 µL |

3. Gently pipet up and down gently 10 times using a pipette set to 20 µL less than the volume of buffer used.
4. Incubate at room temperature for 1 minute.
(Optional) For maximum yield, this incubation can be extended up to 5 minutes total, with gentle tip-mixing after 2 minutes.
5. Place the tube on the magnet for ~1 minute until the beads have formed a tight pellet and the supernatant is clear.
Note: The supernatant contains the purified DNA.
6. Without removing the tube from the magnet, carefully transfer the supernatant containing the DNA to a new, sterile microcentrifuge tube without disturbing the pellet.
7. Discard the used ChargeSwitch® Magnetic Beads.

The purified gDNA is ready for use.

STOPPING POINT (Optional) Store at -20°C. Avoid repeated freezing and thawing of gDNA.

Quantify the DNA concentration

Perform DNA quantitation using UV absorbance at 260 nm or Quant-iT™ Kits.

- **UV absorbance:**
 - a. Prepare a dilution of the DNA solution in 10 mM Tris HCl, pH 8.5, mix well.
 - b. Measure the absorbance at 260 nm (A_{260}) of the dilution in a spectrophotometer (using a cuvette with an optical path length of 1 cm) blanked against 10 mM Tris-HCl, pH 7.5.
 - c. Calculate the concentration of DNA using the formula:
$$\text{DNA } (\mu\text{g/mL}) = A_{260} \times 50 \times \text{dilution factor}$$

For DNA, $A_{260} = 1$ for a 50 µg/mL solution measured in a cuvette with an optical path length of 1 cm.
- **Quant-iT™ Kits:** These kits provide a rapid, sensitive, and specific fluorescent method for dsDNA quantitation. Each kit contains a fluorescent quantitation reagent, DNA standards for standard curve, and a pre-made buffer to allow fluorescent DNA quantitation using standard fluorescence microtiter plate readers or fluorometers.



Troubleshooting

| Observation | Possible cause | Recommended action |
|---|---|---|
| Low DNA yield | Incomplete lysis | Reduce the amount of starting material. |
| | | Be sure to add ChargeSwitch® SDS during lysis. |
| | Poor quality of starting material | Use fresh sample and process immediately after collection, or freeze the sample at -80°C or in liquid nitrogen. The yield and quality of DNA isolated depends on the type and age of the starting material. |
| | Incorrect handling of ChargeSwitch® Magnetic Beads | Vortex the tube containing the ChargeSwitch® Magnetic Beads to fully resuspend the beads before adding them to your sample. |
| | Bead pellet was disturbed or lost during binding or washing steps | Keep the sample on the MagnaRack™ when removing supernatant during the binding or washing steps. |
| | | Remove the supernatant without disturbing the bead pellet. See “Handling magnetic beads” on page 8. |
| | Suboptimum elution conditions | After adding ChargeSwitch® Elution Buffer (E5) to the sample, pipet up and down to resuspend the ChargeSwitch® Magnetic Beads before incubation. |
| Do not use water to elute DNA. Use ChargeSwitch® Elution Buffer (E5). | | |
| Preheating Elution Buffer (E5) to 55–60°C may improve yield. | | |
| No DNA recovered | Water used for elution | The elution buffer must be pH 8.5–9.0 or the DNA will remain bound to the ChargeSwitch® Magnetic Beads. Use ChargeSwitch® Elution Buffer (E5). |
| | ChargeSwitch® Magnetic Beads were stored or handled improperly | Store beads at room temperature. Do not freeze the beads, as they will become irreparably damaged. |
| | | Make sure that the beads are in solution at all times and are not dried. Dried beads are non-functional. |
| Eluate containing DNA is discolored | Magnetic pellet disturbed during elution | Place the sample on the MagnaRack™ until the beads form a tight pellet. |
| | | Remove the eluate to a sterile microcentrifuge tube or sterile microtiter plate, without disturbing the bead pellet. |



| Observation | Possible cause | Recommended action |
|----------------------------|------------------------------------|--|
| DNA is sheared or degraded | Lysate mixed too vigorously | Use an appropriate pipette set to a volume that is lower than the total volume of the solution used to mix the sample. Pipet up and down gently to mix. |
| | Bubbles formed during mixing steps | Make sure that the pipette tip is submerged in the solution during mixing. |
| | DNA repeatedly frozen and thawed | Aliquot the eluted DNA and store at 4°C or -20°C. Avoid repeated freezing and thawing. |
| | DNA contaminated with DNases | Maintain a sterile environment while working (for example, wear gloves and use DNase-free reagents). |



Supplemental information

System specifications

Table 1 ChargeSwitch[®] gDNA Rendered Meat Purification Kit specifications

| Feature | Specification |
|-----------------------|--|
| Starting material | <ul style="list-style-type: none">• Animal feed (500 mg)• Animal meal (250 mg)• Heparin (200–250 mg) |
| Bead binding capacity | 1 mg bead binds 5–10 µg gDNA |
| Bead size | 1 µm |
| Bead concentration | 25 mg/ mL |
| Bead storage buffer | 10 mM MES, pH 5.0, 10 mM NaCl, 0.1% Tween [®] 20 |
| Elution volume | <ul style="list-style-type: none">• 100 µL (feed, meal)• 75 µL (heparin) |
| DNA yield | Up to 7 µg |

Documentation and support

Customer and technical support

Visit www.lifetechnologies.com/support for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Food Safety support

Website: www.lifetechnologies.com/foodsafety

Support email: foodsafety@lifetech.com

Phone number in North America: 1-800-500-6855

Phone number outside of North America: Visit www.lifetechnologies.com/support, select the link for phone support, and select the appropriate country from the dropdown menu.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com
lifetechnologies.com

11 December 2014

