









Publication No. MAN0013997 Rev A.0

	Package Contents	Catalog Numbers K210009XP	Amount: 2 preps
	Storage Conditions	<ul style="list-style-type: none"> Store all components at room temperature. 	
	Required Materials	<ul style="list-style-type: none"> Vacuum source equipped with regulator (capable of -600 to -800 mbar) Appropriately sized tubes and bottles 1000-mL Stericup® Receiver flask 250-mL Stericup® Receiver flask Centrifuge and rotor capable of >12,000 x g at 4°C 	
	Timing	<p>Bacterial culture: overnight Purification: 90 minutes</p>	
	Selection Guide	<p>Go online to view related products: PureLink® Nucleic Acid Purification Kits Expi293™ Expression System</p>	
	Product Description	<ul style="list-style-type: none"> The PureLink® HiPure Expi Gigaprep Kit provides users with the ability to isolate large quantities of transfection-grade plasmid DNA using an enhanced anion exchange resin. The kit includes filtration columns to provide bacterial filtration without centrifugation in a single unit. The PureLink® HiPure Expi Gigaprep Kit typically isolates 12 mg of high quality, ultrapure plasmid DNA with inherently low endotoxin levels from 2.5 L of bacterial culture. High Yield – Isolate over 14 mg of high quality plasmid DNA from a single purification using 5 L of bacterial culture volume. Purity – Low endotoxin levels (0.1–1.0 EU/μg), and A260/280 >1.8, making it ideal for mammalian cell transfection. 	
	Important Guidelines	<ul style="list-style-type: none"> Add RNase A to the Resuspension Buffer (R3) and mix well (see instructions on label). Indicate that RNase A has been added on the bottle label. Store at 4°C. If precipitate is observed in the Lysis Buffer (L7), warm the buffer in a 37°C water bath until the solution clears. Swirl contents gently to resuspend. Grow transformed <i>E. coli</i> in LB medium. Use 2.5 L (high copy number plasmid) or 5 L (low copy number plasmid) of an overnight culture. Do not over-dry DNA. If the DNA pellet is difficult to resuspend, allow the pellet to incubate in TE Buffer for a longer period of time. 	
	Online Resources	<p>Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support.</p>	

Gigaprep Plasmid Isolation Protocol

Steps		Procedure Details
1	Harvest	1. Sediment cells by centrifugation at $4,000 \times g$ for 15 min at 4°C . Discard all medium.
2	Resuspend	2. Add 125 mL Resuspension Buffer (R3) with RNase A to the cell pellet and resuspend the pellet until it is homogeneous.
3	Lyse	3. Add 125 mL Lysis Buffer (L7). Mix gently by inverting the capped tube until the mixture is homogeneous. Do not vortex. Incubate at room temperature for 5 minutes.
4	Precipitate	4. Add 125 mL Precipitation Buffer (N3). Mix immediately by inverting the tube until the mixture is homogeneous. Do not vortex.
5	Clarify	5. Pour the lysate into a lysate filtration cartridge attached to a receiver flask . Incubate for 5 minutes. Connect a vacuum source and filter the lysate.
6	Wash	6. Add 50 mL Wash Buffer (W8) to the filtration cartridge and gently stir precipitate with a spatula. Apply vacuum. The clarified lysate contains the plasmid DNA .
7	Equilibrate	7. Add 200 mL Equilibration Buffer (EQ1) to a DNA-binding cartridge attached to a receiver flask . Connect a vacuum source and drain the cartridge.
8	Bind	8. Load the clarified lysate (from step 6) onto the DNA-binding cartridge. Apply vacuum and drain solution.
9	Wash	9. Add 275 mL Wash Buffer (W8) and apply vacuum. Repeat wash step. Attach DNA-binding cartridge to a new receiver flask .
10	Elute	10. Add 100 mL Elution Buffer (E4) to the DNA-binding cartridge. Apply soft vacuum (-100 to -200 mbar) and draw 30–40 mL of solution. Stop the vacuum and incubate for 1 minute. Apply vacuum to all the liquid has passed from the cartridge.
11	Precipitate and Wash	11. Add 0.7 volume of isopropanol to the eluate. Mix well. Centrifuge at $>12,000 \times g$ for 30 minutes at 4°C . Remove and discard the supernatant. Wash the DNA pellet in 20 mL 70% ethanol. Centrifuge at $>12,000 \times g$ for 10 minutes at 4°C . Remove the supernatant.
12	Resuspend	12. Air-dry the pellet for 10 minutes, then resuspend the purified plasmid DNA in TE Buffer (TE). Store plasmid DNA at -20°C .

Disclaimer

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support, visit www.lifetechnologies.com/support.

21 April 2015

ThermoFisher
SCIENTIFIC