



PRODUCT INFORMATION

Thermo Scientific

Phusion U Multiplex PCR Master Mix

Pub. No. MAN0013365
Rev. Date 27 June 2018 (Rev. B.00)

#
Lot
Expiry Date

Store at -20°C

Ordering information

Table with 4 columns: Component, 2X Phusion U Multiplex PCR Master Mix, 2X Phusion U Green Multiplex PCR Master Mix, Water, nuclease-free. Rows include #F-562S, #F-562L, #F-564S, #F-564L.

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

Rev.3

1. Introduction

Thermo Scientific™ Phusion™ U PCR Master Mix is a ready-to-use, 2X end-point PCR master mix designed for simultaneous amplification of multiple targets in a single tube. Over 20 primer pairs may be combined into a single reaction for highly specific and efficient multiplexing over a broad range of primer and template concentrations.

The Phusion U Multiplex PCR Master Mix contains a proprietary buffer with balanced concentrations of all PCR components eliminating the need for tedious optimization. The Phusion U Green Multiplex PCR Master Mix further simplifies the workflow - it includes a density reagent and two electrophoresis tracking dyes for direct loading of PCR products on gels.

2. Storage conditions

Phusion U Multiplex PCR Master Mix should be stored at -20 °C immediately upon arrival. For short term use, the master mix can be stored at 2-8 °C for up to 3 months without compromising performance.

3. Important Notes

- Use 98 °C for denaturation (see 6.1 & 6.2).
The annealing rules are different from many common DNA polymerases (such as Taq DNA polymerases). Read Sections 5.2 and 6.3 carefully for instructions.
Use 15 -30 s/kb for extension (see 6.4).
Phusion U Hot Start DNA Polymerase produces DNA products with blunt ends.

4. Guidelines for using Phusion U Multiplex PCR Master Mix

Carefully mix and centrifuge all tubes before opening to ensure homogeneity and improve recovery. When using Phusion U Multiplex PCR Master Mix, it is not necessary to perform the PCR setup on ice. Due to the unique nature of Phusion U DNA Polymerase, optimal reaction conditions may differ from standard enzyme protocols.

Table 1. Pipetting instructions (add items in this order)

Table with 4 columns: Component, 20 µL rxn, 50 µL rxn, Final conc. Rows include H2O, 2X Phusion U Multiplex PCR Master Mix, Primer Mix, Template DNA.

*Recommended final primer concentration is 0.3 µM, but it can be varied in a range of 0.2-0.4 µM.

Table 2. Cycling instructions

Table with 4 columns: Cycle step, Temp., Time, Cycles. Rows include Initial denaturation, Denaturation, Annealing, Extension, Final Extension.

*Calculate extension time based on the size of the largest amplicon and the level of multiplexing (see 6.4)

**30 cycles give sufficient yield of PCR product in most cases. Larger number of cycles may be required to increase assay sensitivity when low amounts of template DNA are used.

5. Notes about reaction components

5.1. 2X Phusion U Multiplex PCR Master Mix

Phusion U Multiplex PCR Master Mix contains all the necessary reaction components for PCR except for template DNA and primers. The master mix composition is designed to give optimal results in simultaneous amplification of multiple targets ranging in size from 70 to 2,500 bp.

It contains a proprietary buffer and balanced ratio of all components, eliminating the need for tedious optimization. The master mix provides 1.5 mM MgCl2 and 200 µM of each dNTP in the final reaction concentration.

5.2. Primers

The recommended final primer concentration is 0.3 µM of each primer. If required, the primer concentration may be optimized in the range between 0.2 µM and 0.4 µM. Special attention to primer design parameters is critical for a successful multiplex PCR. Usually, primers of 21-34 nt length are used. Optimal GC content of the primer is 40-60%. Ideally, G and C nucleotides should be distributed uniformly along the primer.

Table 3. Recommended differences between adjacent amplicons

Table with 2 columns: Size of amplicons, Minimum difference. Rows include 50 - 200 bp, 200 - 700 bp, 700 - 1000 bp, 1000 - 2500 bp.

Note that primer quality is also an essential factor for good multiplex PCR results. Primers should be purchased desalted or purified, from a reliable oligonucleotide manufacturer. Dissolve the primers in TE buffer and check the concentration by spectrophotometry.

5.3. Template

The recommended amount of DNA template depends on the range of amplicon sizes. To simultaneously amplify DNA fragments up to 1 kb, it is recommended to use 0.1 ng - 1 µg of template DNA. If the size of the largest amplicon is longer than 1 kb, the amount of template DNA should not exceed 250 ng.

(continued on reverse page)

6. Notes about cycling conditions

6.1. Initial denaturation

Denaturation should be performed at 98 °C. Due to the high thermostability of Phusion U Hot Start DNA Polymerase even higher than 98 °C denaturation temperatures can be used. A 30-second initial denaturation at 98 °C is recommended for most templates. Some templates may require longer initial denaturation time, and the length of the initial denaturation time can be extended up to 3 minutes.

6.2. Denaturation

Keep the denaturation as short as possible. Usually 5-10 seconds at 98 °C is enough for most templates. Note: the denaturation time and temperature may vary depending on the ramp rate and temperature control mode of the cycler.

6.3. Primer annealing

The optimal annealing temperature for Phusion U Hot Start DNA Polymerase may be significantly different from the annealing temperature with other DNA polymerases. **Always use the Tm calculator and instructions on website: www.thermofisher.com/tmcalculator to determine the Tm values of your primers and optimal annealing temperature.**

If necessary, use a temperature gradient PCR to find the optimal annealing temperature for each template-primer mix combination.

6.4. Extension

The extension should be performed at 72 °C. The recommended extension time is 15-30 s/kb and depends on the amplicon length and the level of multiplexing. Calculate the extension time based on the size of the largest amplicon. Increase extension time to 30 s/kb for longer amplicons or higher level of multiplexing.

7. References

1. Nord K. *et al.* (1997) *Nature Biotechnol.* 15: 772–777.
2. Wikman M. *et al.* (2004) *Protein Eng. Des. Sel.* 17: 455–462.


8. Troubleshooting

No product at all or low yield	
<ul style="list-style-type: none">• Repeat and make sure that there are no pipetting errors.• Make sure that the cycling protocol was performed as recommended.• Titrate template amount. As a starting point we recommend using 20-50 ng of genomic DNA per 50 µL reaction.• Increase extension time.• Increase the number of cycles.	<ul style="list-style-type: none">• Check the concentration and quality of template DNA. Process the sample carefully to minimize DNA damage.• Check the purity and concentration of the primers. Check primer design.• Optimize annealing temperature using temperature gradient PCR.
Non-specific products	
<ul style="list-style-type: none">• Make sure that the extension time used was not too long. (Recommended extension time is 15-30 s/kb).• Increase annealing temperature or perform a temperature gradient PCR.	<ul style="list-style-type: none">• Titrate template amount.• Reduce the total number of cycles.• Decrease primer concentration.• Check primer design.• Validate each primer pair in singleplex PCR. Use only those primers which produce a single product of correct size.
High amount of primer dimers	
<ul style="list-style-type: none">• Check primer design.• Ensure that each primer pair produces a single band without primer dimers in singleplex PCR reactions.	<ul style="list-style-type: none">• Reduce the number of cycles.• Reduce primer concentration in final reaction mix.
Uneven amplification	
<ul style="list-style-type: none">• Annealing temperature is too low or too high. Test annealing temperatures in increments of 2 °C or perform temperature gradient PCR.• Do not exceed 250 ng of template DNA to achieve efficient amplification of amplicons longer than 1 kb.• Increase primer concentration for lower yield amplicons.	<ul style="list-style-type: none">• Decrease primer concentration for too abundant amplicons.• Make sure that the thermal cycler has been used according manufacturer's recommendations. Appropriate instrument maintenance might be critical for successful and reproducible multiplex PCR assays. Use a high quality instrument with excellent thermal uniformity and accuracy.

CERTIFICATE OF ANALYSIS

Multiplex PCR assay

Performance in multiplex PCR is tested by the amplification of 15 fragments (99 - 1606 bp) from 50 ng human genomic DNA in 50 µL reaction volume.

Quality authorized by:  Jurgita Zilinskiene

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermofisher.com for Material Safety Data Sheet of the product.

© 2018 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.