









Publication No. MAN0015741 Rev A.0

	<b>Package Contents</b>	<b>Catalog Numbers</b> A31233	<b>Amount:</b> 2 preps
	<b>Storage Conditions</b>	<ul style="list-style-type: none"><li>Store all components at room temperature.</li></ul>	
	<b>Required Materials</b>	<ul style="list-style-type: none"><li>Vacuum source capable of -600 to -800 mbar</li><li>Appropriately sized tubes and bottles</li><li>1000-mL Stericup™ Receiver flask or 45-mm thread sterile glass bottle</li><li>250-mL Stericup™ Receiver flask or 45-mm thread sterile glass bottle</li><li>Centrifuge and rotor capable of &gt;12,000 x g at 4°C</li></ul>	
	<b>Timing</b>	Bacterial culture: overnight Purification: 120 minutes	
	<b>Selection Guide</b>	Go online to view related products: <b>PureLink™ Nucleic Acid Purification Kits</b> <b>Expi293™ and ExpiCHO™ Expression Systems</b>	
	<b>Product Description</b>	<ul style="list-style-type: none"><li>The PureLink™ Expi Endotoxin-Free Giga Plasmid Purification Kit enables isolation of highly pure, endotoxin-free (&lt;0.1 EU/μg) advanced transfection grade plasmid DNA for use in your most sensitive downstream applications.</li><li>Isolate up to 15 mg of high quality, endotoxin-free ultrapure plasmid DNA from 2.5–5 L of bacterial culture.</li><li>Endotoxin-free (&lt;0.1 EU/μg) DNA with ultra-low levels of RNA, gDNA or protein contamination that is ideal for mammalian cell transfection or <i>in vivo</i> experiments.</li><li>Lysis Indicator turns blue during lysis and helps to ensure complete lysis and subsequent neutralization.</li></ul>	
	<b>Important Guidelines</b>	<ul style="list-style-type: none"><li>Add RNase A to the Resuspension Buffer (R3) and mix well (see instructions on label). Indicate that RNase A has been added on the bottle label. Store at 4°C.</li><li>Add Lysis Indicator to Resuspension Buffer (R3) at a 1:1000 ratio if desired.</li><li>If precipitate is observed in the Lysis Buffer (L7), warm the buffer in a 37°C water bath until the solution clears. Swirl contents gently to resuspend.</li><li>Grow transformed <i>E. coli</i> in LB medium. Use 2.5 L (high copy number plasmid) or 5 L (low copy number plasmid) of an overnight culture.</li><li>Perform protocol (step 6 onward) using endotoxin-free/pyrogen free components.</li><li>Do not over-dry DNA. If the DNA pellet is difficult to resuspend, allow the pellet to incubate in TE Buffer for a longer period of time.</li></ul>	
	<b>Online Resources</b>	Visit our product page for additional information and protocols. For support, visit <a href="http://thermofisher.com/support">thermofisher.com/support</a> .	

# Gigaprep plasmid isolation protocol

Steps	Procedure Details
<b>1</b>	<b>Harvest</b> Sediment cells by centrifugation at $4,000 \times g$ for 15 min at $4^{\circ}\text{C}$ . Discard all medium.
<b>2</b>	<b>Resuspend</b> Add 125 mL Resuspension Buffer (R3) with RNase A to the cell pellet and resuspend the pellet until it is homogeneous.
<b>3</b>	<b>Lyse</b> Add 125 mL Lysis Buffer (L7). Mix gently by inverting the capped tube until the mixture is homogeneous. Do not vortex. Incubate at room temperature for 5 minutes. <b>Note:</b> If using Lysis Indicator, the solution turns blue.
<b>4</b>	<b>Precipitate</b> Add 125 mL Precipitation Buffer (N3). Mix immediately by inverting the tube until the mixture is homogeneous. Do not vortex. <b>Note:</b> If using Lysis Indicator, the solution turns completely colorless.
<b>5</b>	<b>Clarify</b> Pour the lysate into a <b>lysate filtration cartridge</b> attached to a <b>receiver flask</b> . Incubate for 2 minutes. Connect a vacuum source and filter the lysate.
<b>6</b>	<b>Wash</b> Add 50 mL Wash Buffer (W8) to the filtration cartridge and gently stir precipitate with a spatula. Apply vacuum. The <b>clarified lysate contains the plasmid DNA</b> .
<b>7</b>	<b>Remove Endotoxin</b> Add 30 mL Endotoxin Removal Buffer (ER) to the clarified lysate and mix by inverting the bottle 10 times.
<b>8</b>	<b>Equilibrate</b> Add 200 mL Equilibration Buffer (EQ1) to a <b>DNA-binding cartridge</b> attached to a <b>receiver flask</b> . Connect a vacuum source and drain the cartridge.
<b>9</b>	<b>Bind</b> Load the <b>clarified lysate</b> (from step 7) onto the DNA-binding cartridge. Apply vacuum and drain solution.
<b>10</b>	<b>Wash</b> Add 275 mL Wash Buffer (W8) and apply vacuum. Repeat wash step. Attach <b>DNA-binding cartridge</b> to a <b>new receiver flask</b> .
<b>11</b>	<b>Elute</b> Add 100 mL Elution Buffer (E4) to the DNA-binding cartridge. Apply vacuum and draw 10–20 mL of solution. Stop the vacuum and incubate for 1 minute. Apply vacuum until all the liquid has passed from the cartridge.
<b>12</b>	<b>Precipitate and Wash</b> Make aliquots of the eluate and add 0.7 volume of isopropanol to each tube. Mix well. Centrifuge at $>12,000 \times g$ for 30 minutes at $4^{\circ}\text{C}$ . Remove and discard the supernatant. Wash the DNA pellets in 5 mL 70% ethanol for each tube. Centrifuge at $>12,000 \times g$ for 10 minutes at $4^{\circ}\text{C}$ . Remove the supernatant.
<b>13</b>	<b>Resuspend</b> Air-dry the pellets for 10 minutes, then resuspend the purified plasmid DNA in TE Buffer (TE). Store plasmid DNA in endotoxin-free tubes at $-20^{\circ}\text{C}$ .

**Disclaimer:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit [thermofisher.com/support](http://thermofisher.com/support)

2 March 2016

**ThermoFisher**  
SCIENTIFIC