

PureLink® Quick Gel Extraction and PCR Purification Combo Kit

Catalog number K2200-01

Publication Part Number 7015020

ΜΔΝηηη4376

Revision Date 14 September 2011

Purifying PCR Products



Before Starting. Add 96–100 % ethanol to Wash Buffer (W1) and isopropanol to Binding Buffer (B2) according to the labels on the bottles.

- 1. Combine. Add 4 volumes of Binding Buffer (B2) with isopropanol to 1 volume of a PCR sample (50-100 µL). Mix well.
- 2. Load. Pipet the sample into a PureLink® Clean-up Spin Column in a Wash Tube. Centrifuge the column at $>10,000 \times g$ for 1 minute. Discard the flow-through.



Wash. Re-insert the column into the Wash Tube and add 650 µL Wash Buffer (W1) with ethanol. Centrifuge the column at $>10,000 \times g$ for 1 minute.



Remove ethanol. Discard the flow-through and place the column in the same Wash Tube. Centrifuge the column at maximum speed for 2-3 minutes.




- **Elute.** Place the column into a clean 1.7-mL Elution Tube. Add 50 µL Elution Buffer (E1) to the column. Incubate the column at room temperature for 1 minute. Centrifuge the column at maximum speed for 1 minute.
- **Store.** *The elution tube contains the purified PCR product.* Store the purified DNA at 4°C for immediate use or at -20°C for long-term storage.

Intended Use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.



Purifying DNA from Gels Using a Centrifuge



1 MIN)

Before Starting. Add ethanol to Wash Buffer (W1) according to the label on the bottle. Equilibrate a water bath or heat block to 50°C.

- 1. Excise. Use a clean, sharp razor blade to excise a minimal area of gel containing the DNA fragment of interest.
- **2. Weigh.** Using a scale sensitive to 0.001 g, weigh the gel slice containing the DNA fragment.
- 3. Solubilize. Add Gel Solubilization Buffer (L3) to the excised gel in a tube as indicated in the table. Incubate the tube at 50°C for 10 minutes (or longer for large gel slices and high concentration gels), and invert the tube every 3 minutes. After the gel slice appears dissolved, incubate the tube for an additional 5 minutes. (Optional) Add 1 gel volume of isopropanol to the dissolved gel slice. Mix well.

Gel	Tube	Buffer L3 Volume
≤2% agarose	1.7-mL polypropylene	3:1
>2% agarose	5-mL polypropylene	6:1

4. Load. Pipet the dissolved gel piece onto a column inside a Wash Tube. Centrifuge the column at $>12,000 \times g$ for 1 minute. Discard the flow-through and place the column into the Wash Tube.

Note: The column reservoir capacity is 850 μ L. Use 1 column per 400 mg of agarose gel.

- 5. Wash. Add 500 μ L Wash Buffer (W1) containing ethanol to the column. Centrifuge the column at >12,000 \times g for 1 minute. Discard the flow-through and place the column into the Wash Tube. Centrifuge the column at maximum speed for 2–3 minutes. Discard the flow-through.
- **6. Elute.** Place the column into a Recovery Tube. Add 50 μ L Elution Buffer (E1) to the column. Incubate the tube for 1 minute at room temperature. Centrifuge the tube at >12,000 × g for 1 minute.
- 7. **Store.** The elution tube contains the purified DNA. Store the purified DNA at 4°C for immediate use or at -20°C for long-term storage.

Purifying DNA from Gels Using a Vacuum Manifold



Before Starting. Add ethanol to Wash Buffer (W1) according to the label on the bottle. Equilibrate a water bath or heat block to 50°C.

- **1. Excise**. Use a clean, sharp razor blade to excise a minimal area of gel containing the DNA fragment of interest.
- **2. Weigh.** Using a scale sensitive to 0.001 g, weigh the gel slice containing the DNA fragment.
- 3. Solubilize. Add Gel Solubilization Buffer (L3) to the excised gel in a tube as indicated in the table. Incubate the tube at 50°C for 10 minutes (or longer for large gel slices and high concentration gels), and invert the tube every 3 minutes. After the gel slice appears dissolved, incubate the tube for an additional 5 minutes. (Optional) Add 1 gel volume of isopropanol to the dissolved gel slice. Mix well.

Gel	Tube	Buffer L3 Volume
≤2% agarose	1.7-mL polypropylene	3:1
>2% agarose	5-mL polypropylene	6:1

4. Load. Remove the Wash Tube, and attach a column to the vacuum manifold. Pipet the dissolved gel mixture containing the DNA of interest onto the center of the column. Apply vacuum until all of the liquid passes through the column.

Note: The column reservoir capacity is 850 μ L. Use 1 column per 400 mg of agarose gel.

- 5. Wash. Add 500–700 μL Wash Buffer (W1) containing ethanol to the column. Apply vacuum until all of the liquid passes through the column. Place the column into a Wash Tube. Centrifuge the column at maximum speed for 2–3 minutes. Discard the flowthrough.
- 6. Elute. Place the column into a Recovery Tube. Add 50 μ L Elution Buffer (E1) to the column. Incubate the tube for 1 minute at room temperature. Centrifuge the tube at >12,000 × g for 1 minute.
- 7. **Store.** *The elution tube contains the purified DNA.* Store the purified DNA at 4°C for immediate use or at −20°C for long-term storage.

Troubleshooting

Problem	Solution
Low DNA yield	 PCR Purification: Check the amplicon on the gel to verify the PCR product prior to purification. Always mix 1 volume of PCR (50–100 μL) with 4 volumes of Binding Buffer (B2). Add 100% isopropanol to the Binding Buffer (B2) and 96–100% ethanol to Wash Buffer (W1). Add Elution Buffer (E1) to the center of the column. Incubate the tube with Elution Buffer for 1 minute before centrifugation
	 Gel Extraction: Ensure that the correct volume of Gel Solubilization Buffer (L3) is added. Verify that the temperature of water bath or heat block is 50°C. Cut large gel slices into several pieces. Mix the gel slice in the buffer every 3 minutes during the dissolution step. Increase the incubation time for elution to >10 minutes. Note: This kit is not designed to purify supercoiled plasmid DNA.
Low A _{260/230} ratio	 Gel Extraction: Do not get any buffer solution in the cap area of the tube. Add a second wash step with Wash Buffer (W1): After washing the DNA with Wash Buffer (step 5): 1. Add another 500–700 μL Wash Buffer containing ethanol. 2. Centrifuge the tube at 12,000 × g. Discard the flow-through and return the column into the Wash Tube. 3. Centrifuge at the tube at maximum speed for 2–3 minutes.
Enzymatic reactions are inhibited	For all Procedures: To remove Wash Buffer, discard Wash Buffer flow-through from the Wash Tube. Place the column into the Wash Tube and centrifuge the column at $>12,000 \times g$ for 2–3 minutes to completely dry the column.

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

©2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA

Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.invitrogen.com/support or email techsupport@invitrogen.com www.lifetechnologies.com

