

Quantifiler™ Human DNA Quantification Kits

PCR Amplification


For safety and biohazard guidelines, refer to the “Safety” section in the Preface of the *Quantifiler™ Human DNA Quantification Kit* and *Quantifiler™ Y Human Male DNA Quantification Kit User’s Manual* (PN 4344790). For each chemical in **bold** type below, read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product Overview

The Quantifiler™ Human DNA Quantification Kit and the Quantifiler™ Y Human Male DNA Quantification Kit are designed to quantify the total amount of amplifiable human (and higher primate) DNA or human male DNA in a human identity sample.

This quick reference card covers the protocol for the preparation of DNA standards and sample reactions. Refer to the Quantifiler kit’s user’s manual or the appropriate instrument quick reference card for instructions on setting up the plate document, running PCR, and analyzing data.

Workflow



Materials

| Material | Source |
|---|------------------------------------|
| Quantifiler Human DNA Quantification Kit | Applied Biosystems (PN 4343895) |
| Quantifiler Y Human Male DNA Quantification Kit | Applied Biosystems (PN 4343906) |
| Glycogen, 20 mg (1 mL) | Roche Applied Science (PN 901 393) |
| Equipment | Source |
| 96-Well Optical Reaction Plates | Applied Biosystems (PN 4306737) |
| Optical Adhesive Covers Starter Kit | Applied Biosystems (PN 4313663) |
| Optical Tubes | Applied Biosystems (PN 4316567) |
| Optical Caps | Applied Biosystems (PN 4323032) |



Figure 1 The Quantifiler™ Human DNA Quantification Kit and Quantifiler™ Y Human Male DNA Quantification Kit.

Preparing the DNA Quantification Standards

The table below shows a sample dilution series.

IMPORTANT! Applied Biosystems recommends:

- Three-fold dilution series with eight concentration points.
- Minimum input volume of 10 μL DNA for dilutions.

Table 1-1 Standards dilution series - example

| Standard | Concentration (ng/ μL) | Amounts | Minimum Amounts | Dilution Factor |
|----------|------------------------------------|---|---|-----------------|
| Std. 1 | 50.000 | 50 μL [200 ng/ μL stock] + 150 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [200 ng/ μL stock]+ 30 μL T ₁₀ E _{0.1} buffer | 4X |
| Std. 2 | 16.700 | 50 μL [Std. 1] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 1] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 3 | 5.560 | 50 μL [Std. 2] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 2] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 4 | 1.850 | 50 μL [Std. 3] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 3] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 5 | 0.620 | 50 μL [Std. 4] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 4] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 6 | 0.210 | 50 μL [Std. 5] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 5] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 7 | 0.068 | 50 μL [Std. 6] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 6] + 20 μL T ₁₀ E _{0.1} buffer | 3X |
| Std. 8 | 0.023 | 50 μL [Std. 7] + 100 μL T ₁₀ E _{0.1} /glycogen buffer | 10 μL [Std. 7] + 20 μL T ₁₀ E _{0.1} buffer | 3X |

Standards Preparation Procedure

1. Label 8 microcentrifuge tubes: Std. 1 through Std. 8.
2. Dispense the required amount of diluent (T₁₀E_{0.1} Buffer with or without glycogen) to each tube.

Note: T₁₀E_{0.1} buffer recipe:

- 10 mM Tris-HCl (pH 8.0)
- 0.1 mM Na₂EDTA
- 20 $\mu\text{g}/\text{mL}$ glycogen (optional)

3. Prepare Std. 1:
 - a. Vortex the Quantifiler Human DNA Standard 3 to 5 seconds.
 - b. Using a new pipette tip, add the calculated amount of Quantifiler Human DNA Standard to the tube for Std. 1.
 - c. Mix the dilution thoroughly by vortexing for 5 sec.

4. Prepare Std. 2 through 8:
 - a. Using a new pipette tip, add the calculated amount of the prepared standard to the tube for the next standard.
 - b. Mix the standard thoroughly.
 - c. Repeat steps 4a and 4b until you complete the dilution series.

Prepare Reactions

IMPORTANT! While preparing the reactions, keep the 96-well reaction plate in its base and do not place it on the counter.

1. Calculate the volume of each component needed:

| Component | Volume Per Reaction (µL) |
|---|--------------------------|
| Quantifiler Human Primer Mix or Quantifiler Y Human Male Primer Mix | 10.5 |
| Quantifiler PCR Reaction Mix | 12.5 |

2. Prepare the reagents:
 - Thaw the primer mix completely, then vortex 3 to 5 seconds and centrifuge briefly before opening the tube.
 - Swirl the **Quantifiler PCR Reaction Mix** gently before using. Do not vortex.
3. Pipette the required volumes of components into an appropriately sized polypropylene tube.
4. Vortex the PCR mix 3 to 5 seconds, then centrifuge briefly.
5. Dispense 23 µL of the PCR mix into each reaction well.
6. Add 2 µL of sample, standard, or control to the appropriate wells. For plate setup examples, see the *Quantifiler™ Human DNA Quantification Kit and Quantifiler™ Y Human Male DNA Quantification Kit User's Manual*.

IMPORTANT! Applied Biosystems recommends running duplicates of the 8 DNA quantification standards for each assay.

7. Seal the reaction plate with the Optical Adhesive Cover.

IMPORTANT! Make sure the corners are sealed properly.

8. Centrifuge the plate at 3000 rpm for about 20 seconds in a tabletop centrifuge with plate holders to remove any bubbles.
9. Place the compression pad over the Optical Adhesive Cover with the gray side down, the brown side up, and the holes positioned directly over the reaction wells.
10. Run the reactions on the appropriate Applied Biosystems instrument. Refer to the Quantifiler kit's user's manual or the appropriate instrument quick reference card for the protocol and data analysis.



iScience. To achieve accurate, reproducible results, life scientists are taking advantage of advanced analysis systems that unite technology, informatics, and traditional laboratory research. In partnership with our customers, Applied Biosystems provides the innovative products, services, and knowledge resources that make this new, **Integrated Science** possible.

© Copyright 2004, Applied Biosystems. All rights reserved.

For Research, Forensic or Paternity Use Only. Not for use in diagnostic procedures.

The PCR process and the 5' nuclease process are covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd.

Applied Biosystems is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

AB (Design), Applied, iScience, and Quantifiler are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and/or certain other countries. All other trademarks are the sole property of their respective owners.

Printed in USA, 06/2004

