

GeneAmp® 10X PCR Buffer I

Catalog Number N8080006, N8080129, 4379876

Pub. No. 4376113 Rev. C

Contents	Cat. no.	Quantity	Storage conditions
GeneAmp® 10X PCR Buffer I	N8080006	1.5 mL	When stored at -20°C to -5°C in a constant temperature freezer, this product will maintain performance until the control date on the tube or bottle label.
	N8080129	6 x 1.5 mL	
	4379876	2 x 75 mL	



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **www.lifetechnologies.com/support**.

Description

GeneAmp® 10X PCR Buffer I provides the preferred pH and ionic strength for PCR amplification reactions performed with AmpliTaq® DNA Polymerase.

Specifications

Composition: 100 mM Tris-HCl, pH 8.3 (at 25° C); 500 mM KCl; 15 mM MgCl₂; 0.01% gelatin in autoclaved, deionized, ultrafiltered water.

Directions

Before use, thaw frozen components to room temperature and mix well by vortexing.

For 75-mL bottles only

To prepare-aliquots:

- 1. Autoclave 2-mL polypropylene microcentrifuge screwcapped tubes or 1.5-mL polypropylene microcentrifuge tubes at 121°C for 20 minutes.
- 2. Pipet 1.5 mL buffer into 2.0 mL polypropylene microcentrifuge screwcapped tubes or 1.0 mL into 1.5 mL polypropylene microcentrifuge flip-top tubes using aseptic techniques (with gloves and sterile pipette tips), capping each one as it is filled.

Store the tubes upright in a rack in a –20°C constant temperature freezer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures.

For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

DISCLAIMER

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Life Technologies is a Thermo Fisher Scientific brand. ©2014 Thermo Fisher Scientific Inc. All rights reserved. AmpliTaq is a registered trademarks of Roche Molecular Systems, Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

For support visit **lifetechnologies.com/support** or email **techsupport@lifetech.com**

lifetechnologies.com

