

	Package contents	<table border="1"> <tr> <th>Catalog number</th> <th>Size</th> </tr> <tr> <td>13000-012</td> <td>50 rxns</td> </tr> <tr> <td>13000-013</td> <td>200 rxns</td> </tr> <tr> <td>13000-014</td> <td>1000 rxns</td> </tr> </table>	Catalog number	Size	13000-012	50 rxns	13000-013	200 rxns	13000-014	1000 rxns	Kit contents
Catalog number	Size										
13000-012	50 rxns										
13000-013	200 rxns										
13000-014	1000 rxns										
	Storage conditions	<ul style="list-style-type: none"> Store all contents at -20°C. 									
	Required materials	<ul style="list-style-type: none"> Template: cDNA, gDNA, λDNA Forward and reverse gene-specific primers Invitrogen™ E-Gel™ General Purpose Gels, 1.2% (Cat. no. G5018-01) Invitrogen™ TrackIt™ 1 kb Plus DNA Ladder (Cat. no. 10488-085) 0.2 or 0.5-mL nuclease-free microcentrifuge tubes Gel loading buffer 									
	Timing	Varies depending on amplicon length									
	Selection guide	<p>PCR Enzymes and Master Mixes Go online to view related products.</p> <ul style="list-style-type: none"> Platinum™ Hot Start PCR 2X Master Mix contains Platinum™ <i>Taq</i> DNA Polymerase in an optimized PCR buffer with Mg²⁺ and dNTPs. The master mix retains all features of the Platinum™ <i>Taq</i> DNA Polymerase. 									
	Product description	<ul style="list-style-type: none"> Specific binding of the anti-<i>Taq</i> antibody inhibits polymerase activity at ambient temperatures. Activity is restored after the initial denaturation step in PCR cycling at 94°C, providing an automatic “hot start” and offering increased sensitivity, specificity, and yield, while allowing reaction assembly at room temperature. Platinum™ GC Enhancer improves amplification of GC-rich targets. 									
	Important guidelines	Click here for important PCR guidelines.									
	Online resources	Visit our product page for additional information and protocols. For support, visit thermofisher.com/techresources .									

Enzyme characteristics

Hot-start:	Antibody
Length:	Up to 5 kb
Fidelity vs. <i>Taq</i>:	1X
Format:	Master mix

PCR setup

Use the measurements below to prepare your PCR experiment, or enter your own parameters in the column provided.

Component	25-µL rxn	50-µL rxn	Custom	Final conc. in 50-µL rxn
Water, Nuclease-free	to 25 µL	to 50 µL	to µL	—
Platinum™ Hot Start PCR 2X Master Mix	12.5 µL	25 µL	µL	1X*
10 µM forward primer	0.5 µL	1 µL	µL	0.2 µM
10 µM reverse primer	0.5 µL	1 µL	µL	0.2 µM
Template DNA	varies	varies		<500 ng/rxn
Platinum™ GC Enhancer (optional)**	5 µL	10 µL	µL	20%

* 1X master mix contains a final concentration of 1.5 mM MgCl₂ and 0.2 mM of each dNTP.

** Recommended only for targets with >65% GC sequences.

PCR protocol

See page 2 for instructions to prepare and run your PCR experiment.






Optimization strategies

Click here for guidelines to optimize your PCR experiment.

Purchaser notification

Click here for Limited warranty, Disclaimer, and Licensing information.

The example PCR procedure below shows appropriate volumes for a single **50- μ L** reaction using GC-rich template DNA. For multiple reactions, prepare a master mix of common reaction components, then dispense appropriate volumes into each 0.2–0.5 mL PCR tube prior to adding template DNA and primers.

Steps	Action	Procedure details																						
1 	Thaw reagents	Thaw, mix, and briefly centrifuge each component before use. Avoid generating bubbles when mixing the Master Mix.																						
2 	Prepare PCR master mix	<p>Add the following components to each PCR tube.</p> <p>Note: Consider the volumes for all components listed in steps 2 and 3 to determine the correct amount of water required to reach your final reaction volume.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>50-μL rxn</th> <th>Final conc.</th> </tr> </thead> <tbody> <tr> <td>Water, Nuclease-free</td> <td>to 50 μL</td> <td>—</td> </tr> <tr> <td>Platinum™ Hot Start PCR 2X Master Mix</td> <td>25 μL</td> <td>1X</td> </tr> <tr> <td>Platinum™ GC Enhancer (<i>optional</i>)*</td> <td>10 μL</td> <td>20%</td> </tr> </tbody> </table> <p>*For targets with >65% GC sequences.</p> <p>Mix and then briefly centrifuge the components.</p>	Component	50- μ L rxn	Final conc.	Water, Nuclease-free	to 50 μ L	—	Platinum™ Hot Start PCR 2X Master Mix	25 μ L	1X	Platinum™ GC Enhancer (<i>optional</i>)*	10 μ L	20%										
Component	50- μ L rxn	Final conc.																						
Water, Nuclease-free	to 50 μ L	—																						
Platinum™ Hot Start PCR 2X Master Mix	25 μ L	1X																						
Platinum™ GC Enhancer (<i>optional</i>)*	10 μ L	20%																						
3 	Add template DNA and primers	<p>Add your template DNA and primers to each tube for a final reaction volume of 50 μL.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>50-μL rxn</th> <th>Final conc.</th> </tr> </thead> <tbody> <tr> <td>10 μM forward primer</td> <td>1 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 μM reverse primer</td> <td>1 μL</td> <td>0.2 μM</td> </tr> <tr> <td>Template DNA</td> <td>varies</td> <td><500 ng/rxn</td> </tr> </tbody> </table> <p>Cap each tube, mix, and then briefly centrifuge the contents.</p>	Component	50- μ L rxn	Final conc.	10 μ M forward primer	1 μ L	0.2 μ M	10 μ M reverse primer	1 μ L	0.2 μ M	Template DNA	varies	<500 ng/rxn										
Component	50- μ L rxn	Final conc.																						
10 μ M forward primer	1 μ L	0.2 μ M																						
10 μ M reverse primer	1 μ L	0.2 μ M																						
Template DNA	varies	<500 ng/rxn																						
4 	Incubate reactions in a thermal cycler	<table border="1"> <thead> <tr> <th colspan="2">Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td colspan="2">Initial denaturation</td> <td>94°C</td> <td>2 minutes</td> </tr> <tr> <td rowspan="3">25–35 PCR cycles</td> <td>Denature</td> <td>94°C</td> <td>30 seconds</td> </tr> <tr> <td>Anneal</td> <td>~55°C (depending on primer T_m)</td> <td>30 seconds</td> </tr> <tr> <td>Extend</td> <td>72°C</td> <td>1 minute/kb</td> </tr> <tr> <td colspan="2">Hold</td> <td>4°C</td> <td>indefinitely</td> </tr> </tbody> </table>	Step		Temperature	Time	Initial denaturation		94°C	2 minutes	25–35 PCR cycles	Denature	94°C	30 seconds	Anneal	~55°C (depending on primer T _m)	30 seconds	Extend	72°C	1 minute/kb	Hold		4°C	indefinitely
Step		Temperature	Time																					
Initial denaturation		94°C	2 minutes																					
25–35 PCR cycles	Denature	94°C	30 seconds																					
	Anneal	~55°C (depending on primer T _m)	30 seconds																					
	Extend	72°C	1 minute/kb																					
Hold		4°C	indefinitely																					
5 	Add gel loading buffer and analyze with gel electrophoresis	<p>Add gel loading buffer to 10 μL of PCR sample, mix, and briefly centrifuge the contents. Analyze the sample using agarose gel electrophoresis. Use your PCR product immediately in down-stream applications, or store it at –20°C.</p>																						