



Product Datasheet

Product No: PLMM01 Quantity: 100 reactions



Expiry Date

Supplied with: 4 x 625µl 2X Taq Master Mix*

3ml of Nuclease-free Water 1ml of 50mM MgCl₂

Store at -20°C

*2X Tag Master Mix consists of Tag DNA Polymerase (0.05u/µl), 2X Vibuffer A, 0.4mM dNTPs and 3.0mM MgCl₂.



info@vivantechnologies.com

Description:

2X Taq Master Mix is an optimized ready-to-use 2X concentrated DNA amplification mixture containing Taq DNA Polymerase, reaction buffer, dNTPs and MgCl₂. It contains all the components required for routine DNA amplification except template and primers.

Features:

- Saves time and reduces contamination due to reduced number of pipetting steps.
- Stable at 4°C for 6 months, allowing immediate reaction setup without the time-consuming thawing of reagent.
- Suitable for all routine DNA amplification applications.
- Generates mostly 3'dA overhang PCR products which are suitable for TA cloning.

Storage and Stability:

- 2X Taq Master Mix is stable at -20°C for one year or at 4°C for 6 months if properly stored.
- 2X Taq Master Mix is stable for 20 freeze-thaw cycles. To avoid frequent freeze-thaw, keeping small aliquot at -20°C is recommended.
- For daily use, keeping an aliquot at 4°C is recommended.

Quality Control:

All preparations are assayed for contaminating endonuclease, exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.



Amplification of 5kb DNA fragment from lambda DNA using 2X Taq Master Mix in a 50μl reaction mixture.

Lane M: VC 1kb DNA Ladder.

Lane 1 : DNA amplification product generated with 1.25u of

Taq DNA polymerase.

: DNA amplification product generated with 2X Taq

Master Mix (store at -20°C).

Lane 3 : DNA amplification product generated with 2X Tag

Master Mix (after 20 freeze-thaw cycles)

0.7% TAE agarose gel.

72°C for 7 minutes	al Extension
72°C for 30 seconds	ension / 1kb
50 - 68°C for 30 seconds	nnealing
94°C for 2 seconds	naturation
94°C for 2 minutes	naturation

dding additional MgCl₂

higher

may be achieved by

25 - 35 cycles

This protocol may change depending on the template DNA and primers used

Reagent:	Volume	Final Concentration
2X Taq Master Mix	25µI	*1X
MgCl ₂ (50mM)	Refer to Table (A)	**For more than 1.5mM MgCl ₂
Primers (Fwd / Rev)	Variable	0.1 - 1 μM each
DNA Template	Variable	0.02 - 5μg
Nater, nuclease-free	Adjust final volume to 50µl	lume to 50µl

Add the following components in a 0.2ml thin walled PCR tube on ice

Spin down briefly and keep on ice

Gently mix all solutions after thawing.

For 50µl reaction volume:

_
a
able
0
\supset
П
윽
: For more
ō
ർ
≑
thar
٠.
-
ੜ'
₹
<u></u>
≅.
<u>a</u>
≤
g
Ω
20
ö
2
(A): For more than 1.5mM final MgCl ₂ concentration
≓
<u>න</u>
<u>e</u> .
Š

2.5	2.0	1.5	1.0	0.5	Volume of MgCl ₂ (50mM) stock to add into 50µl reaction mixture (µl)
4.0	3.5	3.0	2.5	2.0	Final MgCl ₂ concentration (mM)

Product Use Limitation This product is for research purpose an in vitro use only V 1 V a n t 1 S | www.vivantechnologies.com