

BigDye™ Terminator v3.1 Matrix Standard Kit

310 Genetic Analyzer

Catalog Number 4336948

Pub. No. 4336886 Rev. C

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *310 GeneScan™ Reference Guide* (Pub. No. 4303189). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The BigDye™ Terminator v3.1 Matrix Standard Kit is used to perform spectral calibrations on the 310 Genetic Analyzer. The kit contains a set of four matrix standards. Each matrix standard contains DNA fragments labeled with a different fluorescent dye.

Contents and storage

Contents	Amount	Storage
BigDye™ Terminator v3.1 Matrix Standard (G-Dye 1)	1 tube	Store at 2–8°C, protected from light. ^[1] Do not freeze.
BigDye™ Terminator v3.1 Matrix Standard (A-Dye 2)	1 tube	
BigDye™ Terminator v3.1 Matrix Standard (T-Dye 3)	1 tube	
BigDye™ Terminator v3.1 Matrix Standard (C-Dye 4)	1 tube	

^[1] The kit is stable for six months when stored at 2–8°C.

Guidelines for use

- For more information on the use of matrix standards, see the instrument user guide or *DNA Fragment Analysis by Capillary Electrophoresis User Guide* (Pub. No. 4474504).
- IMPORTANT!** Thoroughly mix the contents of the matrix standard tubes, then briefly centrifuge before use.
- To prepare the matrix standard dilution, combine the appropriate volumes of matrix standard and Hi-Di™ Formamide (Cat. No. 4311320). Dilution volumes vary depending on the specific application and instrument.
- Do not prepare the matrix standard more than 2 hours in advance.
- Do not add size standard to the matrix standard.
- IMPORTANT!** Discard any unused reagent that has been diluted in Hi-Di™ Formamide.

Prepare the standard for the 310 Genetic Analyzer

- Combine the following components for each of the four matrix standards:

Component	Volume
Matrix standard	1 µL
Hi-Di™ Formamide	12 µL
Total volume	13 µL

2. Mix thoroughly, then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
3. To denature the DNA fragments, incubate the mixture at 95°C for 2 minutes. Immediately place the mixture on ice.
4. Load each matrix standard in a separate injection.
5. Open the GeneScan™ injection list in the 310 Genetic Analyzer Data Collection Software.
6. Import the 4-dye sample sheet, then select the appropriate Run Module:
 - Seq POP6 (1 mL) E
 - Seq POP6 Rapid (1 mL) E
 - P4StdSeq (1 mL) E
 - P4StdSeq (1 mL) E
7. Start the run.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4336886

Revision	Date	Description
C	05 November 2018	<ul style="list-style-type: none"> • Update instrumentation • Update licensing, trademarks, general style and format.
B	13 August 2002	Baseline for this revision

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.