

GeneScan™ 500 ROX™ Size Standard

Catalog Number 401734

Pub. no. 4340060 Rev. F

Contents [†]	Size	Storage conditions	
16 ROX™ dye-labeled, single-stranded DNA fragments	2 X 200 µL (800 reactions‡)	2–8°C for 1 year from date of manufacture. Do not freeze.	

[†] With the retirement of the ABI PRISM® 377 DNA Sequencer this kit no longer contains the gel loading buffer. If you need to continue to purchase the loading buffer, it is available as a separate part (Cat. no. 402055).

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **www.lifetechnologies.com/support**.

Product description

GeneScan[™] 500 ROX[™] Size Standard is an internal lane size standard developed for use with fluorescence-based DNA electrophoresis systems from Life Technologies. The use of an internal lane size standard enables automated data analysis and is also essential for achieving high run-to-run precision in sizing DNA fragments by electrophoresis. GeneScan[™] 500 ROX[™] Size Standard is designed for sizing DNA fragments in the 35–500 bp range and provides 16 single-stranded, ROX[™] dye-labeled fragments of 35, 50, 75, 100, 139, 150, 160, 200, 250, 300, 340, 350, 400, 450, 490, and 500 bases. Each of the DNA fragments is labeled with a proprietary fluorophore, which results in a single peak when run under both denaturing and non-denaturing conditions.

Instructions for use

Sample preparation

1. Before use, mix the contents of each tube thoroughly and centrifuge briefly to collect the liquid at the bottom of the tube. Typical loading cocktails are as follows:

	Electrophoresis system			
Component	310	3130/3100 Series	3500 Series	3730 Series
Sample	0.5 μL	0.5 µL	0.5 µL	0.5 μL
GeneScan™ 500 ROX™ Size Standard	0.5 µL	0.25 μL	0.25 μL	0.5 μL
Hi-Di™ Formamide (Cat. no. 4311320)†	9.0 µL	9.25 μL	9.25 μL	9.0 µL

[†] Not included in this kit.

Note: We highly recommend using the ratios of DNA sample (PCR product) and size standard presented in the table as a starting point only. Optimize the ratios if necessary, based on your experimental results. Also, see guidelines in the following section.

IMPORTANT! For HID applications, please follow protocol volumes provided with those application kits.

- 2. Heat the loading cocktail for 3 minutes at 95°C.
- 3. Immediately chill on ice for a few minutes and load samples.

Guidelines



Optimize your analysis based on the following information about use of this standard on capillary electrophoresis instruments:

- The 250-bp peak is sensitive to small temperature variations. The 250-bp fragment should not be used when defining the size standard in GeneMapper[®] Software.
- The 340-bp peak is subject to large temperature variations.
- Fragment analysis primer peaks can often interfere with the detection of the 35-bp peak.

Note: Discard any unused reagent that has been diluted in Hi-Di™ Formamide.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

The information in this guide is subject to change without notice.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

© 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

life technologies™