PRODUCT INFORMATION SHEET

Pub. No. MAN0000826

Rev. A.0



Contents

 Catalog No.
 Size

 10416014
 50 μg



(1) Kit contents




Storage

- Product is shipped at ambient temperature.
- Store at –20°C.




Product description

- The Invitrogen[™] 50 bp DNA Ladder is designed for sizing and quantification of double stranded DNA on 2% agarose gels.
- The 50 bp DNA Ladder consists of 17 individual chromatography-purified DNA fragments ranging in size from 50 bp to 2,500 bp.
- Three reference bands at 2,500 bp, 800 bp, and 350 bp are included for easy orientation.
- The ladder is supplied with 6X TrackIt[™] Cyan/ Orange Loading Buffer for sample DNA.





Required materials





Online resources

- Visit our product pages for additional information and protocols.
- Go online to view related DNA ladders and markers.
- For support, visit thermofisher.com/support.



Important guidelines

- Do not heat the 50 bp DNA Ladder before loading.
- Load the same volume of DNA sample and DNA ladder.
- For quantification, adjust the concentration of the sample to equalize it approximately with the amount of DNA in the nearest band of the ladder.
- For DNA bands visualization with GelRed[™] use gel staining after electrophoresis to avoid aberrant DNA migration.



Guidelines for agarose gel preparation

 Determine the required agarose concentration for your gel based on the size of DNA fragments to be separated.

Fragment size	Recommended agarose gel %		
	1X TAE	1X TBE	
800-10,000	0.8	0.7	
400-8,000	1.0	0.85	
300-7,000	1.2	1.0	

- Prepare agarose in a flask with 2-4 times the volume of the agarose solution.
- Exercise caution when handling microwaved agarose. The solution may become superheated and foam over when agitated.
- Refer to the product insert for UltraPure[™] Agarose for detailed instructions on agarose preparation.
- Guidelines for staining gels
- Troubleshooting
- Limited product warranty and disclaimer details

Prepare DNA ladders and samples for electrophoresis

Step			Action	
1	THE THE PARTY OF T	Cast agarose gel	a. Prepare agarose solution (w/v) for the gel percentage appropriate for separating your DNA fragments.b. Microwave agarose solution.c. Cast agarose gel.	
2		Prepare DNA ladder	a. Thaw, mix and briefly centrifuge each component before use. b. Add the following components to prepare enough ladder for a single 5 mm well Component DNA ladder [1] 1 μL (500 ng) 6X TrackIt™ Cyan/Orange Loading Buffer Water, nuclease free 4 μL [1] Scale components up or down depending upon width of wells. Modify volume by 0.2 μL (0.1 μg of DNA) for each 1 mm of width. c. Mix gently. d. Load DNA ladder on gel.	
3		a. Dilute your sample with 6X TrackIt™ Cyan/Orange Loading Buffer (Cat. no. 10482028): mix 1 volume of loading dye with 5 volumes of the DNA samples. b. Mix gently. c. Load DNA ladder on gel.		
4		Perform electrophoresis	a. Add appropriate amount of UltraPure TAE or UltraPure TBE buffer to chamber. b. Set appropriate voltage and perform electrophoresis of samples. DNA size Voltage Buffer	
5		Stain agarose gel	 a. Incubate gel in staining buffer for 30 minutes. b. Visualize DNA ladder and samples. - Use UV transilluminator to detect DNA bands stained with ethidium bromide. - Use blue light transilluminator to detect DNA bands stained with SYBR™ stains. 	