

Mouse (monoclonal) Anti-Monkey CD3

PRODUCT ANALYSIS SHEET

Catalog Number: APS0301

Lot Number: See product label

Clone: FN-18

Size: 0.5 mg, lyophilized

Formulation: Lyophilized purified immunoglobulin. Lyophilized from 0.5 mL phosphate buffered saline

containing 125 mM trehalose.

Purification: Purified from tissue culture supernatant by ion exchange chromatography. Membrane filtered

 $(0.22 \mu m)$.

Myeloma/Fusion

Partners: BALB/c mouse spleen cells were fused with SP2/O Ag14 mouse myeloma cells.

Immunogen: Rhesus monkey peripheral blood leukocytes.

Isotype: IgG1

Specificity: The antigen recognized by this monoclonal antibody is the rhesus monkey RhT3 antigen

(homologue of the human CD3 antigen) with a molecular mass of 20-25 kDa. This monoclonal antibody reacts with CD3 from *Macaca mulatta* (rhesus monkey), *Macaca arctoides* (stumptail monkey), *Macaca fascicularis* (cynomolgus monkey), *Cercocebus atys* (sooty mangabey), *Cercopithecus aethiops* (African green monkey), and *Papio anubis* (baboon). This antibody does not react with *Callithrix jacobus* (marmoset), nor does this antibody react with any of the

Hominoidea tested (human, chimpanzee, gorilla, and orangutan).

Applications: This antibody is suitable for use in flow cytometry. Use approximately 50 ng to label

1 x 10⁶ cells. In mitogenic studies, use approximately 1 μg/mL. For inhibition of

Con A/PMA induced proliferation, use approximately 10 µg/mL.

Reconstitution: Reconstitute with 0.5 mL sterile, distilled water. Further dilutions should be made in phosphate

buffered saline with 1% BSA.

Storage: Store the lyophilized preparation at 2 to 8°C. Reconstituted antibody should be stored at ≤ -20 °C in

aliquots. Avoid repeated freeze-thaw cycles.

Expiration Date: Expires one year from date of receipt when stored as instructed.

This product is for research use only. Not for use in diagnostic procedures.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PI APS0301

(Rev 11/08) DCC-08-1089

References:

Amara, R.R., Villinger, F., et al. (2001) Control of a mucosal challenge and prevention of AIDS by a multiprotein DNA/MVA vaccine. Science 292(5514):69-74.

Ansari, A.A., et al. (2002) Administration of recombinant rhesus interleukin-12 during acute simian immunodeficiency virus (SIV) infection leads to decreased viral loads associated with prolonged survival in SIVmac251-infected rhesus macaques. J. Virology 76(4):1731-1743.

Bostik, P., et al. (2001) Relative resistance in the development of T cell anergy in CD4(+) T cells from simian immunodeficiency virus disease-resistant sooty mangabeys. J. Immunol. 166(1):506-516.

Brice, G.T., et al. (2000) A novel role for tumor necrosis factor-α in regulating susceptibility of activated CD4+ T cells from human and nonhuman primates for distinct coreceptor using lentiviruses. Journal of Acquired Immune Deficiency Syndromes 24(1):10-22 (cites the use of monoclonal antibody FN-18 coated onto tosyl activated coated beads for cell stimulation).

Chakrabarti, L.A., et al. (2000) Normal T-cell turnover in sooty mangabeys harboring active Simian Immunodeficiency Virus infection. J. Virol. February 2000:11209-11223.

Chakrabarti, L.A., et al. (2000) Age-dependent changes in T-cell homeostasis and SIV load in sooty mangabeys. J. Med. Primatol. 29:158-165.

Croix, D.A., et al. (2000) Effect of myobacterial infection on virus loads and disease progression in Simian Immunodeficiency Virus-infected rhesus monkey. AIDS Research and Human Retroviruses 16(17):1895-1908.

Croix, D.A., et al. (2002) Alterations in T lymphocyte profiles of bronchoalveloar lavage fluid from SIV- and Pneumocystis carinii-coinfected rhesus macaques. AIDS Research and Human Retroviruses 18(5):391-401.

Earl, P.L., et al. (2002) Comparison of vaccine strategies using recombinant env-gag-pol MVA with or without an oligomeric env protein boost in the SHIV rhesus macaque model. Virology 294(2):270-281.

Giavedoni, L.D., et al. (2001) Expression of the interleukin-18 gene from rhesus macaque by the Simian Immunodeficiency Virus does not result in increased viral replication. J. Interferon Cytokine Res. 21:173-180.

Harouse, J.M., et al. (2001) Mucosal transmission and induction of simian AIDS by CCR5-specific simian/human immunodeficiency virus SHIVSF162P3. J. Virol. 75(4):1990-1995.

Keller, E.T., et al. (2001) The impact of chronic estrogen deprivation on immunologic parameters in the ovariectomized rhesus monkey (Macaca mulatta) model of menopause. J. Reprod. Immunol. 50(1):41-55.

Kou, Z.-C., et al. (1998) In vivo effects of a bacterial superantigen on macaque TCR repertoires. J. Immunol. 160:5170-5180.

Lafont, B.A.P., et al. (2000) One-round determination of seven leukocyte subsets in rhesus macaque blood by flow cytometry. Cytometry 41:193-202.

Locher, C.P., et al. (2003) Expression patterns of phenotypic markers on lymphocytes from human immunodeficiency virus type 2-infected baboons. AIDS Research and Human Retroviruses 19(1):31-40.

Ma, S., et al. (1997) Genetic construction and characterization of an anti-monkey CD3 single-chain immunotoxin with a truncated diphtheria toxin. Bioconjug. Chem. 8(5):695-701.

Mohri, H., et al. (1998) Rapid turnover of T lymphocytes in SIV-infected rhesus macagues. Science 279(5354):1223-1227.

Neville, D.M., et al. (1996) A new reagent for the induction of T-cell depletion, anti-CD3-CRM9. J. Immunother. Emphasis Tumor Immunol. 19(2):85-92.

Nooij, F.J., et al. (1986) Differentiation antigens on rhesus monkey lymphocytes. II. Characterization of RhT3, a CD-3-like antigen on T cells. Eur. J. Immunol. 16(8):981-984.

Nooij, F.J., et al. (1986) Polymorphism for RhT3, a CD3-like cell surface antigen, expressed on rhesus monkey T lymphocytes. Immunology 59(4):611-620.

Schramm, B., et al. (2000) Cytopathicity of human immunodeficiency virus type 2 (HIV-2) in human lymphoid tissue is coreceptor dependent and comparable to that of HV-1. J. Virol. Oct.: 9594-9600.

Teranishi, K., et al. (2002) Depletion of anti-gal antibodies in baboons by intravenous therapy with bovine serum albumin conjugated to gal oligosaccharides. Transplantation 73(1):129-139.

Thomas, J.M., et al. (1997) Preclinical studies of allograft tolerance in rhesus monkeys: a novel anti-CD3-immunotoxin given peritransplant with donor bone marrow induces operational tolerance to kidney allografts. Transplantation 64(1):124-135.

Thomas, J.M., et al. (2000) Durable donor-specific T and B cell tolerance in rhesus macaques induced with peritransplantation anti-CD3 immunotoxin and deoxyspergualin. Transplantation 69:2497-2503.

Villinger, F., et al. (2003) Evidence for antibody-mediated enhancement of simian immunodeficiency virus (SIV) Gag antigen processing and cross presentation in SIV-infected rhesus macaques. Journal of Virology 77(1):10-24.

Wilson, L.A., et al. (2000) Identification of SIV env-specific CTL in the jejunal mucosa in vaginally exposed, seronegative rhesus macaques (Macaca mulatta). J. Med. Primatol. 29:173-181.

Xu, H., et al. (2003) Studies investigating pretransplant donor-specific blood transfusion, rapamycin, and the CD154-specific antibody IDEC-131 in a nonhuman primate model of skin allotransplantation. J. Immunol. 170(5):2776-2782.

This product is for research use only. Not for use in diagnostic procedures.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

(Rev 11/08) DCC-08-1089

PI APS0301