

Lot
Expiry Date
Concentration
Supplied with

Store at -  $20^{\circ}$ C

Product No : PL3205 Quantity : 200u

:

: 1u/µl

: 2ml of 10X ViBuffer A 1ml of 10X ViBuffer S

1ml of 50mM MgCl<sub>2</sub>

info@vivantechnologies.com

# Description:

Chromo At *Taq* DNA Polymerase is a complex of specific anti-*Taq* monoclonal antibody with top quality thermostable *Taq* DNA Polymerase for automatic "hot start" amplification, resulting in greatly enhanced amplification specificity, sensitivity and yield. At *Taq* DNA Polymerase catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction in the presence of Mg<sup>2+</sup> and has the 5' to 3' exonuclease activity. The enzyme is supplemented with indicator for ease for visualization of the addition of polymerase to the reaction.

#### Features:

- Ultra pure recombinant protein which is reversibly complex with anti-Taq monoclonal antibody that blocks replication activity of the enzyme at moderate temperatures.
- Carefully selected anti-Taq antibodies have high thermal stability, providing protection againts non-specific primer extension from room temperature to 80°C.
- Formation of complexes between *Taq* DNA Polymerase and an anti-*Taq* antibody forms a basis for automatic "hot start" amplification, which allows for the assembly of amplification reactions at room temperature.
- High stability of the complexes allows for the enormous increase in amplification specificity, sensitivity and yield in comparison to the conventional amplification assembly method.
- Increased specificity as a result of reduced amplification artefacts such as primer-dimer formation and mispriming in mutlplex amplification.

#### **Unit Definition:**

1u is defined as amount of enzyme that required to catalyze the incorporation of 10nmoles of dNTP into acid-insoluble material in 30 minutes at 74°C.

### Reaction Buffer:

### **10X ViBuffer A** (without MgCl<sub>2</sub>):

500mM KCI, 100mM Tris-HCI (pHv.1 at 20°C) and 0.1% Triton™X-100. The buffer is optimized for use with 0.1-0.2mM of each dNTP.

### 10X ViBuffer S:

160mM (NH $_4$ ) $_2$ SO $_4$ , 500mM Tris-HCI (pH 9.2 at 22 $^{\circ}$ C), 17.5mM MgCI $_2$  and 0.1% Triton $^{\text{TM}}$ X-100. The buffer is optimized for use with 0.35mM of each dNTP.

#### Storage Buffer:

20mM Tris-HCl (pH 8.0 at 22°C), 100mM KCl, 0.5% Tween<sup>™</sup> 20, 0.5% Nonidet P-40, 0.1mM EDTA, 1mM DTT and 50% glycerol and color dye.

## **Quality Control:**

All preparation are assayed for contaminating endonuclease, exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.

\* This protocol is subjected to changes depending on the template DNA

# Product Size 8.0 - 20.0kb 5.0 - 8.0kb 0.1 - 5.0kb 2.5 2.0 2.0 2.0 2.0 2.5 2.0 2.0 2.0 2.0

**TABLE** 

(A): RECOMMENDED UNITS

FOR

SPECIFIC VIVANTIS DNA POLYMERASES PER 50

μL REACTION VOLUME

	Implate: Plasmid (0.02 - 2 ng)  Lambda (0.1 - 150 ng)  Genomic (0.05 - 5 µg)  Ultrapure DMSO or formamide  DNA Polymerase							
	DNA Polymerase	Ultrapure DMSO or formamide	ViBuffer (1 X)	dNTP Mix	Product Size			
	Refer to the below Table (A)	-	Α	100 μM	100bp – 5kb			
		3%	А	200 μM	5kb – 8kb			
		3%	S	360 µM	8kb - 20kb			

	_	_	_	_	_	_
Final Extension	Cycles	Extension / 1kb	Annealing*	Denaturation	Denaturation	Product Size
72°C, 7 min	25 - 35	72°C, 30 s	50 - 68°C, 30 s	94°C, 30 s	94°C, 2 min	100bp – 5kb
72°C, 7 min	25 - 35	72°C, 45 s	50 - 68°C, 30 s	94°C, 12s	94°C, 2 min	5kb – 8kb
68°C, 7 min	25 - 35	68°C, 1 min	50 - 68°C, 30 s	94°C, 12s	94°C, 2 min	8kb - 20kb

\*Primer dependent

V i V a n t i S | www.vivantechnologies.com

REACTION MIX (FINAL CONCENTRATION):