

**Chromo Max Taq
DNA Polymerase
(recombinant)**



Lot :
 Expiry Date :
 Concentration : 1u/μl
 Supplied with : 2ml of 10X ViBuffer A
 1ml of 10X ViBuffer S
 1ml of 50mM MgCl₂
 Store at - 20°C

Product No : PL2205
 Quantity : 200u

info@vivantechnologies.com



Description:

Chromo Max Taq DNA Polymerase is a modified and optimized thermostable enzyme blend containing Taq DNA Polymerase, Pfu DNA Polymerase and enhancing factors. It exhibits the 3' to 5' proofreading activity, resulting in considerably higher amplification fidelity than possible with unmodified Taq DNA Polymerase. **Recommended for use in amplification to obtain DNA products up to 20kb.** The enzyme is supplemented with indicators for ease of visualization of the addition of polymerase to the reaction.

Features:

- Ultra pure recombinant protein.
- Excellent for multiplex amplification as it exhibits wider tolerance for Mg²⁺ and salt concentrations, pH, template contaminations and has increased half-life in comparison to unmodified Taq DNA polymerase.
- Improves amplification results with critical templates, such as those containing GC-rich regions, palindromes or multiple repeats.
- Increased amplification product yields and purity.

Unit Definition :

1u is defined as amount of enzyme that required to catalyze the incorporation of 10nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

Reaction Buffer:

10X ViBuffer A (without MgCl₂):

500mM KCl, 100mM Tris-HCl (pH9.1 at 20°C) and 0.1% Triton™X-100. The buffer is optimized for use with 0.1-0.2mM of each dNTP.

10X ViBuffer S:


160mM (NH₄)₂SO₄, 500mM Tris-HCl (pH 9.2 at 22°C), 17.5mM MgCl₂ and 0.1% Triton™X-100. The buffer is optimized for use with 0.35mM of each dNTP.

Storage Buffer:

20mM Tris-HCl (pH 8.0 at 22°C), 100mM KCl, 0.5% Tween™ 20, 0.5% Nonidet P-40, 0.1mM EDTA, 1mM DTT and 50% glycerol.

Quality Control:

All preparation are assayed for contaminating endonuclease, 3'-exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.



Amplification Using Vivantis Max Taq DNA Polymerase

Lane M1 : VC Lambda/ Hind III Marker
 Lane 8kb : 8kb amplification products generated using 0.25mM dNTPs, 2u Vivantis Max Taq DNA Polymerase and 3% formamide.
 Lane 10kb-20kb : 10kb, 12kb, 15kb, and 20kb amplification products generated using 0.36mM dNTPs, 2u Vivantis Max Taq DNA Polymerase and 3% formamide.
 Lane 30kb and 40kb : 30kb and 40kb amplification products generated using 0.36mM dNTPs, 2u Vivantis Max Taq DNA Polymerase and 3% formamide.
 Lane M2 : Lambda DNA (Indicateds 48kb)

SUGGESTED INITIAL PCR CONDITIONS FOR VARIOUS PCR PRODUCT SIZES WITH VIVANTIS DNA POLYMERASE (#PL1201 - 06 / #PL2201 - 06 / #PL3201 - 06 / #PL4201 - 06)
 REACTION MIX (FINAL CONCENTRATION):

| | | | | |
|---|-----------------------------|------------------------------|-----------|------------|
| Primers: 0.2 - 1μM | Product Size | 100bp - 5kb | 5kb - 8kb | 8kb - 20kb |
| Template: Plasmid (0.02 - 2 ng) Lambda (0.1 - 150 ng) Genomic (0.05 - 5 μg) | dNTP Mix | 100 μM | 200 μM | 360 μM |
| | ViBuffer (1 X) | A | A | S |
| | Ultrapure DMSO or formamide | - | 3% | 3% |
| | DNA Polymerase | Refer to the below Table (A) | | |

| | | | |
|-----------------|-----------------|-----------------|-----------------|
| Product Size | 100bp - 5kb | 5kb - 8kb | 8kb - 20kb |
| Denaturation | 94°C, 2 min | 94°C, 2 min | 94°C, 2 min |
| Denaturation | 94°C, 30 s | 94°C, 12s | 94°C, 12s |
| Annealing* | 50 - 68°C, 30 s | 50 - 68°C, 30 s | 50 - 68°C, 30 s |
| Extension / 1kb | 72°C, 30 s | 72°C, 45 s | 68°C, 1 min |
| Cycles | 25 - 35 | 25 - 35 | 25 - 35 |
| Final Extension | 72°C, 7 min | 72°C, 7 min | 68°C, 7 min |

TABLE (A): RECOMMENDED UNITS FOR SPECIFIC VIVANTIS DNA POLYMERASES PER 50 μL REACTION VOLUME:

| Product Size | Taq (#PL1201 - 06) | Max Taq (#PL2201 - 06) | AtTaq (#PL3201 - 06) | AtMax Taq (#PL4201 - 06) |
|--------------|--------------------|------------------------|----------------------|--------------------------|
| 0.1 - 5.0kb | 2.0 | 2.0 | 2.0 | 2.0 |
| 5.0 - 8.0kb | 2.5 | 2.0 | 2.5 | 2.0 |
| 8.0 - 20.0kb | 2.5 | 2.0 | 2.5 | 2.0 |
| +20.0kb | -- | 2.0 | -- | 2.0 |

* This protocol is subjected to changes depending on the template DNA.

*Primer dependent

Product Use Limitation
 This product is for research purposes an *in vitro* use only.