MARROWMAXTM Bone Marrow Medium

RODUCT NEW

Figure 1. Chromosome spread from normal bone marrow cells. Cells were cultured in MARROWMAX[™] Medium for 24 h.

- Optimized for cytogenetic analysis
- High mitotic index
- Excellent chromosomal morphology
- Consistent performance

Analysis of human tumors and hematopoietic cells for diagnosis of malignancies is a rapidly growing area of clinical cytogenetics. For the short-term culture of bone marrow, peripheral blood, and hematopoietic cells required for these analyses, many labs use Giant Cell Tumor (GCT) conditioned media containing a variety of hematopoietic growth factors to supplement serum-containing cultures. However, it is difficult to achieve consistent high levels of analyzable cells with commercially available, or homemade formulations supplemented with GCT.

GIBCO[™] MARROWMAX[™] Bone Marrow Medium is a fully supplemented medium developed specifically to support bone marrow and peripheral blood cell culture for *in vitro* cytogenetic analysis of hematological disease.

MARROWMAX[™] Bone Marrow Medium contains a novel human stromal cell conditioned medium. This conditioned medium is composed of a unique blend of hematopoietic growth factors for optimal cell growth. The medium is manufactured under strict controls ensuring consistent performance and superior chromosomal morphology (*figure 1*). Using MARROWMAX[™] Medium results in cultures with a high mitotic index and an increased number of analyzable cells.

Superior Performance

- Outperforms commercially available media containing GCT (*figure 2, reverse*).
- Higher mitotic index and superior chromosomal morphology.
- Consistent lot-to-lot performance (*figure 3, reverse*).


Convenience

- Complete, ready-to-use medium.
- Fully supplemented with serum, antibiotics and L-glutamine.
- Store either frozen or refrigerated.

Quality Assurance

- Manufactured in compliance with the FDA's Quality System Regulation (cGMP) and the current requirements of BS EN ISO 9001.
- Application tested by an independent certified cytogenetics laboratory using human bone marrow cells.
- Extended shelf life of nine months when stored at -20°C and 60 days stored at 4°C.








Figure 2. Comparison of media for

stimulation of mitotic cells. Cells were cultured in basal medium without conditioned medium (Basal), MARROWMAXTM Medium which contains stromal cell-conditioned medium (MARROWMAXTM), Supplier 1 Medium containing GCT-conditioned medium (GCT 1), and Supplier 2 Medium containing GCT-conditioned medium (GCT 2). Mitotic cells were assayed 24 h after plating. Results are mean ± SEM for N = 10 with up to 30 donors.

Figure 3. MARROWMAXTM Medium consistency. Normal bone marrow mononuclear cells were seeded at 1×10^5 cells/ml in 4 different lots of MARROWMAXTM Medium. Results are mean \pm SEM for N = 10.

Ordering Information

Description	Cat. No.	Size	Quantity	Price
MARROWMAX [™] Bone Marrow Medium**	12260-014	100 ml	1–19 20+	\$74.00 66.00
Related Products				
AMNIOMAX [™] II Complete Medium	11269-016	100 ml	1–19 20+	\$33.50 30.80
AMNIOMAX ^{*-} -C100 Complete Medium (system) — The system contains both the medium (90 ml) and the supplement (15 ml)	27000-025	Set	1–19 20+	\$32.55 29.90
AmnioMAX [™] -C100 Basal Medium	17001-082	90 ml	1–19 20+	\$10.10 9.30
	17001-074	450 ml	1–9 10+	\$40.20 36.95
AMNIOMAX [™] -C100 Supplement	17002-080	15 ml	1–19 20+	\$22.45 20.60
	17002-056	75 ml	1–9 10+	\$95.30 87.75
PB-MAX [™] Karyotyping Medium	10386-035	100 ml	1–19 20+	\$20.60 19.00
	10386-043	500 ml	1–9 10+	\$86.70 79.80
KARYOMAX° COLCEMID°, Solution liquid (10 μ g/ml), in HBSS	15210-040	10 ml	1–19 20+	\$9.90 9.20
KARYOMAX® COLCEMID®, Solution liquid (10 μ g/ml), in PBS	15212-012	10 ml	1–19 20+	\$9.90 9.30
KARYOMAX® Giemsa Stain Stock Solution	10092-013	100 ml	1–19 20+	\$24.15 22.20
Gurr Buffer Tablets (pH 6.8)*	10582-013	$50 \times 1 L$	1–19 20+	\$70.90 65.30
Phytohemagglutinin (M Form) (PHA), lyophilized*	10576-015	10 ml	1-9 10+	\$19.55 17.95

United States Headquarters

Invitrogen Corporation 1600 Faraday Avenue • Carlsbad • CA 92008 Phone: 1 760 603 7229 Toll Free: 1 800 955 6288 • Email: techline@invitrogen.com

European Headquarters Invitrogen Ltd

3 Fountain Drive • Inchinnan Business Park • Paisley PA4 9RF • UK Phone (General Enquires): 0800 5345 5345 • Fax: +44 (0) 141 814 6287 Phone (Free Phone Orders): 0800 269 210 • Email: eurotech@invitrogen.com GIBCO[™] Invitrogen Corporation

For specific country contact information, please see our website: www.invitrogen.com

** This product is subject to Limited Label License No. 31. These products are for *in vitro* diagnostic use and are not intended for human or animal therapeutic use. Uses other than the labled intended use may be a violation of federal law. * The noted products are for laboratory research use only and not for diagnostic use. The safety and efficacy of these products in diagnostic or other clinical uses has not been established. MARROWMAX^{III} Bone Marrow medium is subject to Limited Label License 31. COLCEMED[®] is a registered trademark of CIBA-GEIGY Corporation. © 2001 Invitogen Corporation PAS01-024MS Part No. E-08A00-0703

