Lipofectamine[®] LTX & PLUS[™] Reagent

S	Package Contents	Catalog Number A12621 15338-030 15338-100 15338-500* 	Size: LTX/PLUS™ 0.1 mL/40 μL 0.3 mL/0.25 mL 1.0 mL/0.85 mL 15 mL/-	*PLUS™ Reagent (Cat. no. 11514-015) available separately.			
	Storage Conditions	 Store at 4°C (do not freeze). 					
	Required Materials	 Plasmid DNA (0.5–5 µg/µL stock) Opti-MEM[®] Reduced Serum Medium Eppendorf tubes 					
	Timing	Preparation: 10 minutes Incubation: 5 minutes Final Incubation: 1-3 days					
R	Selection Guide	Lipofectamine [®] Reagents Go online to view related products.					
Ç	Product Description	 Lipofectamine[®] LTX and PLUSTM Reagent are proprietary formulations for transfecting nucleic acids into a wide range of eukaryotic cells. 					
	Important Guidelines	 DNA-Lipofectamine[®] LTX and PLUS[™] complexes must be made in serum-free medium such as Opti-MEM[®] Reduced Serum Medium and can be added directly to cells in culture medium, in the presence or absence of serum/antibiotic. 					
		 It is not necessary to remove complexes or change/add medium after transfection. 					
		 The amount of Lipofectamine[®] LTX Reagent required for successful transfection varies depending on the cell type and passage number. Start any new transfection by testing the recommended four concentrations of Lipofectamine[®] LTX Reagent to determine an optimum amount. 					
	Online	Visit our product page f information and protoc					

Resources visit www.lifetechnologies.com/support.



For Research Use Only. Not for use in diagnostic procedures.

Protocol Outline

- A. Plate cells so they will be 70–90% confluent at the time of transfection.
- B. Prepare plasmid DNA-lipid complexes.
- C. Add DNA-lipid complexes to cells.

Lipofectamine® LTX DNA Transfection Reagent Protocol

See page 2 to view a typical plasmid transfection procedure.

Transfection Amounts

Component	96-well	24-well	6-well
DNA per well	100 ng	500 ng	2500 ng
PLUS™ Reagent per well	0.1 µL	0.5 µL	2.5 μL
Lipofectamine [®] LTX Reagent per well	0.2–0.5 μL	1–2.5 µL	5–12.5 μL

Scaling Up or Down Transfections

① Limited Product Warranty and Disclaimer Details



Lipofectamine® LTX DNA Transfection Reagents Protocol

Transfect cells according to the following table. Use the indicated volume of DNA and PLUSTM Reagent with each of the four volumes of Lipofectamine[®] LTX. Each reaction mix is sufficient for triplicate (96-well), duplicate (24-well), and single well (6-well) transfections, and accounts for pipetting variations. For even less toxicity, reduce the amount of DNA-lipid complex to the cells, or reduce the amount of DNA used to make complexes.

Timeline		Timeline	Steps	F	Procedure Details			
Day 0			Seed cells to be 70-90% confluent at transfection	Component	96-well	24-well	6-well	
	1			Adherent cells	$1 - 4 \times 10^{4}$	$0.5 - 2 \times 10^{5}$	$0.25 - 1 \times 10^{6}$	
			Dilute Lipofectamine® LTX Reagent in Opti-MEM® Medium	Opti-MEM [®] Medium	$25~\mu L \times 4$	$50 \ \mu L \times 4$	$150 \ \mu L \times 4$	
	2			Lipofectamine® LTX Reagent	1, 1.5, 2, 2.5 μL	2, 3, 4, 5 µL	6, 9, 12, 15 μL	
			Dilute DNA in Opti-MEM® Medium, then add PLUS™ Reagent	Opti-MEM [®] Medium	125 µL	250 μL	700 µL	
	_	200		DNA (0.5–5 μg/μL)	2.5 µg	5 µg	14 µg	
	3			PLUS TM Reagent	2.5 µL	5 µL	14 µL	
Day 1		0	Add diluted DNA to diluted Lipofectamine® LTX Reagent (1:1 ratio)	Diluted DNA (with PLUS™ Reagent) Total	25 µL	50 µL	150 µL	
	4			Diluted Lipofectamine [®] LTX Reagent	25 µL	50 µL	150 µL	
Õ	5		Incubate	Incubate for 5 minutes at room temperature.				
	5	5	Add DNA-lipid complex to cells	Component	96-well	24-well	6-well	
				DNA-lipid complex per well	10 µL	50 µL	250 µL	
				Final DNA used per well	100 ng	500 ng	2500 ng	
	6			Final Lipofectamine [®] LTX Reagent used per well	0.2–0.5 μL	1.0 –2 .5 μL	5.0–12.5 μL	
Day 2-4	7		Visualize/analyze transfected cells	Incubate cells for 1–3 days at 37°C. Then, analyze transfected cells.				