




# Micromanipulator System NT-88-V3 Series






# ***Fusion of Utmost Precision***

ICSI (Intracytoplasmic sperm injection) and IMSI (Intracytoplasmic morphologically selected sperm injection), IVF (In vitro fertilization) technologies that have received a lot of attention in recent years as a remedy for male infertility, require a highly precise manipulator to respond to delicate movement of the operator's hands.

Stable, smooth manipulation is also required for injection to the stem cell, transgenics and electrophysiology.

The NT-88-V3 series, jointly developed by Nikon and Narishige, ensures easy operation and an unsurpassed level of precision. When used in combination with Nikon's high-performance microscopes, it provides strong support for cutting-edge biomedical research.

- Size is greatly reduced as the 3D motor-driven coarse control manipulator and 3D hydraulic fine micromanipulator are combined in a single unit.
- Assembly of the micromanipulator is fast and easy due to the extremely stable, one-piece mounting adapter.
- Wide varieties of microinjectors for holding or injecting specimens depending on use.





# and Fine Mechanism



# Oil-hydraulic Micromanipulators

## 3D Hydraulic Micromanipulator NT-88-V3

Combined 3D motor-driven coarse control manipulator and the 3D hydraulic fine micromanipulator realize remarkably compact and stable design.

- Less than half the size of a conventional model—minimal length tubes and cords—the manipulator fits easily into the incubator.
- Closer-to-the-center settings are possible, while reduced distance from manipulator mount to microelectrode tip increases stability.
- X-axis control knobs are located symmetrically on opposite sides of the joystick to eliminate left side/right side positioning problems.
- 10mm X-Y operating range provides versatility in manipulation at lower magnifications.
- A “resume” feature returns the pipette to the original position after the pipette has been raised for quicker, easier Petri dish changes.
- The universal joint is outfitted with an angle gauge, which ensures accurate-angle operation. The forward-mounted return mechanism facilitates Petri dish changing.
- Permits one-touch installation on an inverted microscope.



4


## 3D Hydraulic Micromanipulator (Single-hand type) NT-88-V3SH

Suitable for experiments that do not require the holding side's pipette, such as the microinjection or microdissection to an adherent cell



## 3D Hydraulic Micromanipulator (Manual Single-hand type) NT-88-V3MSH

A manual 3D coarse positioning manipulator allows you to quickly move the micro-pipette to the ideal position. And the compact design allows you to position it closer to you with the control knob in easy reach.



		NT-88-V3	NT-88-V3SH	NT-88-V3MSH
<b>Hydraulic fine micromanipulator</b>				
Working range	Control knob operation		Max. 10mm	
	Joystick operation		Max. 2mm	
Displacement by control knob rotation	Full turn of knob		250µm	
	Minimum graduation		2µm	
<b>Coarse manipulator</b>				
Working range		Max. 22mm (motorized)		Max. 30mm (manual)

## Universal Joint UJ-V3

Universal Joint UJ-V3 allows the pipette holder to be moved in a single direction. The angle can be adjusted for both vertical and horizontal directions via two separate control knobs. The forward-mounted return mechanism facilitates changing the petri dish while the built-in angle gauge ensures accurate-angle operation. Moreover, because the UJ-V3 is mountable from above, unlike the conventional type, attachment and removal is much easier.



## Pneumatic Microinjector IM-11-2\*

The IM-11-2 allows suction and injection of a microsample with precision and operability comparable to oil hydraulic injectors. It is suitable for both injecting and holding during ICSI. Air bubbles do not enter the micropipette when the pneumatic microinjector is fitted.

Movement range: 40mm

\* In the U.S.A., use IM-11-2A for clinical use.

## Oil-type Microinjector IM-9B Pneumatic Injector IM-9C

The IM-9B, which is suitable for injection, and the IM-9C, which is suitable for holding, provide smoother and more stable maneuverability by applying a new mechanism to minimize backlashes. Also, their 53mm working distance, large-volume syringe and just-the-right-size control knob with graduated scales all contribute to easy control of the injection. The IM-9B has a glass syringe for easy confirmation of internal air bubbles while its new tube connector with a three-way stopcock makes for easier filling of oil without generating air bubbles.



	IM-9B	IM-9C
<b>Movement distance</b>		
Movement range	53mm	53mm
Full rotation of knob	500µm	6mm
<b>Control value</b>		
Full rotation of knob	Approx. 10µl	Approx. 480µl
Syringe	1,060µl glass type	4,240µl metal type

### Pipette Holder Clamper HIK-5



The HIK-5 is used to keep the pipette holder at the desired position, thereby simplifying pipette exchange. The magnetic clamper is attached to a bench or iron plate. Alternatively, it can be attached to the included round metal plate, which can be affixed at the desired location to a microscope stage using double-sided adhesive tape.

### Microscope Adapters for NT-88-V3

Various adapters are available for use with Nikon microscopes. As both adapters and manipulators are compact, attaching them is very easy.



**NS-TEV3**  
Mounting Adapter to use with inverted microscopes Ti and TE2000



**NS-T1**  
Mounting Adapter to use with inverted microscope TS100



# Micromanipulators for electro-physiology applications

## 3D Water-Hydraulic Micromanipulator MHW-3

The MHW-3 is a hydraulic micromanipulator for patch clamping that reduces unwanted drift due to temperature change. A cartridge with a one-to-five moving ratio enables fine control manipulation equivalent to an oil hydraulic manipulator. The drum type control unit enables high linear repeatability.

Movement range	Fine 2mm (X, Y and Z axes), coarse 30mm (X, Y and Z axes)
Full rotation of knob	50 $\mu$ m
Minimum graduation	Fine 0.2 $\mu$ m



MHW-3 with NN-L3 adapter, configured with patch amplifier

## 3D Motorized Micromanipulator MM-80

The MM-80 has been developed to satisfy the competing goals of "stability" and "compactness" that are essential for research work in electrophysiology. The motorized micromanipulator is constructed in a simplified design with its controller employing simple-design push button switches. The drive unit embodies user-friendly convenient features which help reduce stress during experiments.




Movement range	Fine 20mm (X, Y and Z axes), coarse 30mm (X, Y and Z axes)	
Driving speed (Motorized fine drive)	Minimum	Approx. 10 $\mu$ m/s*
	Maximum	Approx. 150 $\mu$ m/s*

\* DC motor drive, theoretical value

## Manual 3D Micromanipulators NMN-21/25

The NMN-21 and -25 are designed for drift-free operation, making them ideal for delicate, patch-clamp research. The unit has no hydraulic system nor spring mechanism and its unique double-slider structure minimizes vibrations transmitted from the operator's hand to the tip of the pipette via the control knob. The unit's motor-less and actuator-less design also eliminates electromagnetic disturbances. The NMN-25 features a compact, space-saving design without a scale.

	NMN-21	NMN-25
Movement range	Fine 6mm (X, Y and Z axes), coarse 15mm (X, Y and Z axes)	Fine 6mm (X and Z axes), coarse 15mm (X and Z axes), combined coarse and fine working range 15mm (Y axis)
Full rotation of knob	Fine 250 $\mu$ m, coarse approx. 4mm	Fine 250 $\mu$ m (X and Z axes), approx. 250 $\mu$ m (Y axis)
Minimum graduation	Fine 1 $\mu$ m	—



NMN-21




NMN-25

# Recommended combination for various applications

## For ICSI (Intracytoplasmic Sperm Injection), ES cell injection and nucleus transfer

The microinjection system is dedicated to floating cell applications. The inner diameter of the pipette tip is more than 5µm. For the injection of the sperm, embryo-stem cell and nucleus, a manual injector is suitable as it accurately transmits the movement of user's hand via the handle.

Product	Unit	Model name
<b>Components</b>		
Inverted microscope	1	Ti-U NAMC set
3D hydraulic micromanipulator	1	NT-88-V3
Microinjector (injection side)	1	IM-9B
Oil-type microinjector (injection side)	1	IM-9C
Thermal plate warmer: ThermoPlate TP series	1	TP-108R05
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Pipette grinder	1	EG-400
Microforge	1	MF-900
Glass capillaries	1	G-100



## For injection into adherent cells

The microinjection system for adherent cells incorporates a one-axis oil hydraulic micromanipulator—it moves in the direction of a pipette—only on the injection side, since it does not require an injector on the holding side. A motorized injector that can accurately control the injection amount of DNA and fluorescence dyes is suitable.

Product	Unit	Model name
<b>Components</b>		
Inverted microscope	1	Ti-U DIC set
3D hydraulic micromanipulator (single-hand type)	1	NT-88-V3SH
Mounting adapter	1	NS-TEV3
1D hydraulic micromanipulator	1	MMO-220A
Motorized microinjector	1	IM-300
Thermal plate warmer: ThermoPlate TP series	1	TP-108R05
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Glass capillaries	1	GD-1



Product	Unit	Model name
<b>Components</b>		
Stereoscopic microscope	1	SMZ18 FL set
Mechanical micromanipulator	1	MN-153
Mounting adapter	1	NN-V, NR-2
Motorized microinjector	1	IM-300
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Glass capillaries	1	GD-1



## Single-cell patch clamping with inverted microscope

This system is suitable for single-cell patch clamp experiments using an inverted microscope. Two units, one each for recording and stimulating, are included. A water-hydraulic micromanipulator effective for general manipulations is suitable.

Product	Unit	Model name
<b>Components</b>		
Inverted microscope	1	Ti-U DIC set
Mounting adapter	1	NN-L3
3D water-hydraulic micromanipulator	2	MHW-3
Patch amplifier adapter*	1	AP-12A/L/N or AP-13-2
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Microforge	1	MF-830
Glass capillaries	1	GD-1.2
Glass capillaries	1	GD-1.5

\* Please select an adapter that suits the brand of patch amplifier used.



## Single-cell patch clamping with upright microscope

This system is suitable for single-cell patch clamp experiments using an upright microscope. Two units, one each for recording and stimulating, are included. A water-hydraulic micromanipulator effective for general manipulations is suitable.

Product	Unit	Model name
<b>Components</b>		
Upright microscope	1	FN1-FL-DIC set
Mounting adapter	1	NN-R
3D water-hydraulic micromanipulator	2	MHW-3
Patch amplifier adapter*	1	AP-12A/L/N or AP-13-2
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Microforge	1	MF-830
Glass capillaries	1	GD-1.2
Glass capillaries	1	GD-1.5

\* Please select an adapter that suits the brand of patch amplifier used.





## Manual patch clamping

This system is suitable for patch clamp experiments on slice cells using an upright microscope. The two stages that incorporate new functions were developed especially for the FN1. One is independent from the microscope and the other, a manual XY stage, enables the microscope's observation field to be changed. Though three water-hydraulic micromanipulators are used in the picture, any number can be selected.



Product	Unit	Model name
<b>Components</b>		
Upright microscope	1	FN1-FL-DIC set
Stage	1	ITS-FN1
3D water-hydraulic micromanipulator	3	MHW-3
Patch amplifier adapter*	2**	AP-12A/L/N, AP-13-2
<b>Options for making pipettes</b>		
Micropipette puller	1	PC-10
Microforge	1	MF-830
Glass capillaries	1	GD-1.2
Glass capillaries	1	GD-1.5

\* Please select an adapter that suits the brand of patch amplifier used.

\*\* Any one of the following combinations is possible: AP-12A/L/N x2, AP-13-2 x1, AP-12A/L/N x1 + AP-13-2 x2

## Compatibility of patch amplifier adapters

Headstage maker and model name		Manipulator/Patch amplifier adapter				
		MHW-3	MM-80	NMN-21 NMN-25	MWS-1B	MHW-4
AXON	HS-2, HS-4, VG-2, VG-2A, CV-4, CV-5, CV-202 (CV-201AU)	AP-12A			AP-14A	
	CV-203BU, CV-7A/B, HS-9A	AP-13-2			AP-13-3	
HEKA	EPC-7, EPC-8, EPC-9	AP-12L			AP-14L	
	EPC-10	AP-13-2			AP-13-3	
Nihon Kohden	JZ-230J, JZ-240J, JZ-245J	AP-12N			AP-14N	



Patch amplifier adapter

# Accessories

## Micropipette Puller PC-10

The PC-10 is used to pull the glass capillary vertically, utilizing the gravitational force of its own weight. It has two modes: a single pull, pulling capillary at one stretch, and a double pull in which the setting is changed in mid-process. The versatile vertical pull type is suitable for both injection purposes and patch clamping.



## Micropipette Puller PN-31

The PN-31 uses a pulling method by which a glass capillary is fine-drawn using electromagnetic force, enabling the production of long, sharp microneedles. The pulling tension can be adjusted over a wide range, from weak tension that is required at the beginning of pulling to strong tension that is required to make a long, thin microneedle.



## Microforge MF-900

The MF-900 is designed to produce injection and holding pipettes. In addition to the temperature and illumination, the position of the heating element can be adjusted.



## Microforge MF-830

A powerful microscope with 525x total magnification is incorporated to make electrodes with tips less than 2µm for patch clamping.



## Pipette Grinder EG-400

The grinder with minimized irregular movements of the grinding plane and the microscope are combined. Precise grinding while confirming needle contact with grinding plane is facilitated.



## Glass Capillaries G-100, G-1, GD-1/1.2/1.5

G-100 is a thin-walled glass capillary ideal for ICSI, while G-1 is for general purpose.

The GD series is commonly known as double tubing. These glass capillaries contain an internal glass fiber (approx. 100 µm).

The glass capillaries have been prewashed in an ultrasonic washing machine.



### Thermal plate warmer ThermoPlate TP series

A temperature-controllable stage ring with a glass-heating plate ensures more accurate and reliable thermal control of specimens. The temperature can be set at between room temperature +5°C and 50°C in 0.1°C increments. A sterilized sensor allows measurement of the actual temperature of dish contents. Management software and continuous current control provide solutions to a wide range of requirements.

Manufactured by Tokai Hit Co., Ltd.



For inverted microscope




For stereoscopic microscope

### Simple Head Holder Plate (for Mice) MAG-1

After affixing the chamber frame to the head of a mouse, the MAG-1 allows observations/experiments under the microscope or other experimental instruments. Designed to be compact in size, the MAG-1 can be placed under a two-photon microscope or a stereomicroscope.



MAG-1



\* Parts for conventional three-point fixation can be customized.



Configuration with Stereo Microscope SMZ18



Photos courtesy of: Masumi Hirabayashi, Associate professor, National Institute for Physiological Sciences

Specimen 1: Injection of rat's ES cell into 8-cell embryo

Specimen 2: Microinsemination of rat egg

Micromanipulators and related equipment other than microscopes are manufactured by:



**NARISHIGE CO., LTD.**

27-9, Minamikarasuyama 4-chome, Setagaya-ku, Tokyo 157-0062, Japan

Tel: +81(0)3 3308 8233 Fax: +81(0)3 3308 2005

e-mail: sales@narishige.co.jp

http://narishige-group.com

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. September 2013  
©1997-2013 NIKON CORPORATION

 <b>WARNING</b>	TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.
---	--

N.B. Export of the products\* in this brochure is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedure shall be required in case of export from Japan. \*Products: Hardware and its technical information (including software)



**NIKON CORPORATION**

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan  
phone: +81-3-6433-3705 fax: +81-3-6433-3785  
<http://www.nikon.com/instruments/>

**NIKON INSTRUMENTS INC.**

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A.  
phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only)  
fax: +1-631-547-0306  
<http://www.nikoninstruments.com/>

**NIKON INSTRUMENTS EUROPE B.V.**

Tripolis 100, Burgerveshuispad 101, 1076 ER Amsterdam, The Netherlands  
phone: +31-20-7099-000 fax: +31-20-7099-298  
<http://www.nikoninstruments.eu/>

**NIKON INSTRUMENTS (SHANGHAI) CO., LTD.**

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060  
(Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026  
(Guangzhou branch) phone: +86-20-3882-0552 fax: +86-20-3882-0580

**NIKON SINGAPORE PTE LTD**

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

**NIKON INSTRUMENTS KOREA CO., LTD.**

KOREA phone: +82-2-2186-8400 fax: +82-2-555-4415

**NIKON CANADA INC.**

CANADA phone: +1-905-602-9676 fax: +1-905-602-9953

**NIKON FRANCE S.A.S.**

FRANCE phone: +33-1-4516-45-16 fax: +33-1-4516-45-55

**NIKON GMBH**

GERMANY phone: +49-211-941-42-20 fax: +49-211-941-43-22

**NIKON INSTRUMENTS S.p.A.**

ITALY phone: +39-55-300-96-01 fax: +39-55-30-09-93

**NIKON AG**

SWITZERLAND phone: +41-43-277-28-67 fax: +41-43-277-28-61



ISO 9001 Certified  
for NIKON CORPORATION  
Microscope Solutions Business Unit  
Industrial Metrology Business Unit



ISO 14001 Certified  
for NIKON CORPORATION

