# E-Gel<sup>™</sup> 1 Kb Plus DNA Ladder


**PRODUCT INFORMATION SHEET** 

Pub. No. MAN0000772

K	Contents	<b>Catalog No.</b> 10488090	<b>Amount</b> 100 applications	() Kit contents
	Storage	<ul> <li>Product is shipped at ambient temperature.</li> <li>Store at room temperature or at 4°C for up to 6 months, or at -20°C for long term storage.</li> </ul>		


#### **Product description**

- The Invitrogen<sup>™</sup> E-Gel<sup>™</sup> 1 Kb Plus DNA Ladder is designed for sizing and quantification of double stranded DNA on 1.2% E-Gel<sup>™</sup> agarose gels.
- The E-Gel<sup>™</sup> 1 Kb Plus DNA Ladder consists of 18 individual chromatography-purified DNA fragments ranging in size from 100 bp to 15,000 bp.
- A reference band at 1,500 bp is included for easy orientation.
- The ladder is supplied with 1X E-Gel<sup>™</sup> Sample Loading Buffer for sample DNA.



**Rev.** A.0

 Visit our product pages for additional information and protocols.



- Go online to view related DNA ladders and markers.
- For support, visit thermofisher.com/support.

### **Required materials**

- E-Gel<sup>™</sup> EX or E-Gel<sup>™</sup> Agarose Gel with SYBR<sup>™</sup> Safe (See **Choosing the right** DNA ladder for your E-Gel<sup>™</sup> agarose gel)
- TE Buffer (Cat. No. AM9858)
- Ultrapure<sup>™</sup> DNase/RNase-Free Distillated Water (Cat. No. 10977023)



#### Important guidelines

- Do not heat the E-Gel<sup>™</sup> 1 Kb Plus DNA Ladder before loading.
- Load the same volume of DNA sample and DNA ladder.
- For quantification, adjust the concentration of the sample to equalize it approximately with the amount of DNA in the nearest band of the ladder.
- Dilute sample DNA in TE buffer to avoid degradation of DNA sample.
- ⑦ Choosing the right DNA ladder for your E-Gel™ agarose gel
- Troubleshooting
- U Limited product warranty and disclaimer details





invitrogen

## Prepare DNA ladders and samples for electrophoresis

This protocol provides a brief description of how to use the DNA ladder with E-Gel<sup>™</sup> agarose gels. For detailed instructions on using specific types of E-Gel<sup>™</sup> agarose gels, go to thermofisher.com or contact Technical Support.

Step			Action		
1		Prepare DNA ladder	<ul> <li>a. Thaw, mix and briefly centrifuge DNA ladder before use.</li> <li>b. Prepare DNA ladder.</li> <li>For E-Gel<sup>™</sup> EX Agarose Gels, mix 2 μL of DNA ladder with 18 μL of water.</li> <li>For E-Gel<sup>™</sup> Agarose Gels, mix and use the ladder without dilution.</li> <li>For E-Gel<sup>™</sup> 48 Agarose Gels, mix 2 μL of DNA ladder with 13 μL of water.</li> </ul>		
2		Prepare samples	a. Dilute your sample 2- to 10-fold with TE Buffer (Cat. No. AM9858), 1X E-Gel <sup>™</sup> Sample Loading Buffer (Cat No. 10482055), or water. b. Mix gently.		
3		Load samples and DNA ladders	<ul> <li>a. Load DNA ladders and DNA samples into the appropriate wells of the E-Gel<sup>™</sup> agarose gel.</li> <li>Add 20 μL for E-Gel<sup>™</sup> and E-Gel<sup>™</sup> EX Agarose Gels.</li> <li>Add 15 μL for E-Gel<sup>™</sup> 48 Agarose Gels.</li> <li>b. Add water to any empty wells, so that all wells contain an equal volume of liquid.</li> </ul>		
		Perform electrophoresis	a. Choose the appropriate E-Gel <sup>™</sup> run protocol for your gel type based on the electrophoresis device being used.		
4			Gel type Program Recommended run time		
			E-Gel <sup>™</sup> Power Snap Electrophoresis Device (Cat. No. G8100)		
			E-Gel <sup>™</sup> EX Agarose Gel (1%) E-Gel EX 4 1-2% 15 min (20 min max)		
			E-Gel <sup>™</sup> Agarose Gel (0.8%, 1.2%, 2%) E-Gel 0.8-2% 26 min (40 min max)		
			E-Gel <sup>™</sup> E-Base <sup>™</sup> Device		
			E-Gel <sup>™</sup> 48 Agarose Gel (4%) EG 20 min		
			b. Run the program to start electrophoresis.		
		Visualize agarose gel	Visualize DNA ladder and samples.		
5			<ul> <li>Use the E-Gel<sup>™</sup> Power Snap Camera (Cat. No. G8200), E-Gel<sup>™</sup> Imager (Cat. No. 466612), or other blue light imager to detect DNA bands stained with SYBR<sup>™</sup> stains.</li> </ul>		
			<ul> <li>UV transilluminator to detect DNA bands stained with ethidium bromide.</li> </ul>		

For support, visit thermofisher.com/support.

**Thermo Fisher** 

SCIENTIFIC