

Ion ReproSeq™ PGS Kits – Ion S5™ /Ion GeneStudio™ S5 Systems

Catalog Numbers A34899, A34900, A34901

Pub. No. MAN0016714 Rev. B.0

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *Ion ReproSeq™ PGS Kits – Ion S5™ /Ion GeneStudio™ S5 Systems User Guide* (Pub. No. MAN0016712). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The Ion ReproSeq™ PGS Kits – Ion S5™/Ion GeneStudio™ S5 Systems provide reagents and materials for high-throughput whole-genome amplification and sequencing to detect chromosomal aneuploidies, chromosome arm events (>48 Mbp), and copy number variations starting from a single cell. The kits are used with the Ion Chef™ System, and the Ion S5™/Ion S5™ XL System or one of the Ion GeneStudio™ S5 Systems, which enable automated template preparation, chip loading, and sequencing to minimize hands-on time.

The Ion SingleSeq™ Kit, a component of each Ion ReproSeq™ PGS Kit, includes reagents to extract, amplify, and barcode genomic DNA for preparing libraries from up to 24 or 96 samples. The Ion ReproSeq™ PGS Kits also include reagents to prepare templates and load four chips on the Ion Chef™ Instrument, and sequence the chips on an Ion S5™/Ion S5™ XL Sequencer or one of the Ion GeneStudio™ S5 Series Sequencers. Depending on the kit that is used, up to 384 samples can be sequenced per kit.

Kit configurations

Ion ReproSeq™ PGS Kit	Amount per kit	Samples per kit
Ion ReproSeq™ PGS Kit with Ion 510™ Chips (16 samples/run); Cat. No. A34899		
Ion SingleSeq™ Kit (24 reactions)	3 boxes	64
Ion S5™ ExT Chef Supplies	4 boxes	
Ion S5™ Chef Solutions	4 cartridges	
Ion S5™ ExT Chef PGS Reagents	4 cartridges	
Ion S5™ ExT Sequencing Solutions	1 box	
Ion S5™ ExT Sequencing Reagents	4 cartridges	
Ion 510™ Chip Kit	4 chips	
Ion ReproSeq™ PGS Kit with Ion 520™ Chips (24 samples/run); Cat. No. A34900		
Ion SingleSeq™ 96 Kit (96 reactions)	1 box	96
Ion S5™ ExT Chef Supplies	4 boxes	
Ion S5™ Chef Solutions	4 cartridges	
Ion S5™ ExT Chef PGS Reagents	4 cartridges	
Ion S5™ ExT Sequencing Solutions	1 box	
Ion S5™ ExT Sequencing Reagents	4 cartridges	
Ion 520™ Chip Kit	4 chips	

Ion ReproSeq™ PGS Kit	Amount per kit	Samples per kit
Ion ReproSeq™ PGS Kit with Ion 530™ Chips (96 samples/run); Cat. No. A34901		
Ion SingleSeq™ 96 Kit (96 reactions)	4 boxes	384
Ion S5™ ExT Chef Supplies	4 boxes	
Ion S5™ Chef Solutions	4 cartridges	
Ion S5™ ExT Chef PGS Reagents	4 cartridges	
Ion S5™ ExT Sequencing Solutions	1 box	
Ion S5™ ExT Sequencing Reagents	4 cartridges	
Ion 530™ Chip Kit	4 chips	

Kit contents and storage

Ion SingleSeq™ Kit

Component	Amount		Storage
	Ion SingleSeq™ Kit	Ion SingleSeq™ 96 Kit	
Cell Extraction Buffer (green cap)	120 µL	480 µL	-30°C to -10°C
Extraction Enzyme Dilution Buffer (violet cap)	115.2 µL	460.8 µL	
Cell Extraction Enzyme (yellow cap)	4.8 µL	19.2 µL	
Pre-Amplification Buffer (red cap)	115.2 µL	460.8 µL	
Pre-Amplification Enzyme (white cap)	4.8 µL	19.2 µL	
Amplification Buffer (orange cap)	648 µL	3 × 864 µL	
Amplification Enzyme (blue cap)	12 µL	48 µL	
Nuclease-free Water (clear cap)	432 µL	432 µL	
Barcode Plate (Ion SingleSeq™ Barcode Adapters loaded in wells A1–H3 or wells A1–H12)	1 plate; 5 µL/well (Barcodes 1–24)	1 plate; 5 µL/well (Barcodes 1–96)	

Ion S5™ ExT Kit components

Contents	Amount / box	Storage
Ion S5™ ExT Chef Supplies (Part No. A30673)		
Chip Adapter	1	15°C to 30°C
Enrichment Cartridge v2	1	
Ion S5™ ExT Tip Cartridge	1	
PCR Plate and Frame Seal v2	1 each	
Recovery Station Disposable Lid v2	2	
Recovery Tube	12	
Ion S5™ Chef Solutions (Part No. A27754)		
Ion S5™ Chef Solutions	4 cartridges	15°C to 30°C
Ion S5™ ExT Chef PGS Reagents (Part No. A34676)		
Ion S5™ ExT Chef PGS Reagents	4 cartridges	-30°C to -10°C
Ion S5™ ExT Sequencing Solutions (Part No. A30521)		
Ion S5™ ExT Wash Solution	4 × 1.5 L	15°C to 30°C
Ion S5™ Cleaning Solution	250 mL	
Ion S5™ ExT Sequencing Reagents (Part No. A30522)		
Ion S5™ ExT Sequencing Reagents	4 cartridges	-30°C to -10°C

Ion Chip Kits

One of the following Ion Chip kits is included in each Ion ReproSeq™ PGS Kit:

Description	Included in Ion ReproSeq™ PGS Kit (Cat. No.)	Quantity	Storage
Ion 510™ Chip Kit (Cat. No. A34292)	A34899	4 chips	15°C to 30°C
Ion 520™ Chip Kit (Cat. No. A27761)	A34900	4 chips	
Ion 530™ Chip Kit (Cat. No. A27763)	A34901	4 chips	

Precautions

Handle biological samples and all other materials according to local laboratory guidelines and regulations. Avoid cross-contamination of samples and reagents to ensure quality results.

If equipment or reagents are used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

IMPORTANT! Do not substitute components in this kit with components from other Ion Torrent™ library, template, or sequencing kits. This protocol has been verified using these specific materials. Substitution may adversely affect system performance.

Note: Refer to the expiration date on the package labels. Do not use kits after the expiration date.

Customer and technical support

Visit thermofisher.com/support for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Manufacturer: Life Technologies Corporation | 5781 Van Allen Way | Carlsbad, CA 92008

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0016714

Revision	Date	Description
B.0	21 March 2018	Updated for the Ion GeneStudio™ S5 Systems
A.0	7 June 2017	New Product Information Sheet for the Ion ReproSeq™ PGS Kits – Ion S5™ System

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.