Applied Ambion Ambion

X-Gal

5-Bromo-4-Chloro-3-Indolyl**β-D-Galactopyranoside**

Store at -20°C.

AM9944 Catalog #: Amount: 1 g

Product Description: X-Gal (5-Bromo-4-chloro-3-indolyl-β-D-galactoside) is a chromogenic substrate for β-galactosidase that turns blue

upon hydrolysis. It is commonly used for blue/white screening of recombinants when using vectors that support alpha

complementation of β-galactosidase.

Appearance: White crystalline powder

408 6 Molecular Weight:

Molecular Formula: C₁₄H₁₅BrCINO₆

Chemical Analysis: Purity (HPLC) >99.0%

> Solubility 2% in DMF (N-N' dimethyl formamide)

Structure confirmed using ¹H-NMR

Caution: Avoid contact with eyes, skin, and clothing. Wash thoroughly with water after handling.

Storage Conditions: Store at -20°C, with dessication. Warm to room temperature before opening.

USER INFORMATION

General Information: Ambion® X-Gal is normally used as a 20–50 mg/mL stock solution in N, N'-dimethyl formamide. Prepare a stock

solution of 20-50 mg/mL. Wrap the container in aluminum foil to prevent damage by light and store at -20°C. The

solvent will degrade many plastics, so glass containers are recommended

Applications: X-Gal should be added to agar (at 50°C) prior to pouring plates. Plates or overlays typically contain 40 µg/mL of

X-Gal.

An alternative method for preparing plates containing X-Gal is to spread 50 µL of a 20 mg/mL stock solution onto a

100 mm plate and allow it to dry for 15-30 min prior to use.

After an overnight incubation, blue colonies should be readily visible; refrigerating the plates for 2-5 hr may enhance

the intensity of the blue color.

QUALITY CONTROL

Functional Testing: X-Gal is assayed for efficiency of color development of lac+ colonies and the absence of undesirable inhibitors of cell

growth. A range of concentrations (0 µg/mL to 240 µg/mL; 40 µg/mL recommended) is tested by plating a constant number of E. coli cells harboring pUC19 grown on LB/carbenicillin plates. The plates are evaluated for color intensity,

size and number of colonies after an overnight incubation at 37°C.

OTHER INFORMATION

Material Safety Data Sheets: Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following

address: www.ambion.com/techlib/msds. Alternatively, e-mail your request to MSDS_Inquiry_CCRM@appliedbiosystems.com. Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by

telephone or postal mail. (Requests for postal delivery require 1-2 weeks for processing.)

Warranty and Liability: For research use only. Not for use in diagnostic procedures.

Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service

and technical support teams. Warranty information for the accompanying consumable product is available at www.ambion.com/info/warranty in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about

post-warranty support.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that

may appear in this document.

Applied Biosystems disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Applied Biosystems be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection

with or arising from this document, including but not limited to the use thereof.

| Trademarks, Patents, and Licensing: | Applied Biosystems, AB (Design), and Ambion are registered trademarks of Applera Corporation or its subsidiaries in the US and/or certain other countries. All other trademarks are the sole property of their respective owners. |
|-------------------------------------|---|
| | © 2008 Ambion, Inc. All rights reserved. 4386662A |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |