thermo scientific

6X TriTrack DNA Loading Dye

Catalog Number R1161

Pub. No. MAN0013168 Rev

Rev. B.00

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Contents and storage

Contents	Amount	Storage
6X TriTrack DNA Loading Dye	5 x 1 mL	room temperature or at 4 °C for periods up to 12 months. For longer periods, store at -20 °C

Note

Usage Recommendations

fragments, i.e. ~50 bp.

volumes of DNA sample.

O Mix well, spin down and load.

● Add 1 volume of 6X TriTrack[™] DNA Loading Dye to 5

In 1 % agarose gels bromophenol blue co-migrates

with ~300 bp DNA, while xylene cyanol FF co-migrates

with ~4000 bp DNA. Orange G co-migrates with smallest

Description

The 6X Thermo Scientific TriTrack DNA Loading Dye is used to prepare DNA markers and samples for loading on agarose or polyacrylamide gels. The optimized solution contains three dyes: bromophenol blue, xylene cyanol FF and orange G for visual three-color tracking of DNA migration during electrophoresis. The presence of glycerol in the solution ensures that the DNA in the ladder and sample forms a layer at the bottom of the well; the EDTA included in the solution binds divalent metal ions and inhibits metal dependent nucleases.

Composition

10 mM Tris-HCI (pH 7.6) 0.03 % bromophenol blue 0.03 % xylene cyanol FF 0.15 % orange G 60 % glycerol 60 mM EDTA.

Limited product warranty

Life Technologies Corporation and/or it affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

Thermo Fisher Scientific Baltics UAB V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania
For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For Research Use Only. Not for use in diagnostic procedures.

thermofisher.com/support | thermofisher.com/askaquestion

thermofisher.com



USER GUIDE

07 May 2021