Quant-iT™ dsDNA Assay Kit, High Sensitivity

Catalog Number Q33120

Pub. no. MAN0002497 Rev. A.0

Detailed protocol is available online at www.thermofisher.com.

Protocol summary

- 1. Equilibrate the assay components to room temperature.
- Prepare the working solution by diluting the Quant-iT[™] dsDNA HS reagent 1:200 in Quant-iT[™] dsDNA HS buffer.
- 3. Load 200 μ L of working solution in each microplate well.
- 4. Add 10 μL of each of the Quant-iT™ dsDNA HS standards to separate wells and mix well.
- 5. Add 1–20 μ L of each unknown DNA sample to separate wells and mix well.
- 6. Measure the fluorescence using a microplate reader (excitation/emission maxima ~502/523 nm).
- 7. Use a standard curve to determine the DNA amounts.

For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/support

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1,800,955,6288

