

USER GUIDE

applied
biosystems®
by *life* technologies™

Applied Biosystems® 3730/3730xl DNA Analyzer

GETTING STARTED GUIDE

Publication Part Number 4478016 Rev. A

Revision Date May 2012



life
technologies™

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

The information in this guide is subject to change without notice.

DISCLAIMER

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Limited Use Label License No: 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com.


HITACHI This product includes patented technology licensed from Hitachi, Ltd.

TRADEMARKS

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners. Microsoft and Windows are registered trademarks of the Microsoft Corporation. Oracle is a registered trademark of the Oracle Corporation.



© 2012 Life Technologies Corporation. All rights reserved.

Chapter 1





1

Chapter 2




2

Chapter 3




3

Chapter 4




4

Chapter 5



5

Chapter 6



6

Contents

Chapter 1	Preparing the Instrument	9
	Instrument and Parts	10
	Polymer Delivery Pump Detail	11
	Overview	12
	Troubleshooting Instrument Status Lights	14
	Starting the 3730 Series Data Collection Software	17
	Installing the Capillary Array	18
	Replenishing or Changing Polymer Type	20
	Preparing Buffer and Filling the Reservoirs	22
	Placing Reservoirs into the Instrument	27
Chapter 2	Performing Spatial Calibration	31
	Overview	32
	Performing Spatial Calibration	32
	Evaluating the Calibration Data	34
Chapter 3	Performing Spectral Calibration For Sequencing and Fragment Analysis	39
	Overview	40
	Preparing the Spectral Calibration Chemistry	43
	Creating a Spectral Instrument Protocol	47
	Creating a Spectral Calibration Plate Record	50
	Loading the Plate into the Instrument	53
	Running the Spectral Calibration Plate	54
	Evaluating the Spectral Calibration Data	57
	Examples of Passing Sequencing Spectral Calibrations	61
	Example of a Passing Fragment Analysis Spectral Calibration	63
	Spectral Viewer	64
	Troubleshooting	67

Chapter 4	Setting Up the Software for DNA Sequencing	69
	Plate Records and Sequencing Analysis	70
	Creating Required Settings for Automated Sequencing Analysis	74
	Creating and Completing a Sequencing Analysis Plate Record	93
	Fill Down Special	96
Chapter 5	Setting Up the Software for Fragment Analysis	101
	3730/3730xl Analyzer Data Collection and GeneMapper Software	102
	GeneMapper® Software Plate Records	105
	Components of a GeneMapper® Software Plate Record	106
	Creating Required Settings for Fragment Analysis	109
	Creating and Completing a GeneMapper® Software Plate Record	121
	Filling Down the Plate Record	124
Chapter 6	Running the Instrument	127
	Working with Plate Assemblies	128
	Placing Plate Assemblies into the Instrument	132
	Scheduling Runs	133
	Default Load Maps	137
	Barcode Readers	140
	Running the Instrument: Manual vs Auto Mode	142
	Starting the Run	145
	Controlling the Run	147
	Monitoring the Status of the Run	148
	Viewing Real-Time Electrophoresis Data	150
	Viewing Event History	151
	Viewing Electropherogram Data	152
	Viewing the Run History Data	154
	Viewing the Results of Autoextraction	156
Appendix A	Catalog List	163
Appendix B	Dye Sets: G5, G5-RCT, Any4Dye, Any4Dye-HDR, and Any5Dye	165
	Supported Dye Sets	165
	Dye Sets G5 and G5-RCT For Fragment Analysis	166
	Creating a Spectral Calibration for Dye Sets Any4Dye, Any4Dye-HDR, or Any5Dye	168
	Regular Runs Using Any4Dye or Any5Dye Dye Sets	172


Appendix C	KB™ Basecaller Software v1.4.1	175
	Executive Summary	175
	Benefits of using the KB™ Basecaller	176
	Future support of ABI and KB™ Basecallers	179
	Features in KB™ Basecaller Software v1.4.1	179
	Comparison of the ABI and KB™ Basecallers	180
	Differences between the ABI and KB™ Basecallers	181
	FAQs: Processing data with Phred software and .phd.1 Files	184
	FAQs: Quality values	185
	Miscellaneous FAQs	187
	Conference posters and reference	188
Appendix D	Managing Data Collection Software Licenses	189
	Manage software licenses	189
	Obtain and activate a software license	189
	Renew a software license	191
Appendix E	Safety	195
	Symbols on Instruments	195
	Safety Alerts on Instruments	197
	Instrument Safety	197
	Safety and Electromagnetic Compatibility (EMC) Standards	200
	Documentation and Support	203
	Obtaining SDSs	203
	Obtaining Support	204
	Computer Configuration	204
	Limited Product Warranty	204
	Index	205




Preparing the Instrument

1

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.





Notes _____




Instrument and Parts


Polymer Delivery Pump (PDP)




Notes _____



Polymer Delivery Pump Detail



Notes _____




Overview


This chapter explains how to prepare the instrument for a run by installing the capillary array, buffer, and reservoirs.

Powering On the Computer and 3730/3730x/ Analyzer Instrument

1. Press the power button on the monitor to power it on.
2. Press the power button on the computer to power it on.
3. In the **Log On to Windows** dialog box:
 - a. In the **User Name** field, enter your user name.
 - b. In the **Password** field, enter your password.
 - c. Click .
4. Close the oven door.
5. Close the stacker drawer.




6. Close the instrument door.








7. Wait until the monitor displays the desktop of the Windows® operating system.
8. Press the power button on the 3730/3730x/ Analyzer instrument to power it on.


Notes _____



The Status Lights


Status	Status Light	Action
<ul style="list-style-type: none"> The instrument is ready An automated wizard operation is in progress with the instrument door closed 	Solid Green 	Go to page 17.
<ul style="list-style-type: none"> A run is in progress 	Flashing Green 	—
<ul style="list-style-type: none"> The instrument cannot communicate with the computer 	Solid Yellow 	Go to page 15.
<ul style="list-style-type: none"> The instrument is downloading firmware The instrument is performing diagnostics The oven door is open The instrument door is open The buffer reservoir is not installed The capillary array is not installed An automated wizard operation is in progress with the instrument door open 	Flashing Yellow 	Go to page 14.
<ul style="list-style-type: none"> The instrument has detected a problem 	Solid Red 	Go to page 15.

Notes _____





Troubleshooting Instrument Status Lights

Flashing Yellow




To determine the source of the problem:


1. Press on the instrument door to ensure that it is closed. If the 3730/3730xI Analyzer instrument displays the green status light, then the instrument door was open. Go to [page 17](#).
2. If the 3730/3730xI Analyzer instrument continues to display the flashing yellow light:
 - a. Open the instrument door.
 - b. Press on the oven door to verify that it is closed.
 - c. Close the instrument door.
 - d. If the 3730/3730xI Analyzer instrument displays the green status light, then the oven door was open. Go to [page 17](#).




3. If the 3730/3730xI Analyzer instrument continues to display the flashing yellow light:
 - a. Open the instrument door.
 - b. Open the oven door.
 - c. Check that the buffer reservoir and capillary array are installed.
 - d. Close the oven door.
 - e. Close the instrument door.



Notes _____



Solid Yellow Light




To determine the source of the problem, verify that the:

1. Monitor displays the desktop of the Windows operating system.
2. Ethernet cable is connected to the back of the 3730/3730x/ Analyzer instrument.
3. Other end of the Ethernet cable is connected to the computer.
4. Instrument door is closed.
5. Buffer, water, and waste reservoirs are in place.
6. 3730 Analyzer User account password is functional.

If the instrument continues to display the solid yellow light, contact Life Technologies technical support or your service representative for further assistance.


Solid Red Light



To determine the source of the problem:


1. If the instrument continues to display the solid red light:
 - a. Power off the instrument.
 - b. Wait for 30 seconds.
 - c. Power on the instrument.
2. If the instrument continues to display the solid red light:
 - a. Start the 3730 Series Data Collection Software as explained [page 17](#).
 - b. In the navigation pane of the Data Collection Software, double-click **GA Instruments > ga3730 > instrument name > Instrument Status > Event Log**.

Notes _____




Chapter 1 Preparing the Instrument

Troubleshooting Instrument Status Lights



- c. In the Event Log view, find the last message in the log file.
 - d. Using the error code, perform the required tasks to fix the problem.
3. If the instrument continues to display the solid red light, contact Life Technologies technical support or your service representative for further assistance.

Notes _____




Starting the 3730 Series Data Collection Software

1. Select  > **All Programs** > **Applied Biosystems** > **Unified Data Collection** > **Run Unified Data Collection 4.**


The Data Collection Software opens the Service Console dialog box.

Note: The 3730 Series Data Collection Software 4 requires a license to run. Refer to [Appendix D, Managing Data Collection Software Licenses](#) on page 189 for more details.




Red circles indicate that applications of the Data Collection software are not running.

2. Wait for the Service Console dialog box to open the applications of the Data Collection Software.




When the Data Service component displays the yellow triangle, do not press Start All or Stop All; if you do press either, you will need to re-boot the computer.




3. When all applications are running (green squares), the Data Collection Software opens the Data Collection Viewer.

Note: Ensure that all Data Collection Services are running before you launch the AB Navigator tool for security, audit trail and electronic signature features described in the AB Navigator Software Administrator Guide (Part no. 4477853). All services are running when the Service Console contains four green squares.



Applications of the Data Collection Software are running

Notes



Installing the Capillary Array

WARNING CHEMICAL HAZARD.

POP 7™ polymer may cause eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Use for research and development purposes only.

WARNING CHEMICAL HAZARD.

Running Buffer with EDTA causes eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Required Materials

- Capillary array, 96- or 48-capillary
- Lab wipes, lint-free
- Gloves


Guidelines for Capillary Use

- Do not bend the capillaries
- Store capillary arrays using a buffer reservoir and the header shipping cover. For storage information refer to the *Maintenance and Troubleshooting Guide* (Part no. 4477797).

Installing a New or Used Capillary Array


IMPORTANT! Wear gloves when you handle the capillary array.




 **CAUTION** Failure to use the Install Array wizard when changing capillary arrays can result in degraded analysis data.


1. Close the instrument door.
2. In the Data Collection software, select **GA Instruments > ga3730 > instrument name >**.

Notes






3. On the toolbar, select **Wizards > Install Array Wizard**.
4. Install the array as instructed by the Array wizard.
5. Perform a spatial calibration (see [page 32](#)).




Notes _____



Replenishing or Changing Polymer Type

IMPORTANT! Always replace polymer that has been on the instrument longer than one week.



If polymer on the instrument...	Then ...
Has been on less than one week and is in sufficient quantity to complete your runs	Remove all bubbles, and then proceed with instrument preparation.
Has been on less than one week, and insufficient in quantity to complete your runs	Add fresh polymer to the polymer supply by following the Replenish Polymer Wizard.
Has been on longer than one week	
Is the wrong type (a change between POP-4®, POP-6™, and/or POP-7™ polymers is required)	Replace the installed polymer type with a different type by following the Change Polymer Type Wizard.

Before Using the Polymer

1. Remove the polymer from 4°C storage.
2. Loosen the cap and bring the polymer to room temperature.
3. To dissolve deposits, tighten the cap and gently swirl the polymer.

Replenishing the Polymer


IMPORTANT! Wear gloves while handling polymer, the capillary array, septa, or buffer reservoirs.



CAUTION CHEMICAL HAZARD. POP polymer may cause eye, skin, and respiratory tract irritation. Please read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Use for research and development purposes only.

1. Click <Instrument Name> in the tree pane.

Notes _____




2. Select **Wizards > Replenish Polymer Wizard** to replenish polymer.

IMPORTANT! The polymer type defined in the wizard must match the polymer type that you are using.



Changing Polymer Type

IMPORTANT! Wear gloves while handling polymer, the capillary array, septa, or buffer reservoirs.




 **CAUTION** **CHEMICAL HAZARD. POP polymer** may cause eye, skin, and respiratory tract irritation. Please read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Use for research and development purposes only.

1. Click **<Instrument Name>** in the tree pane.
2. Select **Wizards > Change Polymer Type Wizard** to change to a different polymer.



Notes _____



Preparing Buffer and Filling the Reservoirs



WARNING CHEMICAL HAZARD.

Running Buffer with EDTA causes eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Required Materials


- Retainer, buffer/water/waste
- Septa
- Reservoir caps
- Reservoir, buffer/water/waste
- Plate base, water/waste
- Plate base, buffer
- Water, deionized, 180 mL plus, 160 mL for water and waste reservoirs
- 10× Genetic Analyzer Running Buffer with EDTA, 20 mL
- Graduated cylinder, 250-mL
- Gloves, silicone-free, powder-free

Buffer Storage

The 1× run buffer can be stored at:

- 2–8°C for up to 1 month
- Room temperature for 1 week

Notes _____



When to Change the Buffer


Replace the buffer in the reservoirs every 48 hours, or before each batch of runs.

Note: When replacing all liquids, you should not simply ‘top off’. Replacement is critical.

Note: Clean the reservoirs weekly in warm water followed by a rinse with deionized water.

IMPORTANT! Failure to replace buffer may lead to loss of resolution and data quality.

Notes _____



Preparing the 1× Run Buffer


IMPORTANT! Wear gloves when you handle running buffer with EDTA.



WARNING CHEMICAL HAZARD.

Running Buffer with EDTA causes eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

1. Pour 20 mL 10× running buffer with EDTA into a graduated cylinder.
2. Add 180 mL deionized water to bring the total volume to 200 mL.
3. Mix well and set aside.




Filling the Water and Buffer Reservoirs


IMPORTANT! Wear gloves when you handle the reservoir.




1. Close the instrument door.
2. Press the Tray button to bring the autosampler to the forward position.
3. Wait for the autosampler to stop moving and for the green status light to illuminate before you open the instrument door.



Notes _____



4. Unplug the buffer reservoir. Remove the buffer, water, and waste reservoir assemblies from the instrument.
5. Disassemble each reservoir assembly then empty the contents of the reservoirs into an aqueous waste container.



6. Rinse each reservoir using deionized water.

Note: Be sure to clean the buffer jar, as well as water, waste, and buffer reservoirs weekly in warm water followed by a rinse with deionized water.



7. Dry the reservoirs using lint-free wipes.
8. Fill then assemble the reservoirs.


Buffer Reservoir Assembly

- a. Add 80 mL 1X run buffer to the Buffer reservoir.
- b. Assemble the reservoir assembly as shown below:

Water and Waste Reservoir Assemblies

- a. Add 80 mL high-quality deionized water to each reservoir.
- b. Assemble each reservoir assembly as shown below:

Notes _____





9. To prevent damage to the capillary array, inspect each reservoir assembly and verify that the:

- Septa fit snugly and flush on the reservoir cap

Note: Inspect septa weekly and replace any that are worn or discolored.


- Rubber gasket around the edge of the reservoir cap is seated correctly
- Holes of the plate retainer and the septa strip are aligned



10. Dry the reservoirs using lint-free wipes.



Notes _____




Placing Reservoirs into the Instrument


WARNING **CHEMICAL HAZARD.**

Running Buffer with EDTA causes eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.


1. Connect the Buffer reservoir plate base cable into the heater outlet within the instrument.



2. Move the buffer reservoir to the Buffer position (left) making sure the cable is out of the way of the autosampler.




3. Place the Water and Waste reservoirs into the instrument. The reservoirs must be in the following order from left to right:
 - a. Buffer reservoir
 - b. Water reservoir
 - c. Waste reservoir



4. Close the instrument door.



Notes



5. Press the Tray button to return the autosampler to the array position.



Filling the Anode Buffer Jar



WARNING CHEMICAL HAZARD.

Running Buffer with EDTA causes eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Replace the anode buffer:


- Before each group of scheduled runs, or at least every 24–48 hours
- Every time you fill the polymer block with new polymer
- Every time you change the buffer reservoir

Note: Complete replacement of all liquids is critical; do not simply ‘top off’ liquids. Be sure to clean the buffer jar, as well as water, waste, and buffer reservoirs weekly in warm water followed by a rinse with deionized water.


IMPORTANT! Wear gloves when you handle the anode buffer jar.




1. Remove the anode buffer jar by pulling it down and twisting it slowly.
2. Empty the anode buffer jar into an aqueous waste container.
3. Rinse the anode buffer jar using deionized water.



Notes _____



4. Rinse the anode buffer jar using 1× run buffer:
 - a. Add 5 mL 1× run buffer to the anode buffer jar.
 - b. Tilt the anode buffer jar 90°.




- c. Rotate the jar to rinse the interior with buffer.
 - d. Empty the anode buffer jar into an aqueous waste container.

5. Add 67 mL 1× run buffer to the jar.
6. Put the anode buffer jar on the instrument with the overflow hole facing you.


Note: The meniscus should line up just under the red fill line when installed on the instrument.

7. Verify that the electrode is immersed in the buffer.
8. If the reservoir fills completely as polymer is added, perform [steps 1 through 7](#) of this procedure to discard and replace the running buffer.

IMPORTANT! Replace buffer if excess polymer is expelled into the anode jar.



Notes _____




Chapter 1 Preparing the Instrument


Placing Reservoirs into the Instrument

Notes _____

Performing Spatial Calibration



Notes



Overview

What a Spatial Calibration Tells You

The 3730 Series Data Collection Software uses images collected during spatial calibration to establish a relationship between the signal emitted by each capillary and the position where is detected by the CCD camera.

When to Perform a Spatial Calibration


For all dye sets, perform a spatial calibration after you:

- Install a new or used capillary array
- Remove the capillary array from the detection cell block (even to adjust it)
- Move the instrument (even if the instrument was moved on a table with wheels)
- Move the array detection cell


Note: Failure to perform a new spatial calibration can result in poor data quality.

Performing Spatial Calibration

1. In the navigation pane of the Data Collection Software, double-click **GA Instruments** > **ga3730** > **instrument name** > **Spatial Run Scheduler**.




Notes _____



2. In the Spatial Run Scheduler view, do one of the following:

- If the capillaries contain fresh polymer, select **Protocol > SpatialNoFill**.
- Otherwise, select **Protocol > SpatialFill**.


Note: You do not need to fill the capillaries each time you perform a spatial calibration.



3. Click  .

The approximate calibration run times are:

- 48-cap/36cm array with fill, 4 minutes.
- 96-cap/36cm array with fill, 3 minutes.
- No fill, 2 minutes.



4. Evaluate the calibration as explained on [page 34](#).

Notes _____



Evaluating the Calibration Data

Note: Examples of passing spatial calibration profiles start on [page 37](#).

1. Verify that the peaks of the spatial calibration are approximately the same height.

Are the peaks in the profile approximately the same height?



Yes – Go to [step 2 on page 35](#).

No – How does the peak height vary?

- If the peak height increases at the beginning and the end of the spatial profile, then the variation in peak height is acceptable.

Go to [step 2 on page 35](#).

Irregular – If the peak heights are irregular, go to [“If the Calibration Fails”](#) in the *Maintenance and Troubleshooting Guide* (Part no. 4477797).





Magnifying the Spatial Profile

- a. Click and drag the cursor to create a box around the area of interest.


- b. Release the mouse button.

The Data Collection Software displays the selected region.

- c. Press **R** to reset the view.



Notes



2. Verify that an orange cross appears at the top of each peak in the profile.


Does a cross appear at the top of each peak?

Yes – Go to [step 3](#).

No – Where in the profile is the peak located?

- Left side of the profile:
If using a 96-capillary array, a small peak may appear in the left side of the profile.
The peak is normal, go to [step 3](#).

- After the first peak:
The Data Collection Software did not locate the peak correctly.
Move an orange cross to cover the peak. See, “[To move an orange cross](#)” in the *Maintenance and Troubleshooting Guide* (Part no. 4477797).




3. Check the profile for irregular peaks.

Does the profile contain any irregular peaks?


Yes – The calibration run has failed. Go to “[If the Calibration Fails](#)” in the *Maintenance and Troubleshooting Guide* (Part no. 4477797).

No – Go to [step 4](#).

Elements of a poor spatial



Notes _____



4. Examine each row of the 96 Capillary Position table. Typical values for the **Left spacing** and **Right spacing** columns are:

- 4–8 pixels for a 96-capillary array
- 9–11 pixels for a 48-capillary array

Note: Values greater than those stated above are acceptable if you are able to see a corresponding gap in the capillaries in the detection cell.

Be sure to account for all capillaries (e.g., 96 capillary positions for 96 capillary array).

- If *not*, verify that all peaks have crosses. If each peak does not each have a cross, see the Troubleshooting table below.
- If *yes*, go to step 5.


5. Accept or reject the spatial calibration as follows:

If the calibration:


- Passed, click **Accept** writes the calibration data to the database.
- Failed, click **Reject**, then go to **“If the Calibration Fails”** in the *Maintenance and Troubleshooting Guide* (Part no. 4477797).

Capillary	Position (pixels)	Left spacing	Right spacing
1	11	0	6
2	17	6	5
3	22	5	5
4	27	5	5
5	32	5	5
6	37	5	6
7	43	6	5
8	48	5	5
9	53	5	5
10	58	5	5

Left spacing and Right spacing columns




Notes _____




Examples of Passing Spatial Profiles

IMPORTANT! Improper peak identification may lead to sample mistracking on the instrument, and potential sample misnaming.




Passing Profile #1

This example shows a typical passing profile.



Passing Profile #2




Passing Profile #3

Background artifact

This example shows a passing profile with high artifactual background at the left margin.


Notes




Chapter 2 Performing Spatial Calibration

Evaluating the Calibration Data


Notes _____



Performing Spectral Calibration For Sequencing and Fragment Analysis



Notes _____



Overview

A spectral calibration creates a matrix that is used during a run to reduce raw data from the instrument to the 4- or 5-dye data stored in the sample files. Performing a spectral calibration is similar to performing a sample run, except that calibration standards are run in place of samples, and a spectral calibration module is used in place of a run module.

IMPORTANT! Do not run your computer's Internet Connection wizard during a spectral calibration.

Note: A spectral calibration algorithm checks dye order. If the algorithm determines that the dyes are not in the correct order, the error message is "failed calibration due to bad data: Bad dye order detected". It is possible for the major peaks of the matrix standard to appear in the correct order and still receive this error message.

Spectral calibrations are performed with a specific combination of:

- Dye set (G5, G5-RCT, Any4Dye, Any4Dye-HDR, Any5Dye, E or Z). For further information see, "[Preparing the Spectral Calibration Chemistry](#)" on page 43 and, "[Dye Sets: G5, G5-RCT, Any4Dye, Any4Dye-HDR, and Any5Dye](#)" on page 165.
- Array type (48- or 96-capillary)
- Array length (36- or 50-cm)

IMPORTANT! Spectral calibration must be calibrated for dye set, array type, and array length.

When to Perform the Calibration

Perform a spectral calibration:


- Whenever you use a new dye set on the instrument
- After the laser or CCD camera has been realigned/replaced by a service engineer
- If you see a decrease in spectral separation (pull-up and/or pull-down peaks)
- If you alter any condition (dye set, array type, array length, or polymer type)

Note: Life Technologies recommends that you run a spectral calibration each time a new capillary array is installed. In the 3730 Series Data Collection Software, if you install an array that is the same length as the previously installed array, the active spectral calibration still persists. For optimal data quality, perform a new spectral calibration before you perform regular runs.

Changing Capillary Array Lengths

For each dye set, a single spectral calibration cannot be used for all capillary array lengths.

Notes _____



- For every sequencing dye set, you must create a separate spectral calibration for each capillary array length and array type.
- For every fragment analysis dye set, you must create a separate spectral calibration for each capillary array length and array type.

Refer to [page 65](#) for information on how to switch calibrations.

Required Materials

Catalog numbers are located in [Appendix A on page 155](#).

Description


- BigDye® Terminator v3.1 or v1.1 Sequencing Standard or, DS-33 Matrix Standard
- 384- or 96-Well Reaction Plate w/ Barcode
- Multichannel pipettor
- Plate retainer
 - Plate septum with black plate base
- or*
- Heat-seal with gray plate base
- Hi-Di™ Formamide
- Heated block or thermal cycler
- Container with ice
- Centrifuge with microplate adapter
- Microcentrifuge
- Vortex
- Gloves

Two Types of Calibration Standards

Two types of calibration standards are used to create a matrix:

- For Fragment Analysis – Matrix standards are four or five fragments of varying size that are individually labeled with one of the four or five dyes of a set.
- For Sequencing – Sequencing Standards are standard sequencing reaction fragments of varying size that are individually labeled with one of the four dyes.

Notes _____




Select Dye Sets and Calibration Standards

Use the following tables to determine the correct dye set and calibration standard for the application you are using.

Sequencing Chemistry	Dye Set	Calibration Standards
BigDye® Terminator v3.1 Cycle Sequencing Kit	Z_BigDyeV3	BigDye® v3.1 Terminator Sequencing Standard
BigDye® Direct Cycle Sequencing Kit	Z_BigDyeV3	BigDye® v3.1 Terminator Sequencing Standard
BigDye® Terminator v1.1 Cycle Sequencing Kit	E_BigDyeV1	BigDye® v1.1 Terminator Sequencing Standard

Fragment Analysis Chemistry	Dye Set	Calibration Standards
Fragment Analysis	G5	DS-33
Fragment Analysis	G5-RCT	DS-33
SNaPshot® Multiplex Kit	Any5Dye	DS-02

Notes _____



Preparing the Spectral Calibration Chemistry

WARNING CHEMICAL HAZARD.

Formamide causes eye, skin, and respiratory tract irritation. It is a possible reproductive and birth defect hazard. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.


1. Dilute the spectral calibration standard with Hi-Di™ Formamide according to the insert instructions.

2. Vortex thoroughly.

3. Briefly centrifuge the mixture.

4. Heat the standard tube at 95°C for 5 minutes to denature the DNA.

5. Cool the tubes on ice for 2 minutes.



Notes _____



6. Vortex thoroughly and then briefly centrifuge the mixture.

Sealing and Preparing the Plate Assemblies



WARNING Do not use warped or damaged plates.



1. Add the denatured standard to the wells of a 384- or 96-well reaction plate:

If using a:

- **48-capillary, 96-well plate** – Add 10 μL of denatured standard to each well.
- **384-well plate** – Add 5 μL of denatured standard into alternating wells of the plate. See [page 137](#) for load maps.




2. Seal the plate with a septum or heat-seal:

With a septum:


- a. Inspect the septa and be sure to replace any that are worn or discolored.
- b. Place the plate on a clean, level surface.
- c. Lay the septum flat on the plate.
- d. Align the holes in the septum strip with the wells of the plate, then firmly press downward onto the plate. Ensure that:
 - The septa lie flat against the plate. You should not feel any lumps or raised edges.
 - The septa are inserted straight into the wells. You should not see any bent or crooked duckbills when viewing the plate from above.

With heat-seal:


- a. Follow your thermal sealer instrument instructions.



Notes _____




3. Briefly centrifuge the plate.




4. Remove the plate from the centrifuge and verify that each sample is positioned correctly in the bottom of its well.


If the reagents of any well contain bubbles or are not located at the bottom of the well, repeat steps 3 and 4.



5. Assemble the plate assembly as shown in the following figures (see [Appendix A, “Catalog List,”](#) on page 163 for catalog numbers).




WARNING Use only **black** plate bases with septa-sealed plates. If you are using MicroAmp™ Fast 96-Well Reaction Plates (0.1 ml), use only **blue** plate bases and matching retainer.



WARNING Use only **gray** plate bases with heat-sealed plates. If you are using MicroAmp™ Fast 96-Well Reaction Plates (0.1 ml), use only **dark green** plate base and matching retainer.

Notes



6. Verify that the holes of the plate retainer and the septa are aligned.


IMPORTANT! The plate may damage the array if the retainer and the septum holes are not aligned.

7. Make sure when you assemble a plate that the retainer clip is flush with the plate base. A simple way to ensure that they are flush is to run your finger along the edge.

Important Heat Seal Recommendations


- Use 3-mil Life Technologies heat seal film (Cat. no. 4337570). This film is 3-mil before, and 1-mil after, heating.
- *Do not* use heat seal film thicker than 1-mil, after heating, on the 3730/3730xl DNA Analyzer.
- *Do not* use heat-seal film containing adhesives or metals as these may damage the instrument's piercing needles.

Notes _____




Creating a Spectral Instrument Protocol


1. In the navigation pane of the Data Collection Software, click **GA Instruments** > **ga3730** > **Protocol Manager**.




2. In the Instrument Protocols pane, click **New...**. The Protocol Editor opens.



3. Select **Spectral** from the Run Module dropdown list.



Notes



4. The Protocol Editor now displays additional drop-down lists. Select from the following:

If you are using a *matrix standard* for spectral calibration, you can use a 36-cm or 50-cm array length:

- For a 36-cm capillary array, use:
 - Run Module: **Spect36_MtxStd_1**
 - Chemistry: **matrixStandard**

or

- For a 50-cm capillary array, use:
 - Run module: **Spect50_MtxStd_POP-7™_1**
 - Chemistry: **matrixStandard**

IMPORTANT! The array length you select must match the array length information from the Install Array wizard.

If you are using a *sequencing standard* for spectral calibration, you can use a 36-cm or 50-cm array length:



- For a 36-cm capillary array, use:
 - Run module: **Spect36_SeqStd_1**
 - Chemistry: **sequenceStandard**


or

- For a 50-cm capillary array, use:
 - Run module: **Spect50_SeqStd**
 - Chemistry: **sequenceStandard**

Note: The Chemistry file for fragment analysis dye sets automatically defaults to the matrix standard.

IMPORTANT! The array length you select must match the array length information from the Install Array wizard.






Use the following table to select the correct chemistry file for the spectral calibration samples you use.

Dye Sets, Standards, And Chemistry Files

Dye Set	Standard Type	Chemistry File
Z_BigDyeV3	BigDye® v3.1 Terminator Sequencing Standard	Sequence Standard
E_BigDyeV1	BigDye® v1.1 Terminator Sequencing Standard	Sequence Standard

Dye Set	Matrix Standard Set	Chemistry File
G5	DS-33	Matrix Standard
G5-RCT	DS-33	Matrix Standard

1. (Optional) Click **Edit Param** to display the Spectral Params dialog box.
2. Use this dialog box to edit the selection criteria for passing or failing spectral calibrations.




Valid Data Ranges

Parameters	Valid Data Ranges†
Matrix Condition Number Bounds	Lower: 1–10 Upper: 3–20
Locate Start Point	After Scan: 100–5000 Before Scan: 100–5000
Limit Analysis (scans)	400–20,000
Sensitivity	0–0.9
Minimum Quality Score	0.80–0.99

† These ranges are dye-set independent

IMPORTANT! Default parameter values are optimized and are recommended for most situations


Notes




Creating a Spectral Calibration Plate Record

1. In the navigation pane of the Data Collection Software, double-click

GA Instruments > **ga3730** >
instrument name > **Plate Manager**.




2. Click **New** to create a new plate.




3. Complete the New Plate dialog box:

- a. Enter ID or Barcode number
- b. Enter a name for the plate.
- c. (Optional) Enter a description for the plate record.
- d. In the Application drop-down list, select **Spectral Calibration**.
- e. In the Plate Type drop-down list, select **96-Well** or **384-Well**.
- f. Enter desired scheduling. For more information see, [“Globally Modifying a Run Schedule”](#) on page 135.




Notes _____



- g. In the Plate Sealing drop-down list, select **Septa** or **Heat Seal**.
 - h. Enter a name for the owner.
 - i. Enter a name for the operator.
 - j. Click .
4. In the Spectral Calibration Plate Editor, enter the following information:

Note: This example assumes that you are loading the first quadrant.


- a. In the Sample Name column of row A01, enter a sample name, then click the next cell.
- b. In the Comments column of row A01, enter any additional comments or notations for the sample at the corresponding position of the plate.
- c. In the Instrument Protocol 1 column of row A01, select a protocol from the drop-down list.




5. Select the entire row.
6. Select **Edit > Fill Down Special**.

Based on the plate type (96- or 384-well) and capillary array (48 or 96 capillaries) you are using, select the appropriate fill down option:

- 96 capillary/96-well plate: **Fill Down**
- 48 capillary/96-well plate: **Fill down Special (48 Cap)**
- 96 capillary/384-well plate: **Fill down Special (96 Cap)**
- 48 capillary/384-well plate: **Fill down Special (48 Cap)**



Notes



7. Click .


You have successfully created a plate record for the spectral calibration plate.

Note: If multiple cells are selected for copying, select the same number of corresponding target cells before you execute the Paste command.

Note: The Plate Editor Copy and Paste functionality is supported only within one plate editor. To copy and paste the contents of one plate to another plate, use the “Duplicate...” button on the Plate Manager dialog box.


Note: If you use the duplicate plate function, all the information in the plate to be duplicated must be valid. Otherwise, an empty plate is created.

Notes _____



Loading the Plate into the Instrument


1. The name of the plate record you just created is displayed in the Input Stack window of the Data Collection Software, and is ready to run.
2. Open the stacker drawer.
3. Open the In Stack tower door.



Stacker drawer

4. Place the plate assembly into the stacker.

IMPORTANT! The plate must be oriented so that the notched corner of the plate assembly is at the rear-right corner of the stacker.




Notched corner of the plate assembly

5. Close the In Stack tower door.
6. Close the Stacker drawer.




In Stacker tower door

Notes _____




Running the Spectral Calibration Plate

1. In the navigation pane of the Data Collection Software, double-click
GA Instruments > **ga3730** >
instrument name > **Run Scheduler**.
2. In the Run Scheduler view:
 - In the Add Plate field, scan the bar code of a plate to add it to the input stack.
or
 - Type the plate ID then press **Enter** to add it to the input stack.



3. In the toolbar of the Data Collection Software window, click to begin the run.
4. The Processing Plates dialog box opens.
5. Click .




Note: The instrument may pause before running the plate to raise the oven temperature.

Application	Capillary Array Length (cm)	Approximate Spectral Run Time† (min)
Sequencing	50	120
Sequencing	36	60
Fragment Analysis	36	32

† The Data Collection Software may take up to 30 min to calculate the matrices after the run.

6. When the run is finished, remove the plate from the instrument.

Notes _____



Viewing the Pass/Fail Status After the Run


After the instrument completes the spectral calibration run, the pass or fail status of each capillary is recorded in the Events Messages section of the Instrument Status window.

1. In the navigation pane of the Data Collection Software, select
GA Instruments > ga3730 > instrument name > Instrument Status > Event Log.

Type	Date	Time	Publisher	Description
Info	09/05/03	16:41:03	3730C5	Capillary 32 successfully calibrated : q=0.900 c=5.66
Info	09/05/03	16:41:05	3730C5	Capillary 31 successfully calibrated : q=0.957 c=5.72
Info	09/05/03	16:41:04	3730C5	Capillary 30 failed calibration : Failed quality check: q=0.94484 is less than minQ threshold (0.95000)
Info	09/05/03	16:41:04	3730C5	Capillary 29 successfully calibrated : q=0.965 c=5.55
Info	09/05/03	16:41:04	3730C5	Capillary 28 successfully calibrated : q=0.958 c=5.59
Info	09/05/03	16:41:03	3730C5	Capillary 27 failed calibration : Failed quality check: q=0.93434 is less than minQ threshold (0.95000)
Info	09/05/03	16:41:03	3730C5	Capillary 26 successfully calibrated : q=0.970 c=5.62
Info	09/05/03	16:41:02	3730C5	Capillary 25 successfully calibrated : q=0.964 c=5.57
Info	09/05/03	16:41:02	3730C5	Capillary 24 successfully calibrated : q=0.967 c=5.57
Info	09/05/03	16:41:02	3730C5	Capillary 23 successfully calibrated : q=0.966 c=5.62
Info	09/05/03	16:41:01	3730C5	Capillary 22 successfully calibrated : q=0.976 c=5.67
Info	09/05/03	16:41:01	3730C5	Capillary 21 successfully calibrated : q=0.957 c=5.70
Info	09/05/03	16:41:00	3730C5	Capillary 20 successfully calibrated : q=0.957 c=5.60

3

Notes



2. In the Events Messages section of the window, view the status of each capillary.

Type	Date	Time	Publisher	Description
Info	09/05/03	16:41:03	3730C5	Capillary 32 successfully calibrated : q=0.960 c=5.00
Info	09/05/03	16:41:05	3730C5	Capillary 31 successfully calibrated : q=0.957 c=5.72
Info	09/05/03	16:41:04	3730C5	Capillary 30 failed calibration : Failed quality check: q=0.94484 is less than minQ threshold (0.95000)
Info	09/05/03	16:41:04	3730C5	Capillary 29 successfully calibrated : q=0.965 c=5.55
Info	09/05/03	16:41:04	3730C5	Capillary 28 successfully calibrated : q=0.958 c=5.59
Info	09/05/03	16:41:03	3730C5	Capillary 27 failed calibration : Failed quality check: q=0.93434 is less than minQ threshold (0.95000)
Info	09/05/03	16:41:03	3730C5	Capillary 26 successfully calibrated : q=0.970 c=5.67


Dye set G5 status results

For a good-quality calibration, each capillary should have a:

- Q-value:
 - > 0.95 for matrix standards
 - > above 0.93 for sequence standards
- Condition number range, indicated below, for each dye set:

Dye Set	Default Condition Number Range
Sequencing Analysis	
Z_BigDyeV3	2.5–4.5
E_BigDyeV1	3.0–5
Fragment Analysis	
G5	9.5–14.5
G5-RCT	9.5–14.5

Notes _____



Evaluating the Spectral Calibration Data


IMPORTANT! Review and evaluate the spectral calibration profile for each capillary, even if the Spectral Calibration Results box indicated that they all passed.

Note: Pages 61 and 62 contain examples of passing sequencing spectral calibration profiles, and page 63 contains an example of a passing fragment analysis spectral calibration profile.


1. In the navigation pane of the Data Collection Software, select
GA Instruments > ga3730 >
instrument name > Spectral Viewer.

The screenshot shows the 'Spectral Viewer' window in 'Foundation Data Collection Version 3.0'. The navigation pane on the left shows a tree structure with 'GA Instruments' expanded to 'ga3730' and 'C5', with 'Spectral Viewer' selected. The main window contains three panels. The top panel, 'Spectral profile', shows a line graph of intensity vs pixel number (0-400) with four colored curves. The middle panel, 'Raw data (matrix standards)', shows a bar chart of intensity vs scan number (0-8000) with multiple colored bars. The bottom panel, 'Capillary Data', includes a 'Plate diagram' (a 96-well grid with 'A01' and 'H1' labeled), 'Dye Set' (Z-BigDyeV3), 'Matrix used for Capillary 82: 82', 'Condition: 3.877218', 'Q Value: 0.966001', and 'Active Calibration for Dye Set: Z-BigDyeV3' (Thu Jun 19 19:43:46 PDT 2003). A callout box points to the 'Active Calibration' field with the text 'Rename or set the active spectral calibration here'. The status bar at the bottom shows 'System Status: Ready' and 'No Current Run'.


Notes



2. In the Dye Set drop-down list, select the dye set you just created.



3. Select a well on the plate diagram to view the spectral results of the associated capillary.




4. Evaluate the spectral calibration profile for the selected capillary:


- a. Verify that the order of the peaks in the spectral profile from left to right are:
 - 4-dye–blue-green-yellow-red
 - 5-dye–blue-green-yellow-red-orange

If the peaks in the profile:

- Are in the correct order–go to [step a.](#)
- The calibration run has failed–go to [page 67.](#)




Example of a 4-dye spectral calibration profile



Example of a 5-dye spectral calibration profile

Notes _____



- a. Verify that the peaks in the spectral profile do not contain gross overlaps, dips, or other irregularities (see “Tip: Magnifying the Spectral Profile” on page 60).

If the peaks in the spectral profile are:

- Separate and distinct—the capillary has passed. Go to [step 5](#).
- Not separate and distinct—the calibration run has failed. Go to [page 67](#).

- a. Verify that the order of the peaks in the raw data profile from left to right are:



Fragment Analysis

- 5-dye: orange-red-yellow-green-blue

Are the peaks in the wrong order or are there any extraneous peaks that adversely affect the spectral profile?


Yes: The calibration run has failed. Go to [page 67](#).

No: Go to [step 5](#).




Example of a 4-dye sequencing raw data profile

Left to right: Orange, Red, Yellow, Green, Blue




Example of a 5-dye fragment analysis raw data profile

5. Repeat steps 3 and 4 for each capillary in the array.




Notes _____




6. Rename the spectral run. The spectral file default name is the day, date and time of the run.

- a. Click **Rename**.
- b. (*Optional*) In the Rename Calibration dialog box, enter a descriptive name for the spectral calibration including the dye set, array length and polymer type.
- c. Click **OK**.




Tip: Magnifying the Spectral Profile

1. In the navigation pane of the Data Collection Software, click **GA Instruments > ga3730 > instrument name > Spectral Viewer**.
2. In the profile or raw data display, click-drag the cursor to create a box around the area of interest.
3. Release the mouse button.
The Data Collection Software displays the selected region.
4. Press **R** to reset the view.




Selecting an area to magnify in a spectral profile



Magnified area of that spectral profile

Notes _____



Examples of Passing Sequencing Spectral Calibrations

Dye Set Z Created from a Sequencing Standard

Foundation Data Collection Version 3.0

File View Service Tools Wizards Help

GA Instruments > ga3730 > 3730-09 > Spectral Viewer

Intensity vs Pixel Number

Intensity vs Scan Number

Capillary Data: Thu Jun 19 19:35:50 PDT 2003

A01

H12

Dye Set: Z-BigDyeV3

Matrix used for Capillary 4: 4

Condition: 3.40364

Q Value: 0.968057

Active Calibration for Dye Set: Z-BigDyeV3
Thu Jun 19 19:35:50 PDT 2003

List of Calibrations for Dye Set: Z-BigDyeV3
Thu Jun 19 19:35:50 PDT 2003

Override Spectral

Save

Rename Set

System Status System Status: Ready

3

Notes



Dye Set E Created from a Sequencing Standard

Foundation Data Collection Version 3.0

GA Instruments > ga3730 > 3730-09 > Spectral Viewer

Intensity vs Pixel Number

Intensity vs Scan Number

Capillary Data: Mon Jul 28 18:09:54 PDT 2003

A01

H12

Dye Set: E-BigDyeV1

Matrix: used for Capillary 15: 15

Condition: 3.357109

Q Value: 0.981652


Active Calibration for Dye Set: E-BigDyeV1

List of Calibrations for Dye Set: E-BigDyeV1

Mon Jul 28 18:09:54 PDT 2003

System Status Starting Electrophoresis

Notes




Example of a Passing Fragment Analysis Spectral Calibration

Dye Set G5 Created from Matrix Standard Set DS-33

3

Notes



Spectral Viewer

Selecting Active Spectral Calibrations

For best quality data, Life Technologies suggests that you perform spectral calibrations every time a new array is installed in the instrument. However, you may choose to reuse previous spectral calibrations to apply to new data that will be generated on the instrument. Once data is collected, you cannot reapply a different spectral calibration.

IMPORTANT! It is essential that you perform a spectral calibration any time the capillary array is moved or replaced when using DyeSetG5-RCT.

IMPORTANT! When you install an array that is a different length or type (48 vs. 96) from what you were using previously, all spectral calibrations are voided. If a previous spectral calibration for the new array/new condition does not exist, you must run a new spectral calibration. If a previous calibration exists, go to the Spectral Viewer and choose a past calibration to set as the active spectral calibration before you proceed with regular runs, even though spectral profiles are displayed; to do so, follow the directions described next, in [“To select a previous spectral calibration:”](#) on page 65.

IMPORTANT! You cannot link or run a plate unless the dye set used in the plate has been set in the Spectral Viewer. Furthermore, when a plate is running, the Set Active Spectral Calibration function is inactive. Spectral Calibrations can be set only during the idle or ready mode.

Poor quality data or failed analyses are results of using the wrong spectral calibration.

IMPORTANT! Spectral calibrations must be calibrated for dye set, array type, array length, and polymer type.


When a new *spatial* calibration is saved, the current spectral calibration for DyeSet G5-RCT is deactivated. Dye sets G5, E, and Z are not deactivated. If you wish to continue without a spectral recalibration, you can set an active spectral using the following instructions.

All calibrations for your current dye set are listed in the List of Calibrations drop-down list. Therefore, you can choose a spectral calibration to use from the list before you begin a new run.

Note: An asterisk * precedes failing calibrations.

Note: The most recent spectral for each dye set is automatically chosen as the active calibration.


Notes _____




Because each dye set can have its own active calibration, there is no need to manually set the active calibration if you are performing runs with various dye sets.

To select a previous spectral calibration:


1. Select the dye set of interest.
2. In the Spectral Viewer, click the List of Calibrations drop-menu in the lower right-pane.




3. Select the spectral calibration you want to use for future runs.




Notes _____




4. Click **Set** to display your chosen spectral calibration in the Active Calibration text box.





5. (Optional) Click **Rename** to display the Rename Calibration dialog box, enter a new name, then click **OK**.




Notes _____



Troubleshooting

Troubleshooting spectral calibration		
Observation	Possible Cause	Recommended Action
No signal.	Incorrect sample preparation.	Replace samples with fresh samples prepared with fresh Hi-Di™ Formamide.  WARNING CHEMICAL HAZARD. Formamide causes eye, skin, and respiratory tract irritation. It is a possible reproductive and birth defect hazard. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.
	Air bubbles in sample tray.	Centrifuge samples to remove air bubbles.
If the spectral calibration fails, or if a message displays “No candidate spectral files found”.	Clogged capillary.	Refill the capillaries using manual control. Look for clogged capillaries during capillary fill on the cathode side.
	Insufficient filling of array.	Check for broken capillaries and refill the capillary array.
	Expired spectral standards.	Check the expiration date and storage conditions of the spectral standards. If necessary, replace with a fresh lot.
Spikes in the data.	Expired polymer.	Replace the polymer with a fresh lot using the Change Polymer wizard.  WARNING CHEMICAL HAZARD. POP-7™ polymer cause eye, skin, and respiratory tract irritation. Read the SDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.
	Air bubbles, especially in the polymer.	<ul style="list-style-type: none"> • Refill the capillaries using the Bubble Remove wizard. • Properly bring the polymer to room temperature. • Replace expired polymer.
	Possible contaminant in the polymer.	Replace the polymer using the Change Polymer wizard.

Notes _____





Chapter 3 Performing Spectral Calibration For Sequencing and Fragment Analysis

Troubleshooting

Notes _____

Setting Up the Software for DNA Sequencing



Notes



Plate Records and Sequencing Analysis

Overview A plate record is similar to a sample sheet or an injection list that you may have used with other Applied Biosystems® instruments. Plate records are data tables in the instrument database that store information about the plates and the samples they contain. A plate record contains the following information:

- Plate name, type, and owner
- Position of the sample on the plate (well number)
- Sample
- Name, see page [page 87](#)
- Mobility file (in Analysis Protocol), see page [page 80](#)
- Comments about the plate and about individual samples
- Name of the run module and Dye set information (run modules specify information about how samples are run) (in Instrument Protocol), see [page 74](#)
- Name of the Analysis Protocol (Analysis protocols specify how data is analyzed at the end of the run; see page [page 80](#))

Important Notes

- A unique name must be assigned to the instrument computer before 3730 Series Data Collection Software is installed.
- Do not rename the computer once 3730 Series Data Collection Software has been installed. Doing so *will* cause the 3730 Series Data Collection Software to malfunction.

File-Naming Convention Alphanumeric characters that are not valid for user names or file names are:
spaces

\\ : * ? " < > |


An error message is displayed if you use any of these characters. You must remove the invalid character to continue.

When to Create a Plate Record A plate record must be created for each plate of samples for the following types of runs:

- Spectral calibrations
- Sequencing analysis
- SeqScape analysis (Autoanalysis by SeqScape® is no longer supported)

IMPORTANT! A plate record must be created in advance of the first run. Plate records can be created, and plates added to the stacker, while a run is in progress.

Notes _____




Sequencing Analysis Plate Record

The Plate Editor opens an empty plate record for the application that you select in the New Plate dialog box. The data fields within a given plate record vary, depending on the selected application. This section describes the data fields that are present in a sequencing analysis plate record.

The following table and flow chart describe what each file specifies.

Parameters	Description	See Page
Instrument Protocol	Contains everything needed to run the instrument.	74
Analysis Protocol	Contains everything needed to analyze sequencing data.	79
Results Group	Defines the file type, the file name, file save locations, analysis software and autoanalysis by DNA Sequencing Analysis Software 6.	85




Elements of a sequencing analysis plate record

IMPORTANT! For data collection and autoanalysis to be successful, each run of samples must have an instrument protocol, an analysis protocol, and a results group assigned within a plate record. Autoanalysis by SeqScape® is no longer supported; use autoanalysis with DNA Sequencing Analysis Software 6.

Notes _____



Chapter 4 Setting Up the Software for DNA Sequencing
Plate Records and Sequencing Analysis



Default is one sample run. To add additional runs, see [page 95](#).

Blank sequencing analysis plate record

The following table describes the columns inserted in a Plate Record for a sequencing analysis run.

Name	Description
(1.) Sample Name	Name of the sample
(2.) Comment	<i>(Optional)</i> Comments about the sample
(3.) Results Group	Options are: <ul style="list-style-type: none"> • New – Opens the Results Group Editor dialog box • Edit – Opens the Results Group Editor dialog box for the results group listed in the cell • None – Sets the cell to have no selected results group • Select one of the available results groups from the list Note: You must have a results group selected for each sample entered in the Sample Name column. See, “ Results Groups ” on page 85 .

Notes _____



Name	Description
(4.) Instrument Protocol	<ul style="list-style-type: none"> • New—Opens the Protocol Editor dialog box. • Edit—Opens the Protocol Editor dialog box for the instrument protocol listed in the cell. • None—Sets the cell to have no selected protocol. • List of instrument protocols—In alphanumeric order. <p>Note: You must have an Instrument Protocol selected for each sample entered in the Sample Name column.</p> <p>See, “Creating an Instrument Protocol” on page 74.</p>
(5.) Analysis Protocol	<ul style="list-style-type: none"> • New—Opens the Analysis Protocol Editor dialog box. • Edit—Opens the Analysis Protocol Editor dialog box for the instrument protocol listed in the cell. • None—Sets the cell to have no selected protocol. • List of Analysis Protocols—In alphanumeric order <p>Note: You must have an Analysis Protocol selected for each sample entered in the Sample Name column.</p> <p>See, “Creating an Analysis Protocol” on page 80.</p>

Notes _____



Creating Required Settings for Automated Sequencing Analysis

If Settings Already Exist


If the appropriate instrument protocol, analysis protocol, and results group have been created, proceed to “[Creating and Completing a Sequencing Analysis Plate Record](#)” on page 93.

Instrument Protocols

An instrument protocol contains all the settings necessary to run the instrument. An instrument protocol contains the protocol name, type of run, run module, and dye set.

Creating an Instrument Protocol


1. In the navigation pane of the Data Collection Software, select **GA Instruments** > **ga3730** > **Protocol Manager**.




Create instrument protocols here

Create analysis protocols here


Notes



2. In the Instruments Protocols section, click **New...**. The Protocol Editor opens.



3. Complete the Protocol Editor:
- Type a name for the protocol.
 - (Optional) Type a description for the protocol.
 - Select **Regular** in the Type drop-down list.
 - Using the information in the table below, select the correct run module for your run.



Note: To customize a run module, see “[Tip: Customizing Run Modules](#)” on page 77.

Sequencing Run Modules	Capillary Array Length (cm)	Sequencing Run	Approximate Run Times [†] (min)	KB Basecaller QV20 LOR (Bases) [§]
XLRSeq50_POP-7™	50	Extra long read	180	900
LongSeq50_POP-7™	50	Long read	120	850
FastSeq50_POP-7™	50	Fast read	60	700
StdSeq36_POP-7™	36	Standard read	60	700
RapidSeq36_POP-7™	36	Rapid read	35	550
TargetSeq36_POP-7™	36	Short read	20 [‡]	400 [‡]
LongSeq50_POP-6™	50	Long read	150	600
StdSeq36_POP-6™	36	Standard read	60	500

[†] These approximate run times assume oven temperature has reached run temperature

[‡] Time stated for 400 bases. Module can be customized to run 200-400 bases.

[§] Length of read with 98.5% basecalling accuracy, and less than 2% N's, using pGEM-32f (+) as template.

Notes



Chapter 4 Setting Up the Software for DNA Sequencing

Creating Required Settings for Automated Sequencing Analysis

Note: If the BigDye Xterminator[®] Purification Kit was used for sequencing reaction clean up, choose the run modules modified for BDx, as marked by 'BDx'. For additional information, refer to Appendix A in the *BigDye Xterminator[®] Purification Kit Protocol* (Part no. 4374408) for the appropriate run modules.

- e. Using the information in the following table, select the correct Dye Set for your run.

Dye Set	Chemistry
E_BigDyeV1	BigDye [®] Terminator v1.1 Cycle Sequencing Kit
Z_BigDyeV3	BigDye [®] Terminator v3.1 Cycle Sequencing Kit
Z_BigDyeV3	BigDye [®] Direct Cycle Sequencing Kit

- f. Click .

Notes _____




Tip: Customizing Run Modules

You can modify default run modules to suit your particular needs.

1. Click **GA Instruments** > **ga3730** > **instrument name** > **Module Manager**.
2. Click **New...**.
The Run Module Editor dialog box opens.
3. Complete the Run Module Editor dialog box:
 - a. Enter a name for your new module.
 - b. In the Type drop-down list, select the type of module (Regular, Spatial, or Spectral).
 - c. In the Template drop-down list, select a template module as a basis for the new module.
 - d. (Optional) Enter a description of your new run module.
 - e. Change to the desired module parameters using the range for the allowable parameters.
 - f. Click **OK**.

Note: You cannot edit a default module installed with 3730/3730xI Analyzer Data Collection software.



Notes



Editable Run Module Parameters

Parameter Name	Range	Comment
Oven_Temperature	18°C–70°C	Temperature setting for main oven throughout run.
PreRun_Voltage	0–15 kV	Pre run voltage setting before sample injection.
PreRun Time	1–1800 sec	Prerun voltage time.
Injection_Voltage	0–15 kV	Injection voltage setting for sample injection.
Injection_Time	1–90 sec	Sample injection time.
First_ReadOut_time	100–16000 millisec	The interval of time for a data point to be produced. First_ReadOut_time should be equal to Second_ReadOut_time.
Second_ReadOut_Time	100–16000 millisec	The interval of time for a data point to be produced. Second_ReadOut_time should be equal to First_ReadOut_time.
Run_Voltage	0–15 kV	Final run voltage.
Voltage_Number_Of_Steps	0–100 steps	Number of voltage ramp steps to reach Run_Voltage. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Voltage_Step_Interval	0–180 sec	Dwell time at each voltage ramp step. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Voltage_Tolerance	0.1–6 kV	Maximum allowed voltage variation. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel. If it goes beyond tolerance and shuts off, contact Life Technologies tech support.
Current_Stability	0–2000 μ A	Maximum allowed electrophoresis current variation. Current fluctuations above this value will be attributed to air bubbles in system and the voltage automatically powered off. We recommend that you do not change this value unless advised otherwise by Life Technologies support personnel.
Ramp_Delay	1–1800 sec	Delay During Voltage Ramp. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Data_Delay	1–1800 sec	Time from the start of separation to the start of sample data collection.
Run_Time	300–14000 sec	Duration data is collected after Ramp_Delay.

Notes _____



Analysis Protocols

An analysis protocol contains all the settings necessary for analysis and post processing:

- Protocol name – The name, description of the analysis protocol, and the sequence file formats to be used.
- Basecalling settings – The basecaller, DyeSet file, and analysis stop point to be used.
- Mixed Bases – (*Optional*): To use mixed base identification, and if so, define the percent value of the second highest to the highest peak.
- Clear Range – The clear range to be used based on base positions, sample quality values, and/or number of ambiguities (Ns) present.

Note: If you create an appropriate analysis protocol in the Sequencing Analysis software, you can use it in Data Collection Software.

IMPORTANT! Do not delete an analysis protocol during a run while it is being used for that run. Autoanalysis by DNA Sequencing Analysis Software 6 will not be performed if you do so.

Notes _____




Creating an Analysis Protocol


Refer to [Appendix C, “KB™ Basecaller Software v1.4.1,”](#) on [page 175](#) and the *DNA Sequencing Analysis Software 6* (Part no. 4474239) for more information regarding analysis protocols

1. In the Analysis Protocol section of the Protocol Manager, click **New...**

If more than one analysis application is installed on the data collection computer, the Analysis Applications dialog box opens.



2. Select **Sequencing Analysis**, then click **OK**.
The Analysis Protocol Editor opens.




Notes _____



3. Select the **General** tab, then:

- a. Enter a unique name and description for the new protocol.
- b. Select the appropriate Sequence File formats settings.

Option	If checked, the software creates...
Write .Seq File check box	a .seq file for printing the sequence as text file or for using the file in other software. <ul style="list-style-type: none"> • ABI format is used with Applied Biosystems® software. • FASTA format is used with other software
Write Standard Chromatogram Format file (.scf)	When selected, the software creates a .scf file that can be used with other software. When created, the .scf extension is not appended to the file name.
Write Phred (.phd.1) File	When selected and the KB basecaller is used, the software creates a .phd.1 file that can be used with other software.




4. Select the **Basecalling** tab, then:

- a. Select the appropriate basecaller and DyeSet/Primer based on the chemistry and capillary array length you are using.

Note: See [Appendix C, “KB™ Basecaller Software v1.4.1,”](#) on page 175 for a comparison of Basecaller options.

Note: Select Sequencing Analysis Software and 3730 Series Data Collection Software 4 filter .mob file choices to match the chosen .bcp file.



Notes



- b. In the Processed Data pane, select **True** or **Flat Profile**.

Option	Function
<input type="radio"/> True Profile	Used to display data as processed traces scaled uniformly so that the average height of peaks in the region of strongest signal is about equal to a fixed value. The profile of the processed traces will be very similar to that of the raw traces.
<input type="radio"/> Flat Profile	Used to display the data as processed traces scaled semi-locally so that the average height of peaks in any region is about equal to a fixed value. The profile of the processed traces is flat on an intermediate scale (> about 40 bases). Note: This option is applied to data that is analyzed with the KB™ basecaller only. If you use the ABI basecaller, the profile option reverts to True Profile.


- c. If desired, select one or more stop points for data analysis.
 d. Select your Threshold Quality option.

Option	Function
<input type="radio"/> Call all bases and assign QV	When using the KB™ basecaller, use this option to assign a base to every position, as well as the QV.
<input type="radio"/> Assign 'N' for bases with QV < 15	When using the KB™ basecaller, use this option to assign Ns to bases with QVs less than the set point. The QV is still displayed.

5. Select the **Mixed Bases** tab.

Note: This function is active with the KB™ Basecaller only.

- a. For data containing any mixed bases, select **Use Mixed Base Identification**.
 b. The User can set the secondary peak threshold, as a percentage of the primary peak, for consideration as a mixed base by the basecalling algorithm. Reaching this threshold is a necessary but not sufficient condition for arriving at a mixed base determination. Set the percentage by entering a value into the “= __%” field or by dragging the horizontal line above or below the 25% default setting.



Note: Do not use less than 15% as your detection limit.

Notes




6. Select the **Clear Range** tab.

Note: The clear range is the region of sequence that remains after excluding the low-quality or error-prone sequence at both the 5' and 3' ends.

Select one or more Clear Range methods. If you apply multiple methods, the smallest clear range results.


7. Click **OK** to save the protocol and close the Sequence Analysis Protocol Editor.



Editing and Deleting Analysis Protocols

Editing an Analysis Protocol

1. In the Analysis Protocols pane in the Analysis Protocol Manager, select the protocol you want to edit.
2. Click **Edit...**
3. Make changes in the General, Basecalling, Mixed Bases, and Clear Range tabs, as appropriate.
4. Click **OK** to save the protocol and close the Analysis Protocol Editor.



Notes




Deleting an Analysis Protocol

IMPORTANT! Do not delete an Analysis Protocol during a run while it is being used for that run. Autoanalysis by DNA Sequencing Analysis Software 6 is not performed if you do so. Also, you must first delete any plate records using the Analysis Protocol before you can delete or modify the Analysis Protocol for these plate records.

1. In the Analysis Protocols pane in the Analysis Protocol Manager, select the protocol you want to delete.
2. Click **Delete**. The Deletion Confirmation dialog box opens.
3. Click **Yes**.


Note: To reuse a plate after deleting the analysis protocol associated with it, either re-create the analysis protocol with the same name or assign the plate a unique plate name.




Exporting and Importing Analysis Protocols

Exporting an Analysis Protocol

1. In the Analysis Protocols pane in the Analysis Protocol Manager, select the protocol you want to export.
2. Click **Export**. A standard file export dialog box opens.
3. Navigate to the destination folder.
4. Click **Save**.




Notes _____



Importing an Analysis Protocol

1. In the Analysis Protocols pane in the Analysis Protocol Manager, select the protocol you want to import.
2. Click **Import** . a standard file export dialog box opens.
3. Click **Save**.




Results Groups


A Results Group is a component within Data Collection that organizes samples and certain user settings under a single name. It is called a Results Group because it is used to analyze, name, sort, and deliver samples that result from a run.

Creating a Results Group

1. In the navigation pane of the Data Collection Software, click **GA Instruments > Results Group**.
2. Click **New...**.
 The Results Group Editor window opens.



3. Select the General tab, then:
 - a. Type a Results Group Name. The name can be used in naming and sorting sample files. It must be unique (see [page 70](#) for a list of accepted characters).
 - b. (Optional) Type a Results Group Owner. The owner name can be used in naming and sorting sample files.
 - c. (Optional) Type a Results Group Comment.



Notes




4. Select the **Analysis** tab, then:

- a. Select **Sequencing Analysis** from the Analysis Type drop-down list.
- b. In the Analysis Actions section, select **Do Autoanalysis**, if you want your data automatically analyzed after a run by DNA Sequencing Analysis Software 6.

Note: Login ID and password are not required for Sequencing Analysis software.


Note: Autoanalysis by SeqScape® is no longer supported.



5. Select the **Destination** tab, then use the default destination or define a new location for data storage.

To use ...	Then ...
default location	skip to step 1
custom location	complete step 5b

Note: The Results Group Destination tab, and Data Collection Software in general, does not recognize remote storage locations unless they have been mapped to a local drive letter using the Map Network Drive feature of the operating system. Specify the mapped drive letter location in the Results Group Destination tab.



- a. Click **Use Custom Location**, then click **Browse...** to navigate to a different save location.
- b. Click **Test** to test the Location path name connection:
 - If it passes, “Path Name test successful” is displayed.

Notes _____



- If it fails, “Could not make the connection. Please check that the Path Name is correct.” is displayed. Click **Browse** then select a different location.

Sample File Destinations

Locations Where Sample Files Are Placed During Extraction:

- Default Destination, default folder naming: Data / instrument type / instrument name / run folder (No ProcessedData folder)
- Default Destination, custom folder naming: Data/top custom folder/subfolders, and so on.
- Custom Destination, default folder naming: Destination/instrument type/instrument name/run folder
- Custom Destination, custom folder naming: Destination/top custom folder/subfolders, and so on.

1. Select the **Naming** tab.

Use the Naming tab to customize sample file and run folder names.



Note: Sample name, run folder name, and path name, *combined*, can total no more than 250 characters. See page [page 70](#) for accepted characters.

The elements of the Naming tab are discussed in the following sections.

Sample File Name Format Pane

Follow the procedure below to complete the Sample File Name Format pane.



1. (Optional) In the Naming tab, select the **Prefix** box to type a prefix for the file name. Anything that you type here is shown in the Example line (see figure below).



Notes




2. Click the **Name Delimiter** list then select the symbol that will separate the Format elements in the file name (see step 3 below). You can select only one delimiter symbol.



3. Click the Format list, then select the components that you want in the sample name.


Note: Generally, all the samples from a single run are placed in the same run or results folder, so the name of every sample from a single run should be different from each other. However, most of the Format options are not different between samples, you need to take care to select at least one of the options that make the sample names unique within a run.




For example, if a unique identifier is not included in the name, a warning message is displayed. The Results Group makes the file name unique. As you select the elements for the file name, they are reflected in the Example line.

As you continue to select elements for the file name, additional elements are displayed.

Notes _____




The names of the Format elements are eventually shortened, but the Example field remains visible (up to 72 characters).



4. (Optional) Select the **Suffix** box, then type the suffix for the file name.

The File Extension field displays the file extension generated from the Analysis Type specified on the Analysis tab (page 86). For example, Sequencing Analysis produces sample files with an .ab1 extension.



Saving a Results Group

Click **OK** in any tab after you select all the elements within the Results Group.

Note: Even if you create a custom run folder location, a separate default run folder is generated that contains the log file.

Format Elements (Unique Identifiers)


Although you can save a results group by selecting a minimum of one Format element, selecting just the minimum may not provide enough information for you to identify the file or folder later.

Note: If you choose a non-unique file name, the software appends numbers (incrementally) before the file extension.

Notes



If you select elements from the Format lists that do not create unique Sample file or Run folder names, a warning message is displayed below the Example line (see next figure).



To remove the warning message and proceed within the Results Group Editor window, simply select a Format element that distinguishes one file from another (for example, the capillary number is unique but the instrument name is not).

Run Folder/Sub-Folder Name Format Pane

Follow the same steps described above for the Sample File Name Format pane (page [page 87](#)) to specify the run folder name within the run folder.

Notes _____






Importing and Exporting a Results Group

Results Groups can be imported from, or exported to, tab-delimited text files. This allows easy sharing of identical Results Groups between instruments.

Note: Importing Excel files is not supported.

Importing a Results Group





1. In the navigation pane of the Data Collection Software, select  **GA Instruments** >  **Results Group**.
2. Click . A standard File Import dialog box opens.
3. Navigate to the file you want to import.

Note: Import file type is .xml (extensible markup language).

4. Click .

Note: When you import or duplicate a Results Group, the software prompts you to type a name for the new Results Group and for the analysis application type.

Exporting a Results Group

1. In the navigation pane of the Data Collection Software, select  **GA Instruments** >  **Results Group**.
2. Click the Results Group name to select it.
3. Click . A standard file export dialog box opens with the chosen Results Group name.
4. Navigate to the location where you want to save the exported file.
5. Click .

Note: A name conflict occurs with a Results Group that already exists at the save location, the Results group can be duplicated to copy the settings into a similar Results Group without the risk of user error when copying it manually (see procedure below).

Notes _____



Duplicating a Results Group

1. Click the Results Group to select it.
2. Click **Duplicate** .


Note: When you import or duplicate a Results Group, the software prompts you to type a name for the new Results Group and for the analysis application type.

Notes _____



Creating and Completing a Sequencing Analysis Plate Record

1. In the navigation pane of the Data Collection Software, select **GA Instruments** > **ga3730** > **Plate Manager**.
2. Click **New...**. The New Plate Dialog dialog box opens.
3. In the New Plate Dialog:
 - a. Type a plate ID or barcode.
 - b. Type a name for the plate.
 - c. (Optional) Type a description for the plate.
 - d. Select your sequencing application in the Application drop-down list.
 - e. Select **96-well** or **384-well** in the Plate Type drop-down list.
 - f. Select **96-well** or **384-well** in the Plate Type drop-down list.
 - g. Select **heat seal** or **septa**.
 - h. Type a name for the owner and operator.
 - i. Click **OK**. The Sequencing Analysis Plate Editor opens.




Notes _____






Completing a Sequencing Analysis Plate Record

Note: Plate records can be imported and exported as tab-delimited files (.txt)

Note: Importing Excel files is not supported.



1. In the Sample Name column of a row, enter a sample name, then click the next cell. The value 100 is automatically displayed in the Priority column.
2. In the Comments column, enter any additional comments or notations for the sample.
3. In the Results Group 1 column, select a group from the drop-down list (see [page 85](#)).
4. In the Instrument Protocol 1 column, select a protocol from the drop-down list (see [page 74](#)).
5. In the Analysis Protocol 1 column, select a protocol from the drop-down list (see [page 80](#)).
6. To complete the rest of the plate record based on the samples loaded in your plate, do one of the following:
 - For the same samples and protocols – Select the entire row, then select **Edit > Fill Down Special** (see “[Fill Down Special](#)” on [page 96](#))
 - Based on the plate type (96- or 384-well) and capillary array (48 or 96 capillaries) you are using, select the appropriate fill down option:




Notes



- 96 capillary/96-well plate: **Fill Down**.
 - 48 capillary/96-well plate: **Fill down Special (48 Cap)**.
 - 96 capillary/384-well plate: **Fill down Special (96 Cap)**.
 - 48 capillary/384-well plate: **Fill down Special (48 Cap)**.
 - For the same samples and protocols – Select the entire row, then select **Edit > Fill Down**.
 - For the different samples and protocols, complete the plate editor manually.
7. If you want to do more than one run, select **Edit > Add Sample Run**.

Additional Results Group, Analysis Protocol, and Instrument Protocol columns are added to the right end of the plate record.

To add additional runs, select **Edit > Add Sample Run** again.



8. Complete the columns for the additional runs.
9. Click **OK**.

IMPORTANT! After clicking OK within the Plate Editor, the completed plate record is stored in the Plate Manager database, then the plate record can be searched for, edited, exported, or deleted.

Note: If multiple cells are selected for copying, select the same number of corresponding target cells before you execute the Paste command.

Note: The Plate Editor Copy and Paste functionality is supported only within one plate editor. To copy and paste the contents of one plate to another plate, use the “Duplicate..”. button on the Plate Manager dialog box.



Note: If you use the duplicate plate function, all the information in the plate to be duplicated must be valid. Otherwise, an empty plate is created.

Notes _____



Fill Down Special

The following table illustrates the Fill Down Special feature.

If You Choose ...	Then ...																																										
<p>Fill Down Special (48 Cap)</p> 	<p>The fill down pattern matches the 48-capillary load pattern.</p> <table border="1" data-bbox="516 506 703 1098"> <thead> <tr> <th>Well</th> <th>Sample Name</th> </tr> </thead> <tbody> <tr><td>A01</td><td>notMJD</td></tr> <tr><td>B01</td><td>notMJD</td></tr> <tr><td>C01</td><td>notMJD</td></tr> <tr><td>D01</td><td>notMJD</td></tr> <tr><td>E01</td><td>notMJD</td></tr> <tr><td>F01</td><td>notMJD</td></tr> <tr><td>G01</td><td>notMJD</td></tr> <tr><td>H01</td><td>notMJD</td></tr> <tr><td>A02</td><td>MJD</td></tr> <tr><td>B02</td><td>MJD</td></tr> <tr><td>C02</td><td>MJD</td></tr> <tr><td>D02</td><td>MJD</td></tr> <tr><td>E02</td><td>MJD</td></tr> <tr><td>F02</td><td>MJD</td></tr> <tr><td>G02</td><td>MJD</td></tr> <tr><td>H02</td><td>MJD</td></tr> <tr><td>A03</td><td>notMJD</td></tr> <tr><td>B03</td><td>notMJD</td></tr> <tr><td>C03</td><td>notMJD</td></tr> <tr><td>D03</td><td>notMJD</td></tr> </tbody> </table> <p>First Quadrant</p> <p>Second Quadrant</p>	Well	Sample Name	A01	notMJD	B01	notMJD	C01	notMJD	D01	notMJD	E01	notMJD	F01	notMJD	G01	notMJD	H01	notMJD	A02	MJD	B02	MJD	C02	MJD	D02	MJD	E02	MJD	F02	MJD	G02	MJD	H02	MJD	A03	notMJD	B03	notMJD	C03	notMJD	D03	notMJD
Well	Sample Name																																										
A01	notMJD																																										
B01	notMJD																																										
C01	notMJD																																										
D01	notMJD																																										
E01	notMJD																																										
F01	notMJD																																										
G01	notMJD																																										
H01	notMJD																																										
A02	MJD																																										
B02	MJD																																										
C02	MJD																																										
D02	MJD																																										
E02	MJD																																										
F02	MJD																																										
G02	MJD																																										
H02	MJD																																										
A03	notMJD																																										
B03	notMJD																																										
C03	notMJD																																										
D03	notMJD																																										
<p>Fill Down Special (96 Cap) *</p>  <p>* Especially useful for 384-well plates</p>	<p>The fill down pattern matches the 96-capillary load pattern.</p> <table border="1" data-bbox="516 1167 703 1766"> <thead> <tr> <th>Well</th> <th>Sample Name</th> </tr> </thead> <tbody> <tr><td>A10</td><td>12345</td></tr> <tr><td>B10</td><td>12345</td></tr> <tr><td>C10</td><td>12345</td></tr> <tr><td>D10</td><td>12345</td></tr> <tr><td>E10</td><td>12345</td></tr> <tr><td>F10</td><td>12345</td></tr> <tr><td>G10</td><td>12345</td></tr> <tr><td>H10</td><td>12345</td></tr> <tr><td>A11</td><td>12345</td></tr> <tr><td>B11</td><td>12345</td></tr> <tr><td>C11</td><td>12345</td></tr> <tr><td>D11</td><td>12345</td></tr> <tr><td>E11</td><td>12345</td></tr> <tr><td>F11</td><td>12345</td></tr> <tr><td>G11</td><td>12345</td></tr> <tr><td>H11</td><td>12345</td></tr> <tr><td>A12</td><td>12345</td></tr> <tr><td>B12</td><td>12345</td></tr> <tr><td>C12</td><td>12345</td></tr> </tbody> </table>	Well	Sample Name	A10	12345	B10	12345	C10	12345	D10	12345	E10	12345	F10	12345	G10	12345	H10	12345	A11	12345	B11	12345	C11	12345	D11	12345	E11	12345	F11	12345	G11	12345	H11	12345	A12	12345	B12	12345	C12	12345		
Well	Sample Name																																										
A10	12345																																										
B10	12345																																										
C10	12345																																										
D10	12345																																										
E10	12345																																										
F10	12345																																										
G10	12345																																										
H10	12345																																										
A11	12345																																										
B11	12345																																										
C11	12345																																										
D11	12345																																										
E11	12345																																										
F11	12345																																										
G11	12345																																										
H11	12345																																										
A12	12345																																										
B12	12345																																										
C12	12345																																										

Notes




Fill Down Special for a 48 Cap/96-Well Plate

The Fill Down Special function allows you to fill the plate record based on the load pattern of the capillary array that you are using.

To use the fill down special function:


1. In the Plate Manager, double-click the plate of interest to open the Plate Editor.
2. Type the sample name, complete all columns, then click-drag the entire row to select it.
3. Select **Edit > Fill Down Special (48 Cap)** to fill the plate record with the first load pattern.




Notes



- Click A02, type the name of sample 2, complete all columns, then click-drag the entire row to select it.



- Select **Edit > Fill Down Special (48 Cap)** to fill the plate record with the second load pattern.




Notes



Fill Down Special for a 96 Cap/384-well Plate

When you use the Fill Down Special (96 Cap) function on a 384-well plate, the fill-down pattern appears as in the adjoining illustration to the right.



Adding a Sample Run


By adding additional sample runs, you can run samples with different variables (different run modules, for example).

To add a sample run **Select Edit > Add Sample Run.**

- Results Group
- Instrument Protocol
- Analysis Protocol (sequencing only)

To run the plate(s), see [“Running the Instrument” on page 127.](#)

Note: When you add another sample run to a processed plate, confirm that all the information in the processed runs is valid. Otherwise, that data will not be validated, and a new sample run cannot be created.



Notes



Chapter 4 Setting Up the Software for DNA Sequencing

Fill Down Special

SequencingAnalysis Plate Editor

File Edit

Plate Name: 384 Operator: sc

Plate ID: 384 Owner: sc

Plate Sealing: Heat Sealing Scheduling: 1234

vWell	Sample Name	Comment	Results Group 1	Instrument Protocol 1	Analysis Protocol 1
A01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
B01					
C01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
D01					
E01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
F01					
G01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
H01					
I01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
J01					
K01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def
L01					
M01	sample		SeqA	RapidSeq	3730BDTv3-KB-Def

SequencingAnalysis Plate Editor

File Edit

Plate Name: Sample_10 Operator: m

Plate ID: Sample_10 Owner: m



Plate Sealing: Septa

vWell	Instrument Protocol 1	Analysis Protocol 1	Results Group 2	Instrument Protocol 2	Analysis Protocol 2
A01					
B01					
C01					
D01					
E01					
F01					
G01					
H01					
A02					
B02					
C02					
D02					
E02					
F02					
G02					
H02					
A03					
B03					
C03					
D03					
E03					
F03					
G03					

Description | Ok | Cancel

Notes _____

Setting Up the Software for Fragment Analysis



Notes



3730/3730xl Analyzer Data Collection and GeneMapper Software

IMPORTANT! Do not rename the computer after 3730 Series Data Collection Software is installed. Doing so causes the 3730 Series Data Collection Software to malfunction.

File-Naming Convention

Some alphanumeric characters are not valid for user names or file names. The invalid characters are below:

spaces \ / : * ? " < > |

IMPORTANT! An error message is displayed if you use any of these characters. You must remove the invalid character to continue.

Note: Autoanalysis by GeneMapper[®] Software is no longer supported.

Data Analysis

For information on data analysis, refer to *GeneMapper[®] Software 5 Online Help* (Part no. 4474202)

Fragment Analysis and Data Collection

When GeneMapper[®] software is installed on a computer that has 3730 Series Data Collection Software, you can access GeneMapper[®] through the Results Group Editor (see [page 114](#)):

- GeneMapper-Generic
- GeneMapper-<Computer Name>

GeneMapper-Generic

GeneMapper-Generic enables you to generate .fsa files. When completing the Sample Sheet, you need to fill in basic information for Data Collection to complete the run; all other GeneMapper[®] software related fields are text entries. This is useful if you are using other software applications for analysis. This is also useful if you choose to analyze your samples in GeneMapper[®] software on another computer, but do not have the same entries in the GeneMapper[®] software database stored on the Data Collection computer. For example, if you have a customized size standard definition on the other GeneMapper[®] software computer, you can type in that size standard name in the size standard text field and it will populate that column in your GeneMapper[®] software project.

Notes _____




**GeneMapper-
<Computer
Name>**

GeneMapper-<Computer Name> permits the Size Standard, Analysis Method, and Panel columns in the Sample Sheet window to be read directly from the GeneMapper® software database. These components must be created in GeneMapper® software prior to setting up the plate record for a run. There is no way to create a new entry for these columns once inside the plate editor dialog box. If you create a new GeneMapper® software component while the plate record dialog box is open, the columns will not update. The plate record must be closed and reopened to update the GeneMapper® software components.

Notes _____



**Workflow with GeneMapper-
<Computer Name> option
Using GeneMapper®
Software**



Notes _____



GeneMapper® Software Plate Records

Overview Plate records are data tables in the instrument database that store information about the plates and the samples they contain. A plate record contains:

- Plate name, type, and owner
- Position of the sample on the plate (well number)
- Comments about the plate and about individual samples
- Dye set information (in instrument protocol)
- Name of the run module. Run modules specify information about how samples are run (in instrument protocol)

A plate record is similar to a sample sheet or an injection list that you may have used with other Applied Biosystems® instruments.

When to Create a Plate Record You must create a plate record for each plate of samples for:

- Spectral calibrations
- Fragment analysis

Note: A plate record must be created in advance of the first run. Then, plate records can be created, and plates added to the stacker, while a run is in progress.

Parameters	Description	See Page
Instrument protocol	Contains everything needed to run the instrument.	99
Results group	Defines the file type, the file name, and file save locations that are linked to sample injections.	104


IMPORTANT! For data collection and analysis to be successful, each run of samples must have an Instrument Protocol and a Results Group assigned within a plate record.

Note: Autoanalysis by GeneMapper® is no longer supported.

Notes _____




Components of a GeneMapper® Software Plate Record



Notes _____



Descriptions for numbers 1–10 are in the following table.



Default is one sample run. To add additional runs, see [page 115](#).

The following table describes columns 1-10 inserted in a plate record for a fragment analysis run (see the preceding figure).

Table 5-1 Components of the plate record

Column	Description
1. Sample Name	Name of the sample
2. Comment	(Optional) Comments about the sample
3. Sample Type	Use to identify the sample as Sample, Positive Control, Allelic Ladder, or Negative Control.
4. Size Standard IMPORTANT! For GeneMapper-<Computer Name> ONLY: Size Standard, Panel, and Analysis Method must be created in GeneMapper® software before creating a new plate in order to make them available in Data Collection Software	<ul style="list-style-type: none"> (Optional) GeneMapper-Generic: Manually enter size standards in the text field GeneMapper-<Computer Name>: Select a saved size standard from the drop-down list

Notes



Table 5-1 Components of the plate record

Column	Description
5. Panel IMPORTANT! GeneMapper-<Computer Name> ONLY: Size standard, panel, and analysis method must be created in GeneMapper software before creating a new plate	For <ul style="list-style-type: none"> • (Optional) GeneMapper-Generic: Manually enter panels in the text field* • GeneMapper-<Computer Name>: Select a saved panel from the drop-down list
6. Analysis Method IMPORTANT! For GeneMapper <Computer Name> ONLY: Size standard, panel, and analysis method must be created in GeneMapper® software before creating a new plate	<ul style="list-style-type: none"> • (Optional) GeneMapper-Generic: Manually enter analysis methods in the text field* • GeneMapper-<Computer Name>: Select a saved analysis method from the drop-down list
7. Snp	Optional field, typically left blank
8. 3 User-defined columns	Optional text entries
9. Results group	Some options: <ul style="list-style-type: none"> • New: Opens the Results Group Editor dialog box • Edit: Opens the Results Group Editor dialog box for the results group listed in the cell • None: Sets the cell to have no selected results group • Select one of the available Results groups from the list <p>Note: You must have a results group selected for each sample entered in the Sample Name column.</p> <p>See, “Results Groups” on page 114.</p>
10. Instrument protocol	<ul style="list-style-type: none"> • New: Opens the Protocol Editor dialog box. • Edit: Opens the Protocol Editor dialog box for the instrument protocol listed in the cell. • None: Sets the cell to have no selected protocol. • List of Instrument Protocols: In alpha-numeric order. <p>Note: You must have an instrument protocol selected for each sample entered in the Sample Name column.</p> <ul style="list-style-type: none"> • See, “Instrument Protocols” on page 109.

Notes _____



Creating Required Settings for Fragment Analysis

If the Settings Already Exist


If the appropriate data collection and fragment analysis files have been created, go to “[Creating and Completing a GeneMapper® Software Plate Record](#)” on page 121.

Instrument Protocols

An instrument protocol contains all the settings needed to run the instrument. An instrument protocol contains the protocol name, type of run, run module, and dye set.

Creating an Instrument Protocol


1. In the navigation pane of the Data Collection Software, select **GA Instruments** > **ga3730** > **Protocol Manager**.




Notes



- In the Instruments Protocols section, click **New...**. The Protocol Editor opens.




- Complete the Protocol Editor:
 - Type a name for the protocol.
 - (Optional) Type a description for the protocol.
 - Select **Regular** in the Type drop-down list.



- Select **GeneMapper36_POP-7™**.
- Select **G5**.
- Click **OK**.

Importing an Instrument Protocol

- In the Protocol Editor window select **Import** in the Instrument Protocols pane, if you want to use an existing instrument protocol.




Notes _____



- Navigate to the protocol you want to import.

Note: Import file type is .xml (extensible markup language).

- Double-click the protocol to import it.
- The imported files are displayed alphabetically in the Instrument Protocol pane.



Fragment Analysis Run Modules

Select one run module:

Run Module	Capillary Length
HTSNP36_POP-7™_V3 (SNaPshot®)	36 cm
HTSNP50_POP-7™_V3 (SNaPshot®)	50 cm
GeneMapper36_POP-7™	36 cm
GeneMapper50_POP-7™	50 cm
GS1200LIZ_36_POP-7™	36 cm
GS1200LIZ_50_POP-7™	50 cm

Notes _____




Customizing Run Modules

If you need to modify default run modules to suit your particular needs:

1. Select **GA Instrument**
 > **ga3730** > **Module Manager**.
2. Click .
3. Select a template module as a basis for the new module.
4. Change to the desired module parameters using the table below as a guide.

Note: You cannot edit a default module installed with 3730/3730xL Analyzer Data Collection Software.



Choose module template from the drop-down menu (step 3).

Notes _____



The Run Module Parameters that you can edit:

Parameter Name	Range	Description
Oven_Temperature	18–70 C	Temperature setting for main oven throughout run.
PreRun_Voltage	0–15 kV	Pre run voltage setting before sample injection.
PreRun Time	1–1800 sec	Prerun voltage time.
Injection_Voltage	0–15 kV	Injection voltage setting for sample injection.
Injection_Time	1–90 sec	Sample injection time.
First_ReadOut_time	100–16000 millisc	The interval of time for a data point to be produced. First_ReadOut_time should be equal to Second_ReadOut_time.
Second_ReadOut_Time	100–16000 millisc	The interval of time for a data point to be produced. Second_ReadOut_time should be equal to First_ReadOut_time.
Run_Voltage	0–15 kV	Final run voltage.
Voltage_Number_Of_Steps	0–100 steps	Number of voltage ramp steps to reach Run_Voltage. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Voltage_Step_Interval	0–180 sec	Dwell time at each voltage ramp step. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Voltage_Tolerance	0.1–6 kV	Maximum allowed voltage variation. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel. If the instrument goes beyond tolerance and shuts off, contact Life Technologies tech support.
Current_Stability	0–2000 microA	Maximum allowed electrophoresis current variation. Current fluctuations above this value will be attributed to air bubbles in system and the voltage automatically powered off. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Ramp_Delay	1–1800 sec	Delay During Voltage Ramp. Life Technologies recommends that you do not change this value unless advised otherwise by Life Technologies support personnel.
Data_Delay	1–1800 sec	Time from the start of separation to the start of data collection.
Run_Time	300–14000 sec	Duration data is collected after Ramp_Delay.

Notes _____




Results Groups

A Results Group is a component within Data Collection that organizes samples and certain user settings under a single name. A Results Group is used to prepare samples for analysis and to name, sort, and deliver samples that result from a run.

Creating a Results Group


1. In the navigation pane of the Data Collection Software, select **GA Instruments > Results Group**.



2. Click **New**. The Results Group Editor window opens.

3. Select the **General** tab:

- a. Type a Results Group Name. The name can be used in naming and sorting sample files. It must be unique (see page for a list of accepted characters).
- b. (Optional) Type a Results Group Owner. The owner name can be used in naming and sorting sample files.
- c. (Optional) Type a Results Group Comment.




4. Skip the **Analysis** tab, because autoanalysis by GeneMapper® is no longer supported.

Notes _____



5. Select the **Destination** tab, then use the default destination or define a new location for data storage. To use a:
- Default location – Skip to step 6.
 - Custom location – Complete step a and step b below.
- a. Click **Use Custom Location**, then click **Browse...** to navigate to a different save location.
- b. Click **Test** to test the Location path name connection:
- If the test passes, “Path Name test successful,” displays.
 - If the test fails, “Could not make the connection. Please check that the Path Name is correct,” displays. Click **Browse**, then select a different location.

Note: The Results Group Destination tab, and Data Collection Software in general, does not recognize remote storage locations unless they have been mapped to a local drive letter using the Map Network Drive feature of the operating system. Specify the mapped drive letter location in the Results Group Destination tab.



Sample File Locations

Locations Where Sample Files Are Placed During Extraction:

- Default Destination, default folder naming: Data / instrument type / instrument name / run folder (No ProcessedData folder)
- Default Destination, custom folder naming: Data/top custom folder/subfolders, and so on.
- Custom Destination, default folder naming: Destination/instrument type/instrument name/run folder
- Custom Destination, custom folder naming: Destination/top custom folder/subfolders, and so on.


Notes



6. Select the **Naming** tab. Use the Naming tab to customize sample file and run folder names.

Note: Sample name, run folder name, and path name, *combined*, can total no more than 250 characters. See [page 102](#) for accepted characters.

The elements of the Naming tab are discussed in the following sections, see [page 117](#).



7. Skip the **Automated Processing** tab, because Autoanalysis by GeneMapper® is no longer supported.

8. Click **OK** to save the Results Group.



Notes _____



Sample File Name Format Pane

To complete the Sample File Name Format pane:


1. (Optional) Select the **Prefix** box then type a prefix for the file name. Anything that you type here is shown in the Example line (see the following graphic).
2. Click the **Name Delimiter** list choose the symbol that will separate the Format elements in the file name (see step 3). You can only choose one delimiter symbol.




3. Click the **Format** list and then select the components that you want in the sample name.

Generally, all the samples from a single run are placed in the same run or results folder, so the name of every sample from a single run should be different. Most of the Format options are not different between samples, so you need take care to select at least one of the options that makes the sample names unique within a run.

For example, if a unique identifier is not included in the name, a warning message displays. The Results Group makes the file name unique. As you select the elements for the file name, they are reflected in the Example line.



Note: An additional drop-down list of formats is displayed after you select a format option.



Notes



The names of the Format elements are eventually shortened, but the Example field remains visible (up to 72 characters).



Note: To view the shortened format elements, place the cursor on the edge of the window until it turns into a double-arrow. Drag the arrow to expand the window horizontally.

4. (Optional) Click the **Suffix** box then type the suffix for the file name.

The File Extension field displays the file extension generated from the Analysis Type specified on the Analysis tab (page 114). For example, fragment analysis produces sample files with an .fsa extension.

Run Folder/Sub-Folder Name Format Pane

Follow the same steps described above for the Sample File Name Format pane (page 117) to change the sub-folder name within the run folder.



Notes _____




Format Elements (Unique Identifiers)

Although you can select a minimum of one Format element for the Sample file and Run folder names to save a Results Group, selecting the minimum may not provide enough information for you to identify the file or folder later.

Note: If you choose a non unique file name, the software automatically appends numbers (incrementally) before the file extension.

If you select elements from the Format lists that do not create unique Sample file or Run folder names, a warning message displays below the Example line (see the following figure).



To remove the warning message and proceed within the Results Group Editor window, select a Format element that distinguishes one file from another (for example, the capillary number is unique but the instrument name is not).

Importing and Exporting a Results Group

Results Groups can be imported from, or exported to, tab-delimited text files to allow easy sharing of identical Results Groups between instruments.

Note: Importing Excel files is not supported.

Importing a Results Group

1. In the navigation pane of the Data Collection Software, select **GA Instruments > Results Group**.
2. Click **Import**. A standard File Import dialog box opens.
3. Navigate to the file you want to import.

Note: Import file type is .xml (extensible markup language).

Notes _____



Chapter 5 Setting Up the Software for Fragment Analysis

Creating Required Settings for Fragment Analysis

4. Click **Open** .

Note: When you duplicate a Results Group, the software prompts you to type a name for the new Results Group and for the analysis application type.

Exporting a Results Group

1. In the navigation pane of the Data Collection Software, select **GA Instruments > Results Group**.
2. Select the Results Group name.
3. Click **Export** . A standard file export dialog box opens, displaying the chosen Results Group name.
4. Navigate to where you want to save the exported file.
5. Click **Save** .

Note: If a results group with the same name already exists at the save location, you can duplicate the results groups to copy settings into a similar results group without the risk of user error.

Duplicating a Results Group

1. Click the results group to select it.
2. Click **Duplicate** .

Note: When you duplicate a results group, the software prompts you to type a name for the new Results Group and for the analysis application type.


Notes _____



Creating and Completing a GeneMapper® Software Plate Record

Creating the GeneMapper® Software Plate Record

1. In the navigation pane of the Data Collection Software, select
GA Instruments > **ga3730** > **Plate Manager**.
2. Click **New...**. The New Plate Dialog dialog box opens.
3. Complete the information in the New Plate Dialog:
 - a. Type a plate ID.
 - b. Type a name for the plate.
 - c. (Optional) Type a description for the plate.
 - d. Select your GeneMapper application in the Application drop-down list.
 - e. Select **96-well** or **384-well** in the Plate Type drop-down list.
 - f. Schedule the plate. For more information, see [“Scheduling Runs” on page 133](#).
 - g. Select **Heat Sealing** or **Septa**.
 - h. Type a name for the owner and the operator.
 - i. Click **OK**. The GeneMapper Software Plate Editor opens.



Completing a GeneMapper Software Plate Record

1. In the Sample Name column of a row, enter a sample name, then click the next cell.
2. In the Comment column, enter any additional comments or notations for the sample.
3. In the Sample Type column, select a sample type from the drop-down list.
4. In the Size Standard column, select a size standard from the drop-down list.

vWell	Sample Name	Comment	Sample Type
A01			
B01			
C01			
D01			
E01			
F01			

Notes




5. In the Panel column, select a panel from the drop-down list.
6. In the Analysis Method column, select a method from the drop-down list.
7. (Optional) In the Snp Set column, select a SNP set from the drop-down list.
8. Enter text for User-Defined columns 1 to 3.
9. In the Results Group 1 column, select a group from the drop-down list.
10. In the Instrument Protocol 1 column, select a protocol from the drop-down list.

4	5	6	7
Size Standard	Panel	Analysis Method	Snp Set

8			9	10
User-Defined 1	User-Defined 2	User-Defined 3	Results Group ▶	Instrument Protocol ▶

11. To complete the rest of the plate record based on the samples loaded in your plate, do one of the following:

- For the same samples and protocols – Select the entire row, then select **Edit > Fill Down Special**. For more information see, “Filling Down the Plate Record” on page 124.
- Based on the plate type (96- or 384-well) and capillary array (48, 50, or 96 capillaries) you use–Select the appropriate fill down option:
 - 96 capillary/96-well plate: **Fill Down**
 - 48 capillary/96-well plate: **Fill down Special (48 Cap)**
 - 96 capillary/384-well plate: **Fill down Special (96 Cap)**
 - 48 capillary/384-well plate: **Fill down Special (48 Cap)**
- For the different samples and protocols, complete the plate editor manually.




Notes _____



12. To do more than one run, select **Edit > Add Sample Run.**

Additional Results Group and Instrument Protocol columns are added to the right end of the plate record.

To add additional runs select **Edit > Add Sample Run**, again (for more information see, [“Adding a Sample Run” on page 126](#)).



13. Complete the columns for the additional runs.

14. Click to save, then close the plate record.

IMPORTANT! After clicking OK within the Plate Editor, the completed plate record is stored in the Plate Manager database. After the plate record is in the Plate Manager database, the plate record can be searched for, edited, exported, or deleted.

Note: If multiple cells are selected for copying, select the same number of corresponding target cells before you execute the Paste command.

Note: The Plate Editor Copy and Paste functionality is supported only within one plate editor. To copy and paste the contents of one plate to another plate, use the “Duplicate...” button on the Plate Manager dialog box.



Note: If you use the duplicate plate function, all the information in the plate to be duplicated must be valid. Otherwise, an empty plate is created.

Notes _____



Filling Down the Plate Record

The Fill Down Special function allows you to fill a plate record based on the load pattern of the capillary array that you use, as shown in the following table.

If You Choose ...	Then ...																																																		
<p>Fill Down Special (48 Cap)</p> 	<p>The fill down pattern matches the 48-capillary load pattern.</p> <table border="1"> <thead> <tr> <th>Well</th> <th>Sample Name</th> </tr> </thead> <tbody> <tr><td>A01</td><td>notMJD</td></tr> <tr><td>B01</td><td>notMJD</td></tr> <tr><td>C01</td><td>notMJD</td></tr> <tr><td>D01</td><td>notMJD</td></tr> <tr><td>E01</td><td>notMJD</td></tr> <tr><td>F01</td><td>notMJD</td></tr> <tr><td>G01</td><td>notMJD</td></tr> <tr><td>H01</td><td>notMJD</td></tr> <tr><td>A02</td><td>MJD</td></tr> <tr><td>B02</td><td>MJD</td></tr> <tr><td>C02</td><td>MJD</td></tr> <tr><td>D02</td><td>MJD</td></tr> <tr><td>E02</td><td>MJD</td></tr> <tr><td>F02</td><td>MJD</td></tr> <tr><td>G02</td><td>MJD</td></tr> <tr><td>H02</td><td>MJD</td></tr> <tr><td>A03</td><td>notMJD</td></tr> <tr><td>B03</td><td>notMJD</td></tr> <tr><td>C03</td><td>notMJD</td></tr> <tr><td>D03</td><td>notMJD</td></tr> <tr><td>E03</td><td>notMJD</td></tr> <tr><td>F03</td><td>notMJD</td></tr> <tr><td>G03</td><td>notMJD</td></tr> <tr><td>H03</td><td>notMJD</td></tr> </tbody> </table> <p>First Quadrant</p> <p>Second Quadrant</p>	Well	Sample Name	A01	notMJD	B01	notMJD	C01	notMJD	D01	notMJD	E01	notMJD	F01	notMJD	G01	notMJD	H01	notMJD	A02	MJD	B02	MJD	C02	MJD	D02	MJD	E02	MJD	F02	MJD	G02	MJD	H02	MJD	A03	notMJD	B03	notMJD	C03	notMJD	D03	notMJD	E03	notMJD	F03	notMJD	G03	notMJD	H03	notMJD
Well	Sample Name																																																		
A01	notMJD																																																		
B01	notMJD																																																		
C01	notMJD																																																		
D01	notMJD																																																		
E01	notMJD																																																		
F01	notMJD																																																		
G01	notMJD																																																		
H01	notMJD																																																		
A02	MJD																																																		
B02	MJD																																																		
C02	MJD																																																		
D02	MJD																																																		
E02	MJD																																																		
F02	MJD																																																		
G02	MJD																																																		
H02	MJD																																																		
A03	notMJD																																																		
B03	notMJD																																																		
C03	notMJD																																																		
D03	notMJD																																																		
E03	notMJD																																																		
F03	notMJD																																																		
G03	notMJD																																																		
H03	notMJD																																																		
<p>Fill Down Special (96 Cap) *</p>  <p>* Especially useful for 384-well plates</p>	<p>The fill down pattern matches the 96-capillary load pattern.</p> <table border="1"> <thead> <tr> <th>Well</th> <th>Sample Name</th> </tr> </thead> <tbody> <tr><td>A10</td><td>12345</td></tr> <tr><td>B10</td><td>12345</td></tr> <tr><td>C10</td><td>12345</td></tr> <tr><td>D10</td><td>12345</td></tr> <tr><td>E10</td><td>12345</td></tr> <tr><td>F10</td><td>12345</td></tr> <tr><td>G10</td><td>12345</td></tr> <tr><td>H10</td><td>12345</td></tr> <tr><td>A11</td><td>12345</td></tr> <tr><td>B11</td><td>12345</td></tr> <tr><td>C11</td><td>12345</td></tr> <tr><td>D11</td><td>12345</td></tr> <tr><td>E11</td><td>12345</td></tr> <tr><td>F11</td><td>12345</td></tr> <tr><td>G11</td><td>12345</td></tr> <tr><td>H11</td><td>12345</td></tr> <tr><td>A12</td><td>12345</td></tr> <tr><td>B12</td><td>12345</td></tr> <tr><td>C12</td><td>12345</td></tr> </tbody> </table>	Well	Sample Name	A10	12345	B10	12345	C10	12345	D10	12345	E10	12345	F10	12345	G10	12345	H10	12345	A11	12345	B11	12345	C11	12345	D11	12345	E11	12345	F11	12345	G11	12345	H11	12345	A12	12345	B12	12345	C12	12345										
Well	Sample Name																																																		
A10	12345																																																		
B10	12345																																																		
C10	12345																																																		
D10	12345																																																		
E10	12345																																																		
F10	12345																																																		
G10	12345																																																		
H10	12345																																																		
A11	12345																																																		
B11	12345																																																		
C11	12345																																																		
D11	12345																																																		
E11	12345																																																		
F11	12345																																																		
G11	12345																																																		
H11	12345																																																		
A12	12345																																																		
B12	12345																																																		
C12	12345																																																		

To use the fill the plate record based on the 48 capillary load pattern:

1. In the Plate Editor, complete the sample information in a row within the first quadrant you want to fill.
2. Select the entire row.
3. Select **Edit > Fill Down Special (48 Cap)** to fill the quadrant.

Notes _____



4. Click position A02, type the sample information, then select the entire row.

GeneMapper Plate Editor

File Edit

Plate Name: GeneMapper Operator: MJD
 Plate ID: GeneMapper Owner: MJD
 Plate Sealing: Septa

Well	Sample Name	Comment	Sample Type	Size Standard	Panel	Analysis Method	Snp Set	User-Defined 1	User-Defined 2	User-Defined 3	Results Group	Instrument Protocol
A01	a										GM	GeneMapper
B01	a										GM	GeneMapper
C01	a										GM	GeneMapper
D01	a										GM	GeneMapper
E01	a										GM	GeneMapper
F01	a										GM	GeneMapper
G01	a										GM	GeneMapper
H01	a										GM	GeneMapper
A02												
B02												
C02												
D02												
E02												
F02												
G02												
H02												
A03	a										GM	GeneMapper
B03	a										GM	GeneMapper
C03	a										GM	GeneMapper
D03	a										GM	GeneMapper
E03	a										GM	GeneMapper
F03	a										GM	GeneMapper
G03	a										GM	GeneMapper
H03	a										GM	GeneMapper
A04												
B04												

Description


5. Select **Edit > Fill Down Special (48 Cap)** to fill the second quadrant (see the preceding figure).

Notes



Filling Down a 96-Cap/384-well Plate Record

When you use the Fill Down Special (96-Cap) feature on a 384-well plate, the fill down pattern appears as shown in the following figure.



Adding a Sample Run

By adding additional sample runs, you can run samples that have different variables (different run modules, for example).


Adding a sample run opens an additional:

- Results group
- Instrument protocol

To add a sample run, select **Edit > Add Sample Run**.

To run the plate(s), see [“Running the Instrument” on page 127](#).


Note: When you add another sample run to a processed plate, confirm that all the information in the processed runs is valid. Otherwise, that data will not be validated, and a new sample run cannot be created.



Notes



Running the Instrument




Notes _____




Working with Plate Assemblies


Plate Assembly Components

 **WARNING** Do not use warped or damaged plates.

Materials Required for Each Septa Assembly:

- Plate retainer
- Plate septa
- Sample plate
- Base plate

 **WARNING** Use only *black* plate bases with septa-sealed plates.





Notes _____



Materials Required for Each Heat-Sealed Assembly

- Plate retainer
- Heat seal film
- Sample plate
- Base plate

 **WARNING** Use only *gray* plate bases with heat-sealed plates.



Heat Seal Film Guidelines



- Use 3-mil heat seal film (Catalog no. 4337570) which is 3-mil before and 1-mil after, heating.
- *Do not* use heat seal film that is thicker than 1-mil, after heating, on the 3730/3730xl DNA Analyzer.
- *Do not* use heat-seal film containing adhesives or metals because they may damage the instrument's piercing needles

Notes _____




Preparing a Septum-Sealed Plate Assembly

1. Seal the plate:
 - a. Place the plate on a clean, level surface.
 - b. Inspect septa weekly and be sure to replace any that are worn or discolored.
 - c. Lay the septum flat on the plate.
 - d. Align the holes in the septa strip with the wells of the plate, then firmly press downward onto the plate.
2. To prevent damage to the capillary array, inspect the plate and septa to verify that the septum fits snugly and flush on the plate.



3. Assemble the plate assembly:
 - a. Place the sample plate into the plate base.
 - b. Snap the plate retainer onto the plate and plate base.
 - c. Make sure when you assemble a plate that the retainer clip is flush with the plate base. A simple way to ensure that they are flush is to run your finger along the edge.




Notes _____



4. Verify that the holes of the plate retainer and the septa strip are aligned. If not, reassemble the plate assembly (see [step 3](#)).

IMPORTANT! Damage to the array tips occurs if the plate retainer and septa strip holes do not align correctly.




Notes _____




Placing Plate Assemblies into the Instrument

1. Open the stacker drawer.
2. Open the door of the In Stack tower.




3. Place the plate assemblies into the stacker in any order, making sure that each plate is oriented so that the notched corner of the plate assembly is at the rear right corner of the stacker.

IMPORTANT! Do not place more than 16 plates in the stacker.



4. Close the metal In Stack tower door.
5. Close the Stacker drawer.




Notes _____




Scheduling Runs

In the navigation pane of the Data Collection Software, select
GA Instruments > ga3730 > instrument name > Run Scheduler.



384-Well Plate Mapping and Default Run Scheduling

Samples within a plate are run in the order of their well designations. For example, a default 384-well injection pattern looks like the following:



- Quadrant 1: wells A1, C1, E1, G1...
- Quadrant 2: wells B1, D1, F1, H1...
- Quadrant 3: wells A2, C2, E2, G2...
- Quadrant 4: wells B2, D2, F2, H2...

Notes



- Plates that contain samples in a single quadrant and with more than one instrument protocol specified run all the protocols in the order in which they appear on the plate record before the next quadrant is run.

Note: The analysis module of a sample does not affect the order in which the sample quadrant runs.

Default Run Priorities and Load Positions

For information on setting up a plate record for:

- Sequencing – See [page 70](#).
- Fragment analysis – See [page 105](#).

The following table indicates the default run priorities and load positions.

Number of Capillaries	Plate Size	Run Priority	Quadrant	First Load Position
96	384-well	1	Q1	Well A1
		2	Q2	Well B1
		3	Q3	Well A2
		4	Q4	Well B2
48	96-well	1	Q1, load 1	Well A1
			Q1, load 2	Well A2
48	384-well	1	Q1 , load 1	Well A1
			Q1 , load 2	Well A3
		2	Q2 , load 1	Well B1
			Q2 , load 2	Well B3
		3	Q3 , load 1	Well A2
			Q3 , load 2	Well A4
		4	Q4 , load 1	Well B2
			Q4 , load 2	Well B4

Note: When using a 384-well plate and a 48-capillary array, you can change the run order of the main quadrant (**bold** numbers above) but not the load numbers.

Notes _____




Globally Modifying a Run Schedule


You can change the run order of quadrants and then apply it to all 384-well plates.

To modify the run order for all 384-well plates:

1. Click your instrument name in the navigation pane.
2. Select **Instrument > Scheduling Preference**.
The Default 384 well scheduling preference dialog box opens.
3. Select the quadrant priority (run order) from the Quadrant list.



You can select any run order. The example to the right shows a 4-3-2-1 quadrant priority (run order). With a 384-well and a 96-capillary array, the samples run in the order B2, A2, B1, A1...



Locally Modifying a Run Schedule

To locally modify the run order of quadrants within a single 384-well plate:

1. In the Plate Manager, click **New Plate**.
Note: For information about the Plate Manager, see [page 93](#) for sequencing, and [page 121](#) for fragment analysis.
2. Select **384-Well** from the Plate Type list.
The Scheduling box is activated.

Notes _____



3. Type the run priority in the Scheduling box.
4. Click **OK**.

Type run priorities here

New Plate Dialog

ID (Barcode): test

Name: test

Description:

Application: GeneMapper-Generic

Plate Type: 384-Well

Scheduling: 1234

Plate Sealing: Heat Sealing

Owner Name: user

Operator Name: user

OK Cancel

	1	2
A	1	3
B	2	4


Notes _____




Default Load Maps

Refer to the following load maps for different sized arrays and sample plates.

96-Well Plate, 48 Capillaries



96-Well Plate, 96 Capillaries




Notes _____




384-Well Plate, 48 Capillaries


First quadrant pickup




Second quadrant pickup



Third quadrant pickup



Fourth quadrant pickup




Notes




384-Well Plate, 96 Capillaries


First quadrant pickup




Second quadrant pickup



Third quadrant pickup



Fourth quadrant pickup



For a 384-well plate, injections are made from every other well and every other row. A full 384-well plate requires 4 runs for a 96-capillary array, and 8 runs for a 48-capillary array, to inject all the samples once.

Notes



Barcode Readers



CAUTION ELECTRICAL HAZARD. Power off the instrument and the computer before connecting an external barcode reader to the instrument.

Internal Barcode Reader

The 3730/3730x/ Analyzer internal barcode reader supports the following formats:


- Code 128
- Code 39
- Code 93
- LOGMARS
- EAN-8

Note: All Applied Biosystems® barcoded plates for the 3730/3730x/ Analyzer use code 128 format.

Note: The barcode reader cannot read spaces or the characters \ / : * ? " < > |.

External Barcode Readers

KEYENCE BL-80VE




An external barcode reader can also be used with the 3730/3730x/ Analyzer. The KEYENCE BL-80VE (see the preceding photo) connects to the instrument computer keyboard. With this reader, you can scan barcodes into any text box in the Data Collection software.

Notes _____



KEYENCE 80RKE



Another option is the KEYENCE 80RKE which you connect to the instrument serial port. With this reader, you can scan barcode information only into specific text boxes within the Data Collection Software.

Note: The 80RE is not supported for the 3730 or 3730x/ DNA Analyzers.

Notes _____



Running the Instrument: Manual vs Auto Mode

Accessing Modes You can schedule a run or runs using either manual mode or auto mode. Both modes are described in the following sections. Access either mode by selecting in the navigation pane:

Run Scheduler > Instrument > Instrument Name > Run mode (Auto or Manual)


Note: You must be in the Run Scheduler view to see the instrument run mode menu.

Manual Mode Features

- Plates can be added to the stacker individually and in order; runs are scheduled in the order the plates are in the stack.
- The internal reader is not necessary to link plates to plate records in the local database.
- Plates do not need to have a barcode.

Scheduling Runs Using Manual Mode (Default)


1. In the navigation pane, select **Instrument > Instrument Name > Manual mode**.
2. Click **Search** in the Run Scheduler to search for plate record(s).




The **Add Plates to In Stack** dialog box opens.

3. Type the name of the plate(s) or scan the plate ID, then click **Search**.

Notes _____



Barcode search



Advanced search

- Select the run(s) to add, then click **Add** to add the plate record(s) to the Input Stack in the order in which you want them to run.



- Click **Done** to close the Add Plates to In Stack dialog box.



- Physically stack the plates in the In Stack in order. The bottom plate runs first.

IMPORTANT! The order of the plate record must match the stack order of the plates in the In Stack. If the order does not match, processed runs have the wrong plate record information.

Note: You can assign more plates in the Run Scheduler than are actually available in the stacker.

- Click (**Run**).

As the plates are retrieved by the autosampler, they are run in the order they were placed in the In Stack.

Notes _____



Auto Mode Features

- Plates must have barcodes.
- an internal barcode reader is necessary to link plates to plate records in the local database.
- You can add plates to the In Stack in any order.
- Plates can be added or removed during instrument operation.

To schedule runs using the Auto mode:


1. Select **Run Scheduler > Instrument Name > Auto mode**.

Notice that the Search, Up, and Down buttons are not available (as they are in Manual mode). Also, the Add Plate (Scan or Type Plate ID) option is not available in Auto mode.

2. Physically place plates in the In Stack in any order. Remember that the bottom plate runs first and the top plate runs last.

3. Click  (Run).

As the plates are retrieved by the autosampler, plate barcodes are scanned and their plate records are associated with those stored in the local data collection database.




Notes _____



Starting the Run

1. Verify that the active spectral calibration matches your dye set and capillary array length.
2. If you want to review the run schedule before beginning the run, click
GA Instruments > ga3730 >
instrument name > Run Scheduler
3. Select the green button in the toolbar.
The Processing Plates dialog box opens.
4. Click **OK**.



5. The software automatically checks the:
 - Capillary array length and polymer type in the Instrument Protocol column of the plate record against the capillary array length and polymer type
 - Available space in the database and in drive E

If the database or drive E is:

- Full – A warning is displayed. Do the following:
 - Delete unneeded files, see “Maintaining Adequate Space for Database and Sample Data Storage” in the *Applied Biosystems® 3730/3730xl DNA Analyzer Maintenance and Troubleshooting Guide* (Part no. 4477797).
 - Click the green button to start the run.
- Not full – The run starts.

Note: A PostBatch Utility, which runs automatically, powers off the oven and the laser at end of a batch of runs.

Notes _____



DNA Sequencing Run Times

The following table lists the approximate run times of common DNA sequencing analysis runs:

Application	Capillary Array Length (cm)	Run Module	Approximate Run Time [†] (min)
Short read DNA Sequencing	36	TargetSeq36_POP-7™	20‡
Rapid read DNA sequencing	36	RapidSeq36_POP-7™	35
Standard read DNA sequencing	36	StdSeq36_POP-7™	60
Fast DNA sequencing	50	FastSeq50_POP-7™	60
Long read DNA sequencing	50	LongSeq50_POP-7™	120
Extra Long DNA sequencing	50	XLRSeq50_POP-7™	180

† Times assume oven is at temperature

‡ Approximate time to run 400 bases. The run module can be customized to run 200-400 bases.

Fragment Analysis Run Times

The following table indicates the approximate run time of a common fragment analysis run:

Application	Capillary Array Length (cm)	Run Module	Approximate Run Time (min)
Fragment Analysis	36	GeneMapper36_POP-7™	32
Fragment Analysis	50	GeneMapper50_POP-7™	43
SNPlex® Genotyping	36	HTSNP36_POP7_V3	15
SNPlex® Genotyping	50	HTSNP50_POP-7™	25







Notes _____



Controlling the Run

You can use the toolbar at the top of the Data Collection Software window to control the run.




To ...	Click ...	Action
Start the run		Starts run(s).
Stop the current run		Stops the current run.
Stop after the current run		Finishes current run and then stops.
Skip to next run		Stops the current run and begins next scheduled run.
Pause after current run		Finishes current run and then waits for resume command to begin next scheduled run.
Resume after pause		Begin the next scheduled run after a pause.

Notes _____




Monitoring the Status of the Run

In the navigation pane of the Data Collection Software, select  (Instrument Status) to view the status of the instrument or the current run.

System Status must be 'Ready' before a run starts

Array and polymer information



The screenshot shows the 'Instrument Status' window with the following data:

- Status Overview:** Instrument ID: M27PT, Run ID: (blank), Plate ID: (blank), System Status: Ready
- Array and Polymer Information:** Array Serial Number: 56789, Number of Capillaries: 96, Array Length: 50 cm, Array Usage: 1
- Sensor States:** Main Door: Closed, Oven Door: Closed, Stacker Drawer: Closed, In Stack: Not Empty, Out Stack: Not Empty, Laser: IDLE, EP: ON, Oven: OFF, Cell Heater: OFF, Buffer Heater: OFF, Instrument: IDLE
- Sensor Values:**
 - EP Voltage: 25.0 KV, EP Current: 2000.0 uA
 - Laser Power: 45.0 mW, Laser Current: 15.0 A
 - Buffer Heater: 100.0 degC, Cell Heater: 100.0 degC
 - Oven Temp: 100.0 degC
- Tray States:** AutoSampler Plate Type: 96-well septa, Buffer Station Plate: Present, Water Station Plate: Present, Waste Station Plate: Present, Parking Station Plate: Empty
- Events:** Spectral calibration has completed, Batch ended, System Status: Postprocessing, Run_M27PT_2003-06-10_11-33_0001 status I, Sample Plate Unloaded, Stacker status has changed, Outstacker error state has changed, System Status: Post-Batch, The number of runs has changed to 0, Run completed, System Status: Idle, System Status: Ready
- Errors:** Expected Raman line was found at pixel 233, Can't start run, spatial calibration was n, Sealing type (from instr) is Heat Sealing, Container and numCap validation failed, Sealing type (from instr) is Heat Sealing, Container and numCap validation failed.

System Status changes from green to flashing red when errors occur.

Notes



Events Box Displays the:

- Recent actions of the instrument
- Status of each capillary (passed or failed) at the end of a spectral calibration
- Calibration data at the end of a spatial calibration

Some of the events listed in the Events box provide information for service engineers.

Errors Box Displays errors that have occurred during the current run

Some of the error messages provide information for service engineers. A “fatal” error usually requires that you restart the Data Collection Software.


Notes _____



Viewing Real-Time Electrophoresis Data


Use the EPT Viewer to view real-time electrophoresis (EP) data during a run.

To access the viewer, in the navigation pane of the Data Collection Software, select **GA Instruments** > **ga3730** > **instrument name** > **Instrument Status** > **EPT Chart**.



GA Instruments > ga3730 > Run History > EPT Viewer

Select a run to view: Run_Huahine_2003-10-04



Parameter	Color	Show
EpVoltage	Blue	<input checked="" type="checkbox"/>
EpCurrent	Blue	<input checked="" type="checkbox"/>
OvenTemp	Blue	<input checked="" type="checkbox"/>
LaserPower	Green	<input checked="" type="checkbox"/>
AmbientTemp	Green	<input checked="" type="checkbox"/>
CCDTemp	Orange	<input checked="" type="checkbox"/>
OvenPeltierTemp	Red	<input checked="" type="checkbox"/>
LaserCurrent	Red	<input checked="" type="checkbox"/>
BufferTemp	Black	<input checked="" type="checkbox"/>
ArrayHolderTemp	Black	<input checked="" type="checkbox"/>

Show Legend

Notes




Viewing Event History

Use the Event log window to view a record of operational events, as shown in the next figure.

To access the Event Log window, in the navigation pane of the Data Collection Software, click **GA Instruments** > **ga3730** > **instrument name** > **Instrument Status** > **Event Log**.

IMPORTANT! To delete error messages, select all error messages, then click **Clear Errors**. The system status light flashes red until all errors are cleared.

Note: Using the Event Log window, you can also verify the capillary-by-capillary processing status during a spectral calibration run.



Note: If an error is generated while using manual control, reboot the instrument then restart the Data Collection Software to recover from the error stage.


Notes



Viewing Electropherogram Data

Viewing Data in the Capillary Viewer

Use the Capillary Viewer to examine the quality of electropherogram data from multiple capillaries during a run. In the navigation pane of the Data Collection Software, select **GA Instruments > ga3730 > instrument name > Instrument Status > Capillary Viewer**.





Electropherogram Displays

An electropherogram is a graph of relative dye concentration as a function of time, plotted for each dye. The displayed data has been corrected for spectral overlap (multicomponented).

How to Zoom

To zoom an area of an electropherogram:

1. Click-drag the mouse over the area of interest.
2. Release the mouse, then click  to expand the view.
3. Click  to return to full view.

Click individual colors to view or hide them.




Notes




Viewing Data in the Array Viewer

Use the Array Viewer during or after a run to examine the quality of your data from all capillaries. You can view all the capillaries (vertical axis) as a function of time/data point (horizontal axis).

To open the Array Viewer window in the navigation pane of the Data Collection Software, select **GA Instruments** > **ga3730** > **instrument name** > **Array Viewer**.



How to Zoom

1. To expand the view, click-drag the mouse over the area of interest.
2. Click  to return to full view.

Displaying or Hiding Color



Click individual colors in the color bar to view or hide the color in the Array View (same in Capillary Viewer).

Notes _____



Viewing the Run History Data

Run History Components

To view the Run History utility can be used only with completed runs stored in the local 3730/3730xI Analyzer Data Collection database. It does not provide real-time viewing of collecting runs.

In the navigation pane, click the icon next to the function to launch it.

Run History Views	Icon
EPT Viewer Note: If Cleanup Database has been used, you cannot view processed data in Run History.	
Spatial Calibration Viewer	
Capillary Viewer Note: If Cleanup Database has been used, you cannot view processed data in Run History.	
Array Viewer Note: If Cleanup Database has been used, you cannot view processed data in Run History.	
Spectral Calibration Viewer	
Reextraction Note: If Cleanup Database has been used, you cannot view processed data in Run History.	


Viewing Data from a Completed Run

There are two formats for viewing data within the 3730/3730xI Analyzer Data Collection Software under the Run History icon:

- In the Array Viewer
- In the Capillary Viewer capillary-by-capillary

1. In the navigation pane of the 3730/3730xI Analyzer Data Collection software, select (**Run History**).

Notes _____



2. Search for the run you want to use by either Barcode or Advanced search.
3. After choosing the run, select the **Array Viewer** or the **Capillary Viewer** in the navigation pane.

Notes




Viewing the Results of Autoextraction

After a run is completed extraction and analysis are performed automatically, according to the settings in the Plate Editor and the Results group. The results of extraction and analysis can be viewed in the Reextraction Panel. Samples can be extracted again with the same settings, or with different Analysis Protocols or different Results Groups. This can be useful for several reasons:

- The destination location may not have been available during extraction.
- Some samples may have failed analysis and a different Analysis Protocol might be more successful.
- Samples might be saved in different locations, or with no analysis at all to save space.
- Sample files are created based on the your destination and folder naming selections.


Runs Stopped Before Complete Autoextraction

Runs that are stopped before completion display the “Completed” status in the Run Scheduler, and the associated plate is moved to the Out Stack. In the Instrument View the status is changed to “Ready”. Successfully extracted and analyzed runs display the “Processed” status in the Run Scheduler.

The auto extractor component of the 3730/3730xI Analyzer Data Collection automatically extracts data from stopped runs. If autoextraction fails, click Reextraction  to extract data.

Selecting and Queuing Samples for Reextraction


You can queue individual samples for reextraction. This is especially useful for experimenting with different analysis protocols for samples that have failed initial extraction.

1. Click  (Run History).
2. Enter the plate ID for a plate that has been run, then click **Search**. All completed runs from that plate appear in the window and can be reextracted. Pending runs from the plate do not appear in the window.
3. Select a run from the list.

Notes _____



Run Name	Plate ID	Plate Name	Type	Size	Operator	Last Modified
Run_Huahine_2002-10-18_04-09_3	DS33InstallPlate	DS33InstallPlate	GeneMapper	96-Well	maf	2002-10-23 22:49:10.0
Run_Huahine_2002-10-18_20-37_7	DS33InstallPlate	DS33InstallPlate	GeneMapper	96-Well	maf	2002-10-23 22:49:10.0
Run_Huahine_2002-10-18_20-37_8	DS33InstallPlate	DS33InstallPlate	GeneMapper	96-Well	maf	2002-10-23 22:49:10.0
Run_Huahine_2002-10-18_20-37_9	DS33InstallPlate	DS33InstallPlate	GeneMapper	96-Well	maf	2002-10-23 22:49:10.0
Run_Huahine_2002-10-18_20-37_10	DS33InstallPlate	DS33InstallPlate	GeneMapper	96-Well	maf	2002-10-23 22:49:10.0
Run_Huahine_2002-10-23-03_1	DS33	DS33Install	GeneMapper	96-Well	install	2002-10-23 22:39:37.0
Run_Huahine_2002-10-24_02-32_2	JaimeTest	Jaime	GeneMapper	96-Well	Jaime	2002-10-24 02:29:28.0
Run_Huahine_2002-10-25_02-08_2	Verification_Plate	Verification_Plate	SequencingAnalysis	96-Well	3730User	2002-10-25 02:06:38.0
Run_Huahine_2002-10-25_04-50_3	LRSPlate	LRSPlate	SequencingAnalysis	96-Well	KK	2002-10-25 04:49:47.0

4. Click  (Reextraction) in the navigation pane. The Reextraction window opens.
5. Select the checkboxes in the Extract column that correspond to the samples to be reextracted.
6. Click **Extract** to start the reextraction.

Note: Reextracted sample files are saved in the original folder that data was extracted to, unless you modify the results group settings.

Notes



Reextraction Window for Sequencing Analysis

Click the boxes to select samples to be reextracted

Select a run

Extraction Result column on the Reextraction window

Extract	Cap	Well	Extraction Result	Results Group	Analysis Protocol	Analysis Result	Score	Sample Name	Extraction Comment
<input checked="" type="checkbox"/>	8	A01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	5.0	s	
<input checked="" type="checkbox"/>	7	B01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	5.0	s	
<input checked="" type="checkbox"/>	6	C01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	48.510754	s	
<input checked="" type="checkbox"/>	5	D01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	4	E01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	43.687626	s	
<input checked="" type="checkbox"/>	3	F01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	1.0	s	
<input checked="" type="checkbox"/>	2	G01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	22.686636	s	
<input checked="" type="checkbox"/>	1	H01	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	6.0	s	
<input checked="" type="checkbox"/>	16	A03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	1.0	s	
<input checked="" type="checkbox"/>	15	B03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	6.0	s	
<input checked="" type="checkbox"/>	14	C03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	1.0	s	
<input checked="" type="checkbox"/>	13	D03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	12	E03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	43.992992	s	
<input checked="" type="checkbox"/>	11	F03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	31.95	s	
<input checked="" type="checkbox"/>	10	G03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	44.98598	s	
<input checked="" type="checkbox"/>	9	H03	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	24	A05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	49.410084	s	
<input checked="" type="checkbox"/>	23	B05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	22	C05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	28.083334	s	
<input checked="" type="checkbox"/>	21	D05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	23.813953	s	
<input checked="" type="checkbox"/>	20	E05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	29.166666	s	
<input checked="" type="checkbox"/>	19	F05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	18	G05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	49.475758	s	
<input checked="" type="checkbox"/>	17	H05	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	2.0	s	
<input checked="" type="checkbox"/>	32	A07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	1.0	s	
<input checked="" type="checkbox"/>	31	B07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	1.0	s	
<input checked="" type="checkbox"/>	30	C07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	28.904762	s	
<input checked="" type="checkbox"/>	29	D07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	ERROR: Analysis Failed	<NA>	s	
<input checked="" type="checkbox"/>	28	E07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	27	F07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	ERROR: Analysis Failed	<NA>	s	
<input checked="" type="checkbox"/>	26	G07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	33.59501	s	
<input checked="" type="checkbox"/>	25	H07	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	46.06404	s	
<input checked="" type="checkbox"/>	40	A09	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	48.61644	s	
<input checked="" type="checkbox"/>	39	B09	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	8.0	s	
<input checked="" type="checkbox"/>	38	C09	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	4.0	s	
<input checked="" type="checkbox"/>	37	D09	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	3.0	s	
<input checked="" type="checkbox"/>	36	E09	SUCCESS: Ext...	SeqA	3730BDTV3-KB-DeNovo...	SUCCESS: Analysis Succ...	16.0	s	

Click here to start extraction

Use these if several samples are highlighted

Notes



Reextraction Window for Fragment Analysis

Click the check boxes to select samples to be reextracted

Select a run — Extraction Result column of the Reextraction window

Extract	Cap	Well	Extraction Result	Results Group	Sample Name	Comment	Sample Type	Size Standard	Plate
<input checked="" type="checkbox"/>	1	A01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	3	B01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	5	C01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	7	D01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	9	E01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	11	F01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	13	G01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	15	H01	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	2	A02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	4	B02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	6	C02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	8	D02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	10	E02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	12	F02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	14	G02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:
<input checked="" type="checkbox"/>	16	H02	SUCCESS: Extr	gm_runbyrun	s		Sample	GSS00LIZ	D:

Click here to start extraction — Use these if several samples are highlighted

Notes _____




Results Column of the Reextraction Window

The results of extraction and analysis are color coded in the Results column of the Reextraction window. The following table indicates the colors and their values.

Color	Value	Notes
Red	Extraction or analysis failed	Descriptive messages can be viewed by resizing the Results column to view all text (click on the arrow)
Yellow *	Warnings for extraction or analysis	
Green	Successful extraction (with no analysis intended), or successful extraction and analysis.	

* Note: The text message for samples that produce yellow is: "FAILURE: Analysis Fail Bad Data; Error Number=nnnnn WARNING..."

The Results column, by default, shows only the beginning of any processing message. The entire message returned from extraction and autoanalysis can be viewed by expanding the cell.



Quality Column of the Reextraction Window

The Quality column represents the quality values for an entire sequence. Quality values are assigned only to analyzed samples when using the KB™ Basecaller. The Quality column is empty (white) if:

- Analysis was not performed
- Analysis failed
- ABI Basecaller was used for analysis. ABI basecaller does not assign quality values

Notes



Results Group and Analysis Protocol Columns

The Results Group and the Analysis Protocol (Analysis Method in the GeneMapper® software) can be edited and the changes used for reextraction.

Note: Select an entire column in the Reextraction window by clicking the column header. For example, clicking the Extract column header selects all samples. Clicking the Uncheck or Check buttons at the bottom of the window, enables or disables the checkboxes for each sample. Additionally, the fill-down command (Ctrl+D) works the same here as in the Plate Editor for easier information input.

Sorting The Samples

The samples can be sorted according to any of the column properties by holding down the Shift key while clicking on the column header. Shift-clicking a column a second time sorts the column contents in the reverse order. This is most useful for sorting by capillary number, by well position, by results, by quality, and by the Extract column. For example, it is often useful to bring all the samples that failed analysis or extraction to the top of the column where they can be examined without having to scroll down to each sample individually.

Reextracting Selected Samples

1. Expand the Results column cells for any yellow or red results, to see a description of the warning or failure.
2. Select a new Results Group, or edit the current one. This allows you to turn off autoanalysis, change the samples and folder naming options, the location where they are placed, the owner of the Results Group, and so on.
3. If desired, change the analysis protocol to experiment with different ways of analyzing the sample, using a different basecaller for example.
4. Select the check box in the Extract column for the samples you wish to extract again.
5. Click **Extract**.

IMPORTANT! Reextraction creates a new sample file and does not replace the previously saved sample file. The presence of a previous sample file has no effect on the creation of a new sample file. If the naming options that are used for reextraction are identical to those used previously, a number is added to the filename. For example, if the first sample is, “sample01.ab1” then the second sample would be, “sample01.2.ab1”.

Notes _____



Chapter 6 Running the Instrument

Viewing the Results of Autoextraction

Notes _____

Catalog List

Item	Cat. and Part no.
3730 36-cm capillary array	4331247
3730 50-cm capillary array	4331250
3730xl 36-cm capillary array	4331244
3730xl 50-cm capillary array	4331246
3700/3730 BigDye Terminator v3.1 Sequencing Std	4336943
3700/3730 BigDye Terminator v1.1 Sequencing Std	4336799
Matrix Standard Set DS-33	4345833
HiDi™ Formamide, 25 mL	4311320
POP-7™ Polymer (1 bottle of 25ml each)	4363929
POP-7™ Polymer (10 bottles of 25ml each)	4363935
POP-7™ Polymer (30 bottles of 25ml each)	4335611
POP-6™ Polymer (1 bottle of 7ml each)†	4352757
POP-6™ Polymer (1 bottle of 3.5ml each)†	4363783
Buffer (10X) with EDTA - 500 mL	4335613
Buffer (10X) with EDTA - 4L	4318976
96-Well sample plates w/barcode	4306737
96-Well sample plates, no bar code	N801-0560
96-Well plate septa	4315933
96-Well plate base (septa sealed)	4334873
96-Well plate base (heat sealed)	4334875
96-Well plate retainer (septa sealed)	4334869
96-Well and 384-well Plate Retainer (heat sealed)	4334865

Notes _____

Item	Cat. and Part no.
FAST (0.1ml) 96-Well Plate Retainer for 3730 (septa-sealed)	4367472
FAST (0.1ml) 96-Well Plate Base for 3730 (septa-sealed)	4367469
FAST (0.1ml) 96-Well Plate Retainer for 3730 (heat-sealed)	4367474
FAST (0.1ml) 96-Well Plate Base for 3730 (heat-sealed)	4367473
384-Well Sample plates with barcode	4309849
384-Well plate septa	4315934
384-Well plate base (septa-sealed)	4334874
384-Well plate base (heat-sealed)	4334877
384-Well plate retainer (septa-sealed)	4334868
Heat seal film, 3-mil	4337570
Applied Biosystems® 3730/3730xI DNA Analyzer Getting Started Guide	4359476
Applied Biosystems® BigDye Xterminator® Purification Kit Protocol	4374408
AB Navigator Software Administrator Guide	4477853

† Call Technical Support for an adaptor to use POP-6 on the 3730 Series.

Notes _____

Dye Sets: G5, G5-RCT, Any4Dye, Any4Dye-HDR, and Any5Dye

Supported Dye Sets

Sequencing Analysis Dye Sets for All Applications

Dye Set	Application Name
E_BigDyeV1	DNA sequencing with BigDye® Terminator v1.1 Cycle Sequencing Kit
Z_BigDyeV3	DNA sequencing with BigDye® Terminator v3.1 Cycle Sequencing Kit
Z_BigDyeV3	DNA sequencing with BigDye® Direct Cycle Sequencing Kit, with combined DNA PCR Amplification/Clean-up/Cycle Sequencing kit

Fragment Analysis Dye Sets for All Applications

Dye Set	Application Name
G5	DNA sizing for 5-dye chemistry
G5-RCT	DNA sizing for 5-dye chemistry
Any5Dye	SNaPshot® Multiplex System

Additional Dye Sets

Dye Set	Application Name
Any4Dye-HDR	DNA sizing and DNA sequencing
Any4Dye	DNA sizing and DNA sequencing

Notes _____

Dye Sets G5 and G5-RCT For Fragment Analysis

Overview Even small levels of crosstalk could be a concern for users of the 3730/3730xl instruments who perform fragment analysis as well as for applications with a high dynamic range. In fragment analysis applications that have few sample peaks and varying peak intensities, a crosstalk peak may appear as a real sample peak and be incorrectly identified as an allele. Crosstalk is not a concern with sequencing applications as there is a constant stream of peaks electrophoresing past the detector.

Dye Set G5-RCT To reduce crosstalk for fragment analysis applications, a new dye set has been created for Data Collection Software v3.0, called dye set G5-RCT. G5-RCT uses the same chemistry as dye set G5 (6-FAM™, VIC® NED™, PET®, LIZ® dyes). This dye set reduces signal, but reduces potential crosstalk to a greater degree, so the reduction in signal-to-noise ratio is less pronounced than the reduction in signal overall. Higher concentration peaks can be used without going offscale, this results in a higher dynamic range for the G5-RCT dye set.

Recommendations for Using G5 or G5-RCT Dye set G5-RCT may be especially useful for users performing fragment analysis with a 96 capillary array, as well as users interested in applications with a high dynamic range (large peaks much higher than small peaks). For most other conditions, users prefer the G5 dye set.

Life Technologies supports:

- Fragment analysis on the 96-capillary array using G5-RCT only
- G5 and G5-RCT on the 48-capillary array.

Notes _____

Refer to the following table for more information about the advantages and issues to consider for each dye set.

Dye Set	Features
Standard Z, E Dye Sets	<p>When to use/Advantages:</p> <ul style="list-style-type: none"> • All DNA sequencing applications using BigDye® Terminators v3.1 and v1.1 and BigDye® Direct. • Higher signal relative to the Any4Dye-HDR dye set • Optimized for the highest signal-to-noise ratio <p>Issues:</p> <ul style="list-style-type: none"> • More susceptible to samples within a plate with large variation in peak height relative to the Any4Dye-HDR dye set
Any4Dye	<p>When to use/Advantages:</p> <ul style="list-style-type: none"> • Use of unsupported dyes. Provides an open platform for system capable applications <p>Issues:</p> <ul style="list-style-type: none"> • Performance of system has not been tested nor can the performance be guaranteed • More susceptible to samples within a plate with large variation in peak height relative to the Any4Dye-HDR dye set
Any4Dye-HDR (High Dynamic Range)	<p>When to use/Advantages:</p> <ul style="list-style-type: none"> • High dynamic range when samples within a plate have a large variation in peak height • Resequencing/Mutational Profiling applications • 4-Dye Fragment Analysis applications • Use of unsupported dyes. Provides an open platform for system capable applications <p>Issues:</p> <ul style="list-style-type: none"> • Signal intensity is reduced by approximately half relative to the standard dye sets, along with a minimal reduction in the noise, resulting in a very slight decrease in the signal/noise ratio when compared to data generated using the standard dye sets • Essential that spectral calibrations are performed each time the capillary array is replaced or moved within the detection cell


Notes _____

Creating a Spectral Calibration for Dye Sets Any4Dye, Any4Dye-HDR, or Any5Dye


The steps to creating and running a customized 4- or 5-Dye Set are similar to running a supported dye set.

The following example illustrates the use of Any4Dye dye set; it works the same for Any5Dye dye set.

1. In the navigation pane of the Data Collection Software, click **GA Instruments** > **ga3730** > **Protocol Manager**.
2. In the Instrument Protocols pane, click **New...**. The Protocol Editor opens.
3. In the Protocol Editor, create a spectral protocol for the 4Dye dye set, specifying the appropriate protocol parameters.
4. Click **OK** to save the spectral protocol.



Note: Customize the Spectral parameters as needed. For more information, see [step 1 on page 49](#).



Notes

5. Click **New** in the Plate Manager to display the New Plate Dialog box.
6. Create a spectral plate for the Any4Dye dye set by completing the New Plate Dialog box.
7. Click **OK**.
8. Create an instrument protocol. For more information, see [page 47](#).
9. In the Plate Editor, select the Instrument Protocol that you just created in the previous steps, then click **OK** to save the plate.

New Plate Dialog

ID (Barcode): Any4Dye_Spectral
Name: Any4Dye_Spectral
Description:

Application: Spectral Calibration
Plate Type: 96-Well
Scheduling: 1234
Plate Sealing: Septa
Owner Name: sc
Operator Name: sc

OK Cancel

Spectral Calibration Plate Editor

File Edit

Plate Name: Any4Dye_Spectral Operator: sc
Plate ID: Any4Dye_Spectral Owner: sc
Plate Sealing: Septa


Well	Sample Name	Comment	Instrument Protocol 1
A01	s		4Dye_Spectral
B01	s		4Dye_Spectral
C01	s		4Dye_Spectral
D01	s		4Dye_Spectral
E01	s		4Dye_Spectral
F01	s		4Dye_Spectral
G01	s		4Dye_Spectral
H01	s		4Dye_Spectral
A02	s		4Dye_Spectral
B02	s		4Dye_Spectral
C02	s		4Dye_Spectral
D02	s		4Dye_Spectral
E02	s		4Dye_Spectral
F02	s		4Dye_Spectral
G02	s		4Dye_Spectral
H02	s		4Dye_Spectral
A03	s		4Dye_Spectral
B03	s		4Dye_Spectral
C03	s		4Dye_Spectral
D03	s		4Dye_Spectral
E03	s		4Dye_Spectral
F03	s		4Dye_Spectral
G03	s		4Dye_Spectral
H03	s		4Dye_Spectral

Description:

OK Cancel


Notes

10. In the Run Scheduler, add the spectral plate to the Input Stack, then run the plate.



Notes _____

11. Verify that spectral matrices for all capillaries meet acceptance criteria (pass). Override individual capillaries and rename calibration as needed.




Notes

Regular Runs Using Any4Dye or Any5Dye Dye Sets

The following example shows the use of Any4Dye dye set. This process works the same for Any5Dye set.


1. In the Protocol Editor, create a regular instrument run protocol for the Any4Dye dye set, then choose the appropriate default run module template. (You can create a customized run module in the Module Editor if desired).



2. In the Plate Manager, create a regular plate, selecting the Any4Dye instrument protocol you created in step 1.


Notes _____

- In the Plate Editor, select the instrument protocol that you created in step 1, then click **OK** to save the plate.



Notes _____

4. In the Run Scheduler, add this plate to the Input Stack, then run the plate.



Notes

KB™ Basecaller Software v1.4.1

April 2012

Contents:

Executive Summary	175
Benefits of using the KB™ Basecaller	176
Future support of ABI and KB™ Basecallers	179
Features in KB™ Basecaller Software v1.4.1	179
Comparison of the ABI and KB™ Basecallers	180
Differences between the ABI and KB™ Basecallers.....	181
FAQs: Processing data with Phred software and .phd.1 Files	184
FAQs: Quality values.....	185
Miscellaneous FAQs	187
Conference posters and reference	188

Executive Summary

Applied Biosystems® KB™ Basecaller Software v1.4.1 reduces manual data review time and increases the read length of high-quality bases in sequences. This algorithm accurately extracts more bases out of the sequencing data generated on Applied Biosystems® DNA Analyzers and Genetic Analyzer Instrument and chemistry platforms. KB™ Basecaller Software v1.4.1 supports all BigDye® Terminator v3.1 and v1.1 and BigDye® Direct chemistries and run modules available on Applied Biosystems® instruments.

- 310 Genetic Analyzer
- 3100/3100-*Avant* Genetic Analyzers
- 3130/3130*xl* Genetic Analyzers
- 3730/3730*xl* DNA Analyzers
- 3500 Dx and 3500 Dx/3500xL Dx Genetic Analyzers
- 3500 and 3500/3500xL Genetic Analyzers.

Notes _____

Software integration

KB™ Basecaller Software v1.4.1 is integrated with:

- Sequencing Analysis Software 6 and v5.4
- SeqScape® Software 3 and v2.7
- Variant Reporter™ Software 2 and v1.1
- 3130 Series and 3730 Series Data Collection Software 4
- 3500 Series Data Collection Software
- MicroSEQ® ID Analysis Software v2.2

KB™ Basecaller Software v1.4.1 is *not* integrated with:

- MicroSeq® ID software versions 2.1 and older
- Any versions of Data Collection Software for the 310 and 3100/3100-*Avant*
- 3130/3130*xl* and 3730/3730*xl* Data Collection Software versions before v3.1
- Sequencing Analysis Software before v5.4
- SeqScape® Software versions before v2.7
- Variant Reporter™ Software versions before v1.1.

During the co-installation of Sequencing Analysis Software 6 and SeqScape® Software 3 with Data Collection Software 4, KB™ Basecaller Software v1.4.1 is installed into your Data Collection Software 4 on the same computer.

Testing on more than 50,000 sequencing samples shows that version 1.4.1 of the algorithm offers many advantages, including longer, accurate read lengths.

Details of the test and validation process are in the poster *Longer Reads and More Robust Assemblies with the KB™ Basecaller*.

IMPORTANT! Life Technologies strongly recommends using the KB™ Basecaller.

Benefits of using the KB™ Basecaller

Some benefits of using the KB™ Basecaller include:

- Increased length of read
- Per-base quality value predictions using an equation that is standardized by Phred software
- Optional detection of mixed-base with quality values
- Analysis of short PCR products
- Accurate start point detection
- Increased accuracy in regions of low signal-to-noise or anomalous signal artifacts
- Detection of failed samples
- Trimming of data using per-base quality value
- Per-sample quality value that helps to determine the quality of each read

Notes _____

- Optional detection of PCR stop
- Optional assignment of Ns
- Optional generation of .phd.1 files

Increased length of read

KB™ Basecaller accurately extracts more bases than ABI Basecaller from the 3' and 5' ends of a sequence. Tests on genomic BAC samples, performed on data generated using 3730/3730xl instruments, indicate an improvement of approximately 100 bases in length-of-read as compared to the same data analyzed by the ABI Basecaller and Phred software (v0.020425.c). The gain in read length varies depending on the run module used to collect the data. The accuracy of start point estimation and the first 50 bases of called sequence is substantially increased. Typically, ~10 more correct calls on average are identified at the 5' end as compared to the ABI Basecaller.

Per-base quality value predictions

The KB™ Basecaller assigns quality values to every basecall. The quality prediction algorithm is calibrated to return Q values that conform to the industry-standard relation established by the Phred software. The KB™ Basecaller and its output are, therefore, interchangeable in processes requiring Phred software for output.

Quality value calibration was performed using a set of correct-sequence annotated sample files, representative of production sequencing data generated on capillary electrophoresis platforms. Over 52.1 million basecalls were used to calibrate KB™ Basecaller Quality Values and over 32.9 million distinct basecalls were used to test the calibration.

Accuracy in start point detection

Improved start point detection contributes to better mobility shift corrections and greater basecalling accuracy in the first 50 bases. Because the KB™ Basecaller detects the start point accurately, you do not need to manually set start points for each sample.

Optional detection of mixed-base with quality values

The KB™ Basecaller can detect mixed base positions, and assign two-base (R, Y, K, M, S, W) IUB codes and quality values to those positions. Quality values are assigned to mixed basecalls using an algorithm similar to that for pure bases.

The definition conforms to the Phred relation. Quality values for mixed bases are inherently lower than those of pure bases due to the higher error risk of interpreting more complex signals. Note that when using the ABI Basecaller or ABI Basecaller and Phred software, a separate analysis stage is required to determine mixed bases.

Increased accuracy in regions of low signal-to-noise or anomalous signal artifacts

The KB™ Basecaller increases the accuracy of sequence reads from low-signal regions or from data that are partially contaminated by a secondary sequence or by other sources of “chemistry noise”.

Basecalling errors caused by anomalous chemistry and/or instrument signals such as dye blobs and fluorescent spikes are substantially reduced. These artifacts often occur in otherwise high-quality “clear-range” data. They result in the loss of high-quality bases that are downstream from the noise region. Tests indicate that KB™ Basecaller distinguishes between target DNA peaks and the most common artifacts better than ABI Basecaller.

Notes

Analysis of short PCR products	The KB™ Basecaller has been tested for accuracy in basecalling and quality value estimates on PCR products as short as 100 bases. Although KB™ Basecaller may be able to basecall products with less than 100 bases, these types of sample files were not tested.
Detection of failed samples	The KB™ Basecaller indicates the gross sample quality of each analysis as “Success without warnings,” “Success with warnings,” or “Failure due to poor data quality”. A common failure mode is no signal—insufficient detection of DNA peaks. For failed samples, the KB™ Basecaller uses “NNNNN” as the sequence, indicating that the sample quality is very low and may need to be omitted from further analysis. Failed samples are flagged in reports in the analysis software. Note that this behavior is different from the ABI Basecaller, which <i>always</i> tries to call bases, resulting in sequences of many Ns.
Option to trim data using per-base quality value	You can use software with KB™ Basecaller to automatically determine the clear range region by trimming the ends using the per-base quality values. The parameters used for trimming are similar to those in other tools used by the genome community.
Per-sample quality value (QV) evaluates quality of reads	Software with the KB™ Basecaller uses the QV from the KB™ Basecaller to trim and determine a sample score. The sample score is the average QV in the clear range, or, if no clear range is determined, in the entire read. This single number value is a measure of the quality of the data. The sample score appears in reports generated by Sequencing Analysis Software, SeqScape® Software, Sequence Scanner Software, Variant Reporter™ Software, and/or MicroSeq® ID Software.
Optional detection of PCR stop	You can set the KB™ Basecaller to end basecalling at a PCR stop. Note that samples with enzymatic failure may have signal properties similar to those in PCR stop conditions. The KB Basecaller may not be able to distinguish between these two conditions.
Optional assignment of Ns	By default, the KB™ Basecaller does not generate Ns. However, you may choose to reassign Ns to bases with QVs below a user-specified threshold for both pure and mixed base positions.
Optional generation of .phd.1 files	.phd.1 files can be generated by autoanalysis or in analysis software. You can use the .phd.1 files for further analysis by downstream software such as Phred software.

Notes

Future support of ABI and KB™ Basecallers

Life Technologies will continue to provide technical support for the ABI Basecaller. However, further development and defect fixes will occur only on the KB™ Basecaller. If you encounter a defect in the ABI Basecaller, please use the KB™ Basecaller instead. In future releases, ABI Basecaller support files are removed from the software wherever they duplicate support in the KB™ Basecaller.

Features in KB™ Basecaller Software v1.4.1

- A basecalling algorithm that supports Applied Biosystems® 310, 3100/3100-*Avant*, 3130/3130*xl*, 3730/3730*xl*, 3500/3500xL, and 3500 Dx/3500xL Dx Genetic Analyzers
- Improvements over all earlier versions of KB™ Basecaller (v1.0, v1.1, v1.1.1, v1.1.2, v1.2, v1.3, and v1.4)

Note: Basecalling results with KB™ Basecaller Software v1.4.1 may differ slightly from results obtained with previous versions of KB™ Basecaller.

Notes _____

Comparison of the ABI and KB™ Basecallers

Question	ABI Basecaller	KB™ Basecaller
What does the software do?	<ul style="list-style-type: none"> Processes raw traces Provides processed traces Provides AGCTN calls 	<ul style="list-style-type: none"> Processes raw traces Provides processed traces Provides pure bases only <i>or</i> Provides pure and mixed calls Provides quality values Generates .phd.1 and .scf files Provides a sample score
What are the resulting basecalls?	<p>One option available: Only mixed bases are assigned as Ns. Further processing (either manual or using additional software) is required to assign IUB codes to the Ns or pure bases.</p>	<p>Four options are available. The software can assign an:</p> <ul style="list-style-type: none"> ACGT and Q value to each peak. ACGT and Q value to each peak. Any peak with a Q value below a defined threshold is reassigned an N. ACGT or a mixed base and a Q value to each peak. ACGT or a mixed base and a Q value to each peak. Any peak with a Q value below a defined threshold is reassigned an N.
How are failed samples handled (for example, no signals, chemistry failure)?	Attempts to call all bases so a sample results with many Ns.	Assigns five Ns to the entire sample to indicate that the sample failed analysis. The analysis report flags these files.
Baseline in processed data	Appears smoother than in KB™ Basecaller.	Appears less smooth than in ABI KB™ Basecaller.
What are the steps to process data?	Calls bases on Windows OS.	Calls bases and estimates QVs on Windows OS.
Data and future support	Supports the 310, 3100, 3100- <i>Avant</i> , Applied Biosystems® 3130/3130 <i>xl</i> and 3730/3730 <i>xl</i> instruments. Further development has stopped.	Applied Biosystems® 310, 3100/3100- <i>Avant</i> , 3130/3130 <i>xl</i> , 3730/3730 <i>xl</i> , 3500/3500 <i>xL</i> , and 3500 Dx/3500 <i>xL</i> Dx Genetic Analyzers. Development is ongoing.

Notes _____

Differences between the ABI and KB™ Basecallers

Question	Answer	
	ABI Basecaller	KB™ Basecaller
Can the KB™ Basecaller basecall short PCR products?		The KB™ Basecaller has been tested for accuracy in basecalling and quality value estimation on PCR products as short as 100 bases. Although it may be possible to basecall products with less than 100 bases, such sample files have not been tested. Samples shorter than 100 bases may not contain enough signal information to basecall the sample file.
Why is the baseline less smooth when the data are analyzed with the KB™ Basecaller?	<p>Processed signals or traces from the ABI Basecaller appear smoother than those from the KB™ Basecaller because each software application uses an algorithm that processes the signals differently.</p> <p>The ABI Basecaller assigns only AGCT and Ns to each peak. Therefore, you must manually search for mixed bases or use a secondary software to complete the task. To facilitate this secondary process, the ABI Basecaller subtracts an aggressive baseline estimate to show a cleaner baseline in the processed signals.</p>	The KB™ Basecaller can determine pure and mixed bases. Therefore, second-stage processing, which allows less aggressive baseline subtraction, is not needed. The processed traces have a higher baseline. If you have mixed bases, turn on the mixed-base detection option and allow KB™ Basecaller to call mixed bases. Use the mixed base calls and the associated QVs to review mixed bases – do not look only at the baseline.
What is the signal to noise value found with data analyzed with the KB™ Basecaller?	<p>The signal-to-noise value is the average of the signal intensity of the A, C, G, or T base divided by the average of the noise for that base.</p> <p>The ABI Basecaller calculates only the signal intensity. The signal-to-noise value is more indicative of data quality than the signal intensity value alone. Both properties are important in determining quality.</p>	KB™ Basecaller calculates the information and presents the data in the Annotation view and analysis report.

Notes _____

Question	Answer	
	ABI Basecaller	KB™ Basecaller
What scaling options are available with the KB™ Basecaller?	The ABI Basecaller uses a scaling method closer to the “True profile” option than the “Flat profile” option.	<p>The KB™ Basecaller can display scaled data in two ways:</p> <ul style="list-style-type: none"> • True profile scaling With this method, the processed traces are scaled uniformly so that the average height of peaks in the region of strongest signal is about equal to a fixed value (for example, 1000). The profile of the processed traces is very similar to that of the raw traces. • Flat profile scaling The processed traces are scaled semi-locally so that the average height of peaks in any region is about equal to a fixed value (for example, 1000). The profile of the processed traces is flat on an intermediate scale (> about 40 bases). <p>You must decide which option is better suited to your circumstances. The sequence and QVs called by the KB™ Basecaller are independent of the selected scaling option.</p>
Does the KB™ Basecaller produce more usable sample files than the ABI Basecaller?		<p>Tests show that medium- and high-quality data result in more usable bases (longer read length) when analyzed by the KB™ Basecaller than by the ABI Basecaller.</p> <p>For very poor-quality data (samples with no, low, or noisy signal), the KB™ Basecaller does not provide more bases but instead fails the samples. By calling a string of “NNNNN” for the failed samples (instead of a sequence containing low QVs), the KB™ Basecaller indicates that the sample is unusable.</p>
Can the KB™ Basecaller analyze data generated on the ABI PRISM® 373, 377, or 3700 instruments?		<p>No, the KB™ Basecaller is calibrated to basecall and estimate the basecall quality for BigDye® Terminator chemistries on 310, 3100, 3100-Avant, and 3130/3130xI Genetic Analyzers, 3730/3730xI DNA Analyzers, and 3500/3500xL and 3500/3500 Dx/3500xL Dx Genetic Analyzers. Life Technologies has stopped support for the 373, 377, and 3700 instruments and data analysis.</p>

Notes _____

Question	Answer	
	ABI Basecaller	KB™ Basecaller
How can I determine which basecaller was used to analyze each sample file?		The Annotation view for each sample file and for the print header displays the basecaller name and version number. When displaying samples files, files analyzed by the KB™ Basecaller have QV value bars displayed above the electropherogram.
Are there any known incompatibilities when a sample file is analyzed with the KB™ Basecaller?		Life Technologies does not know of any incompatibility issues when a sample file (.ab1) is analyzed with the KB™ Basecaller and used in third-party software.

Notes _____

FAQs: Processing data with Phred software and .phd.1 Files

Question	Answer
<p>Can I analyze sample files with the KB™ Basecaller and then reprocess them with Phred software?</p>	<p>In principle, yes, but this is not recommended. The resulting quality values from Phred software are not calibrated—i.e., it is possible that Phred will over or under-predict quality in certain circumstances because it has not been trained on the type of processed electropherogram produced by the KB™ Basecaller. (Phred has been trained using the ABI Basecaller to produce the processed traces.)</p> <p>In addition, Phred replaces (and ignores) the initial called sequence. Reprocessing KB-analyzed samples with Phred, on average, degrades the accuracy of the analysis in terms of actual sequence error. Analysis improvements in KB™ Basecaller outlined above are lost.</p> <p>Studies by Life Technologies indicate that running Phred software on sample files processed by the KB™ Basecaller degrades the quality of the results.</p> <p>Analysis with KB™ Basecaller can generate .phd.1 files, which are interchangeable with any processes that currently depend on Phred.</p>
<p>Which Applied Biosystems® software generates .phd.1 files?</p>	<p>The following software products have KB™ Basecaller (version varies for each software) integrated and can generate .phd.1 files:</p> <ul style="list-style-type: none"> • ABI PRISM® 3100-<i>Avant</i> Data Collection Software v2.0 • ABI PRISM® 3100 Data Collection Software v2.0 • Applied Biosystems® 3130/<i>xl</i> and 3730/<i>xl</i> Data Collection Software v3.0 and later • Sequencing Analysis Software v5.2 and later • SeqScape® Software v2.5 and later • MicroSeq® ID Software v1.0 and later • Variant Reporter™ Software v1.0 and later

Notes _____

FAQs: Quality values

Question	Answer
How do I use quality values to review data?	<p>When analyzing data with pure bases, Life Technologies Corporation recommends that you use the following settings:</p> <p>Pure bases – Low QV = <15, Medium QV= 15–19, High QV= 20+ (default)</p> <p>When reviewing data with pure bases, use the QVs to briefly review bases with high QV(>20). Pay close attention to bases with medium QVs because you may need to make edits. Quickly review low-QV bases, although you will likely discard these bases from further analysis.</p> <p>Mixed base quality values will be lower than pure bases. For mixed bases, review all mixed basecalls. You may want to accept basecalls with quality values as low as 1.</p> <p>Mixed bases – Low QV = <5, Medium QV = 5–10 (investigate to determine the best range for your application)</p> <p>In all cases, keep in mind that, by definition, the predicted probability of error for a particular basecall is $10^{-q/10}$.</p>
What are the differences between quality values of mixed bases and pure bases?	<p>Pure bases and mixed bases have the same probability of error for the associated basecall ($10^{-q/10}$). Note the following:</p> <ul style="list-style-type: none"> • High-quality pure bases typically have QVs of 20 or higher. • The distribution of quality values for mixed bases differs dramatically from that of pure bases. • For mixed bases, quality values greater than 20 are rare. • Accurate mixed basecalls may be assigned quality values as low as 1, because the probability of error with mixed bases is higher. Review all mixed basecalls.
Can I trim my data using quality values?	<p>Yes. When using Data Collection, you can set trimming using QVs in the analysis protocols.</p> <p>When using Sequencing Analysis Software, SeqScape® Software, MicroSeq® ID Software or Variant Reporter™ Software, you can set trimming using QVs in the Analysis settings.</p>

Notes _____

Question	Answer																																				
<p>Is there a table that shows each quality value and its corresponding probability of error?</p>	<p>The following table shows each quality value and its corresponding probability of error. For a more extensive table, look in the Help menu or the Sequencing Analysis Software or the SeqScape® Software user guides.</p> <table border="1" data-bbox="634 422 1373 814"> <thead> <tr> <th>QV</th> <th>Pe</th> <th>QV</th> <th>Pe</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>79.0%</td> <td>35</td> <td>0.032%</td> </tr> <tr> <td>5</td> <td>32.0%</td> <td>40</td> <td>0.010%</td> </tr> <tr> <td>10</td> <td>10.0%</td> <td>41</td> <td>0.0079%</td> </tr> <tr> <td>15</td> <td>3.2%</td> <td>45</td> <td>0.0032%</td> </tr> <tr> <td>20</td> <td>1.0%</td> <td>50</td> <td>0.0010%</td> </tr> <tr> <td>21</td> <td>0.79%</td> <td>60</td> <td>0.00010%</td> </tr> <tr> <td>25</td> <td>0.32%</td> <td>99</td> <td>0.0000000013%</td> </tr> <tr> <td>30</td> <td>0.10%</td> <td></td> <td>–</td> </tr> </tbody> </table>	QV	Pe	QV	Pe	1	79.0%	35	0.032%	5	32.0%	40	0.010%	10	10.0%	41	0.0079%	15	3.2%	45	0.0032%	20	1.0%	50	0.0010%	21	0.79%	60	0.00010%	25	0.32%	99	0.0000000013%	30	0.10%		–
QV	Pe	QV	Pe																																		
1	79.0%	35	0.032%																																		
5	32.0%	40	0.010%																																		
10	10.0%	41	0.0079%																																		
15	3.2%	45	0.0032%																																		
20	1.0%	50	0.0010%																																		
21	0.79%	60	0.00010%																																		
25	0.32%	99	0.0000000013%																																		
30	0.10%		–																																		
<p>Where can I see quality value bars and numbers?</p>	<p>Sequencing Analysis Software, SeqScape® Software, MicroSeq® ID Software, and Variant Reporter™ Software allow you to display or hide quality value (QV) bars in displays and printouts. You can customize the color and range for low, medium, and high quality values. For QVs ≤ 50, the length of a bar is proportional to the corresponding quality value. Quality values above 50 will have the same color and QV bar length as those defined for a QV of 50. To see the quality value for a particular base, place the computer mouse over the QV bar.</p> <p>In SeqScape® Software, MicroSeq® ID Software, and Variant Reporter™ Software, the per-base quality values also appear in the reports corresponding to bases identified as mutations.</p>																																				
<p>Why are the quality value bars displayed in gray?</p>	<p>A quality value is assigned to a specific basecall. When you change a basecall, the quality value does not apply to the new base, and therefore, it is displayed as a gray bar.</p> <p>Also when you reassign Ns to bases below a certain QV, the QV bar does not apply to the N basecall, and therefore it is displayed as a gray bar.</p>																																				
<p>Are quality value bars printed for the Electropherogram or Sequence views?</p>	<p>You can show or hide the QV bars when printing the Electropherogram and Sequence views of the sample file. QV bars are not printed if you print more than seven panels per page (due to space limitations). The quality value numbers cannot be printed.</p>																																				
<p>Which Life Technologies software can display the quality values?</p>	<p>Sequencing Analysis Software v5.X, Sequencing Analysis Software 6, SeqScape® Software v2.X, SeqScape® Software 3, MicroSeq® ID Software v1.X, v2.X, Variant Reporter™ Software v1.X, and Variant Reporter™ 2 can display quality values.</p>																																				
<p>Can I view quality values from KB™ Basecaller with other software?</p>	<p>Quality value graphics from KB™ Basecaller are customized for processing by other Life Technologies software. The KB™ Basecaller allows other Life Technologies software to perform additional functions, such as clear range trimming and more streamlined editing.</p>																																				

Notes

Miscellaneous FAQs

Some frequently asked questions regarding Ns, spacing values, and providing feedback are shown below.

Question	Answer
When do Ns appear in samples analyzed by the KB™ Basecaller Software?	<p>When using the KB™ Basecaller, the sequence “NNNNN” appears in the sample file when the sample fails analysis. Omit this file from further analysis. The Analysis Report in Sequencing Analysis Software will also flag these files.</p> <p>In addition to pure and mixed bases shown with QV bars, N's and gray QV bars are also shown when you reassign Ns to all bases before the user-specified QV threshold. This allows you to view the longer read length and more accurate basecalling of KB™ Basecaller while still viewing data with software that does not display QVs.</p>
Why does the spacing value sometimes appear in red?	When the ABI Basecaller fails to determine a spacing value for a sample file, it uses a default value of 12.00 for all run conditions. This number appears as in red in the Sample Manager, and the Annotation view displays “-12.00”.
Why does the spacing value sometimes have a negative value?	When the KB™ Basecaller fails to determine a spacing value for a sample file, it uses a default value specific to the instrument/polymer/chemistry/run condition used to generate the sample file. This value appears in red in the Sample Manager. The Annotation view displays -1 times this value.
How can I provide feedback to the KB™ Basecaller product team?	Email information to your local Life Technologies applications support representative at www.lifetechnologies.com/support . If applicable, please include sample files and details (including analysis settings) on how to reproduce your observation.

Notes _____

Conference posters and reference

- Posters**
- ABRF 2007 – *Improved Accuracy for Mutation and SNP Detection: Variant Reporter™ Software*, Ming Li et. al.
 - ESHG 2007 – Direct Sequencing Quality Control
 - AGBT 2004 – Longer Reads with the KB™ Basecaller
 - ABRF 2004 – Integrated Sequencing Analysis Solutions using the KB™ Basecaller from Applied Biosystems
 - ESHG 2009 Performance of the KB™ Basecaller for a New Sequencing System

These posters and other literature can be found at:

www.lifetechnologies.com

Click **Support**, then **Products and Technical Literature**. Search with the keyword *KB*.

- Reference** B. Ewing and P. Green, *Genome Research*, 8:186-194, 199.

Notes _____

Managing Data Collection Software Licenses

Manage software licenses

The 3730 Series Data Collection Software 4 requires a license to run.

IMPORTANT! If you replace or add a network card in the computer running the software, or relocate the software to a new computer, contact Life Technologies to update your license for the new network card or computer.

Obtain and activate a software license

The 3730 Series Data Collection Software 4 Software Activation dialog box is displayed when you start the software if no license is installed and activated on your computer.


This task is typically performed by the Life Technologies service representative during installation of the instrument.

1. Ensure that all network cards in the computer are enabled.

IMPORTANT! You can run the 3730 Series Data Collection Software 4 using only the network cards enabled when you activate the software license. For example, if you activate the software when your wireless network card is disabled, you will not be able to run the software when the wireless network card is enabled.

Notes _____

2. Display the Software Activation dialog box by starting the 3730 Series Data Collection Software 4.




3. Obtain the license key. The license key is provided on the 3730 Series Data Collection Software 4 CD case, or in an email from Life Technologies.
4. Request the software license file by performing steps 1a, 1b, and 1c as listed on the activation screen.

IMPORTANT! Keep a record of the email address used to activate the software license. You must use the same email address to renew the software license when it expires.

5. Obtain the software license file from your email.
6. Make a copy of the software license file and keep in a safe location.
7. Copy the software license file to the desktop of the 3730 Series Data Collection Software 4 computer.



Notes _____

8. If the Software Activation dialog box has closed, start the 3730 Series Data Collection Software 4 to open it.
9. Click **Browse**, then navigate to the software license file saved on your computer.
10. Click **Install and Validate License**. A message is displayed when the license is installed and validated.
11. Click **Close**.




Renew a software license

1. Ensure that all network cards in the computer are enabled.
2. Display the Software License Renewal dialog box by doing either of the following:
 - Select **Activate License Now** in the Warning: License Will Expire Soon dialog box that is displayed 8–30 days prior to expiration.
 - Click **Yes** in the Critical Warning: License Will Expire Soon dialog box that is displayed within 7 days of expiration.



Notes _____

3. Choosing to **Activate/Install License Now** will result in the display of 3730 Series Data Collection Software 4 box, shown here for the 3130. Click **Reinstall License** in the Lower Left Corner.



Notes _____

4. Complete the License Renewal dialog box as described below:

5. Enter the email address used to activate the software license.

IMPORTANT! You must use the same email address to activate and renew the software license. If you do not have the activation email address available, enter any email address, click the licensing link in the Software Renewal dialog box, then click **Contact Support** in the License Renewal web page displayed.

6. Request the renewed software license file by performing step **1c** as listed on the renewal screen.

7. Obtain the renewed software license file from your email.

8. Copy the renewed software license file to the desktop of this computer.

9. Click **Browse**, then navigate to the renewed software license file saved on your computer.

10. Click **Install and Validate License**. A message is displayed when the license is installed and validated.

11. Click **Close**.

3xxx Series Data Collection Software 4 Software License Renewal

1. Request license file for Computer ID:
002564ee13a4 002564ee13a5
This ID is unique to this computer and cannot be used to obtain a license file for another computer.

a. Enter the license key (from CD or email):
Dev-a5c5-50f3-c9da-4d1a-a918-0bfa-1667-226f

b. Enter your email address:

c. Is this computer currently connected to the internet?

2. Retrieve the license file from email, then save it to the desktop of this computer.

3. Find the license file:

4. Click

Notes _____

Appendix D Managing Data Collection Software Licenses

Renew a software license

Notes _____



WARNING

GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.













- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, etc). To obtain SDSs, see the “Documentation and Support” section in this document.
- All testing should be performed in accordance with local, regional and national acceptable laboratory accreditation standards and/or regulations.




Symbols on Instruments

Symbols may be found on the instrument to warn against potential hazards or convey important safety information. In this document, the hazard symbol is used along with one of the following user attention words described:

- **CAUTION!** – Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.
- **WARNING!** – Indicates a potentially hazardous situation that, if not avoided, could result in death or serious injury.
- **DANGER!** – Indicates an imminently hazardous situation that, if not avoided, will result in death or serious injury.



Symbol	English	Français
	Caution, risk of danger Consult the manual for further safety information.	Attention, risque de danger Consulter le manuel pour d'autres renseignements de sécurité.
	Caution, hot surface	Attention, surface chaude
	Caution, risk of electrical shock	Attention, risque de choc électrique
	Laser radiation	Rayonnement laser

Symbol	English	Français
	Caution, piercing hazard	Attention, danger de perforation
	Potential biohazard	Danger biologique potentiel
	Ultraviolet light	Rayonnement ultraviolet
	On	On (marche)
	Off	Off (arrêt)
	On/Off	On/Off (marche/arrêt)
	Standby	En attente
	Earth (ground) terminal	Borne de (mise à la) terre
	Protective conductor terminal (main ground)	Borne de conducteur de protection (mise à la terre principale)
	Terminal that can receive or supply alternating current or voltage	Borne pouvant recevoir ou envoyer une tension ou un courant de type alternatif
	Terminal that can receive or supply alternating or direct current or voltage	Borne pouvant recevoir ou envoyer une tension ou un courant continu ou alternatif
	Do not dispose of this product in unsorted municipal waste CAUTION! To minimize negative environmental impact from disposal of electronic waste, do not dispose of electronic waste in unsorted municipal waste. Follow local municipal waste ordinances for proper disposal provision and contact customer service for information about responsible disposal options.	Ne pas éliminer ce produit avec les déchets usuels non soumis au tri sélectif. CAUTION! Pour minimiser les conséquences négatives sur l'environnement à la suite de l'élimination de déchets électroniques, ne pas éliminer ce déchet électronique avec les déchets usuels non soumis au tri sélectif. Se conformer aux ordonnances locales sur les déchets municipaux pour les dispositions d'élimination et communiquer avec le service à la clientèle pour des renseignements sur les options d'élimination responsable.

Conformity mark	Description
	Indicates conformity with safety requirements for Canada and U.S.A.
	Indicates conformity with European Union requirements for safety and electromagnetic compatibility.
	Indicates conformity with Australian standards for electromagnetic compatibility.

Safety Alerts on Instruments

Additional text may be used with one of the symbols described above when more specific information is needed to avoid exposure to a hazard. See the following table for safety alerts found on the instrument.

English	French translation	Location on Instrument
 <p>DANGER! Class 3B (III) visible and/or invisible laser radiation present when open and interlocks defeated. Avoid exposure to beam.</p>	<p>DANGER! Rayonnement laser visible ou invisible de classe 3B (III) présent en position ouverte et avec les dispositifs de sécurité non enclenchés. Éviter toute exposition au faisceau.</p>	<p>Detection cell cover</p> 

Instrument Safety

General



CAUTION

Do not remove instrument protective covers. If you remove the protective instrument panels or disable interlock devices, you may be exposed to serious hazards including, but not limited to, severe electrical shock, laser exposure, crushing, or chemical exposure.



CAUTION

Solvents and Pressurized fluids. Wear eye protection when working with any pressurized fluids. Use caution when working with any polymeric tubing that is under pressure:

- Extinguish any nearby flames if you use flammable solvents.
- Do not use polymeric tubing that has been severely stressed or kinked.
- Do not use polymeric tubing with tetrahydrofuran or nitric and sulfuric acids
- Be aware that methylene chloride and dimethyl sulfoxide cause polymeric tubing to swell and greatly reduce the rupture pressure of the tubing.
- Be aware that high solvent flow rates (~40 mL/min) may cause a static charge to build up on the surface of the tubing and electrical sparks may result.

Physical injury



CAUTION

Moving and Lifting Injury. The instrument is to be moved and positioned only by the personnel or vendor specified in the applicable site preparation guide.


Improper lifting can cause painful and permanent back injury.

Things to consider before lifting or moving the instrument or accessories:


- Depending on the weight, moving or lifting may require two or more persons.
- If you decide to lift or move the instrument after it has been installed, do not attempt to do so without the assistance of others, the use of appropriate moving equipment, and proper lifting techniques.
- Ensure you have a secure, comfortable grip on the instrument or accessory.
- Make sure that the path from where the object is to where it is being moved is clear of


obstructions.


- Do not lift an object and twist your torso at the same time. Keep your spine in a good neutral position while lifting with your legs.
- Participants should coordinate lift and move intentions with each other before lifting and carrying.
- For smaller packages, rather than lifting the object from the packing box, carefully tilt the box on its side and hold it stationary while someone else slides the contents out of the box.


 **CAUTION Moving Parts.** Moving parts can crush, pinch and cut. Keep hands clear of moving parts while operating the instrument. Disconnect power before servicing.

Electrical


 **WARNING Fuse Installation.** Before installing the instrument, verify that the fuses are properly installed and the fuse voltage matches the supply voltage. Replace fuses only with the type and rating specified for the unit. Improper fuses can damage the instrument wiring system and cause a fire.


 **DANGER ELECTRICAL SHOCK HAZARD.** Severe electrical shock can result from operating the Applied Biosystems® 3730/3730xI DNA Analyzer without its instrument panels in place. Do not remove instrument panels. High-voltage contacts are exposed when instrument panels are removed from the instrument.

 **WARNING Voltage Selector Switch.** Before installing the instrument, verify that the voltage selector switch is set for the supply voltage. This will prevent damage to the instrument, reduce risk of fire, and enable proper operation.

 **WARNING Ensure appropriate electrical supply.** For safe operation of the instrument:

- Plug the system into a properly grounded receptacle with adequate current capacity.
- Ensure the electrical supply is of suitable voltage.
- Never operate the instrument with the ground disconnected. Grounding continuity is required for safe operation of the instrument.

 **WARNING Power Supply Line Cords.** Use properly configured and approved line cords for the power supply in your facility.

 **WARNING Disconnecting Power.** To fully disconnect power either detach or unplug the power cord, positioning the instrument such that the power cord is accessible.

Overvoltage Rating

The Applied Biosystems® 3730/3730xI DNA Analyzer has an installation (overvoltage) category of II, and is classified as portable equipment.

Cleaning and decontamination



CAUTION

Cleaning and Decontamination. Using a cleaning or decontamination method not specified by the manufacturer may result in damage to the equipment. For the protection of others, ensure the instrument is properly decontaminated prior to having the instrument serviced at your facility or before sending the instrument for repair, maintenance, trade-in, disposal, or termination of a loan. Decontamination forms may be requested from customer service.

Laser



WARNING

LASER HAZARD. Under normal operating conditions, the Applied Biosystems® 3730/3730x/ DNA Analyzer are categorized as a Class I laser product. However, removing the protective covers and (when applicable) defeating the interlock(s) may result in exposure to the internal Class 3B laser. Lasers can burn the retina, causing permanent blind spots. Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure. To ensure safe laser operation:

- Never look directly into the laser beam.
- Do not remove safety labels, instrument protective panels, or defeat safety interlocks.
- The system must be installed and maintained by an Life Technologies Technical Representative.

Life Technologies Technical Representatives are instructed to:

- Remove jewelry and other items that can reflect a laser beam into your eyes or those of others
- Wear proper eye protection and post a laser warning sign at the entrance to the laboratory if the laser protection is defeated for servicing.

DO NOT operate the laser when it cannot be cooled by its cooling fan; an overheated laser can cause severe burns on contact.

Note the laser warnings provided in [“Safety Alerts on Instruments” on page 197](#).



CAUTION

LASER HAZARD, Bar Code Scanner. The bar code scanner included with the instrument system is a Class 2 laser. To avoid damage to eyes, do not stare directly into the beam or point into another person's eyes.

Laser Classification

The 3730/3730x/ DNA Analyzer uses a laser. Under normal operating conditions, the instrument laser is categorized as a Class I laser. When safety interlocks are disabled during certain servicing procedures, the laser can cause permanent eye damage, and, therefore, is classified under those conditions as a Class 3B laser.

The Applied Biosystems® 3730/3730x/ DNA Analyzer has been tested to and complies with 21 CFR, 1040.10 and 1040.11, as applicable.

The 3730/3730x/ DNA Analyzer laser has been tested to and complies with standard EN60825-1, “Radiation Safety of Laser Products, Equipment Classification, Requirements, and User’s Guide.”

Safety and Electromagnetic Compatibility (EMC) Standards

The instrument design and manufacture complies with the standards and requirements for safety and electromagnetic compatibility as noted in the following table:

Safety

Reference	Description
EU Directive 2006/95/EC	European Union “Low Voltage Directive”
IEC 61010-1 EN 61010-1 CSA C22.2 No. 61010-1	<i>Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements</i>
IEC 61010-2-010 EN 61010-2-010 UL 61010-2-010	<i>Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials</i>
IEC 61010-2-081 EN 61010-2-081 UL 61010-2-081	<i>Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes</i>
IEC 60825-1 EN 60825-1	<i>Safety of laser products – Part 1: Equipment classification and requirements</i>
21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007, as applicable	U.S. FDA Health and Human Services (HHS) “Radiological health performance standards for laser products” and “Radiological health performance standards for specific purpose laser products”


EMC

Reference	Description
Directive 2004/108/EC	European Union “EMC Directive”
EN 61326-1	<i>Electrical Equipment for Measurement, Control and Laboratory Use – EMC Requirements – Part 1: General Requirements</i>
FCC Part 18 (47 CFR)	U.S. Standard “Industrial, Scientific, and Medical Equipment”
AS/NZS 2064	<i>Limits and Methods of Measurement of Electromagnetic Disturbance Characteristics of Industrial, Scientific, and Medical (ISM) Radiofrequency Equipment</i>
ICES-001, Issue 3	<i>Industrial, Scientific and Medical (ISM) Radio Frequency Generators</i>


Environmental design


Reference	Description
Directive 2002/96/EC	European Union “WEEE Directive” – Waste electrical and electronic equipment

Chemical safety

 **WARNING GENERAL CHEMICAL HANDLING.** To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the “Documentation and Support” section in this document.
 - Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
 - Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
 - Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended in the SDS.
 - Handle chemical wastes in a fume hood.
 - Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
 - After emptying a waste container, seal it with the cap provided.
 - Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
 - Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
 - **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.
-

 **WARNING HAZARDOUS WASTE (from instruments).** Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.

 **WARNING 4L Reagent and Waste Bottle Safety.** Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.

Biological hazard safety



WARNING BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Individuals should be trained according to applicable regulatory and company/institution requirements before working with potentially infectious materials. Read and follow the applicable guidelines and/or regulatory requirements in the following:

In the U.S.:

- U.S. Department of Health and Human Services guidelines published in Biosafety in Microbiological and Biomedical Laboratories found at: www.cdc.gov/biosafety
- Occupational Safety and Health Standards, Bloodborne Pathogens (29 CFR§1910.1030), found at: www.access.gpo.gov/nara/cfr/waisidx_01/29cfr1910a_01.html
- Your company's/institution's Biosafety Program protocols for working with/handling potentially infectious materials.
- Additional information about biohazard guidelines is available at: www.cdc.gov

In the EU:

- Check local guidelines and legislation on biohazard and biosafety precaution and refer to the best practices published in the World Health Organization (WHO) Laboratory Biosafety Manual, third edition, found at: www.who.int/csr/resources/publications/biosafety/WHO_CDS_CSR_LYO_2004_11/en/

Documentation and Support

Related documentation

The following related documents are shipped with the system:

Document title	Pub. Part no.
<i>Applied Biosystems® 3730/3730xl DNA Analyzer Maintenance and Troubleshooting Guide</i>	4477797
<i>Applied Biosystems® 3730/3730xl DNA Analyzer Quick Reference Card</i>	4477852
<i>Applied Biosystems 3730/3730xl DNA Analyzer and 3130/3130xl Genetic Analyzers AB Navigator Software Administrator Guide</i>	4477853

Portable document format (PDF) versions of this guide and the documents listed above are also available on the Applied Biosystems® 3730 Series Data Collection Software 4 CD.

Note: To open the user documentation included on the Applied Biosystems® 3730 Series Data Collection Software 4 CD, use the Adobe® Reader® software available from www.adobe.com.

Note: For additional documentation, see “Obtaining Support” on page 204.

Obtaining SDSs

Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Note: For the SDSs of chemicals not distributed by Life Technologies Corporation, contact the chemical manufacturer.

Obtaining Support

For the latest services and support information for all locations, go to:

www.lifetechnologies.com/support

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

Computer Configuration

Life Technologies Corporation supplies or recommends certain configurations of computer hardware, software, and peripherals for use with its instrumentation. Life Technologies Corporation reserves the right to decline support for or impose extra charges for supporting nonstandard computer configurations or components that have not been supplied or recommended by Life Technologies Corporation. Life Technologies Corporation also reserves the right to require that computer hardware and software be restored to the standard configuration prior to providing service or technical support. For systems that have built-in computers or processing units, installing unauthorized hardware or software may void the Warranty or Service Plan.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies Corporation' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies Corporation at www.lifetechnologies.com/support.

Symbols

- .fsa files
 - from GeneMapper-Generic 102

Numerics

- 3730/3730xl Data Collection, starting 17

A

- analysis protocol
 - creating 80
 - creating for autoanalysis 79
 - deleting 84
 - editing 83
 - exporting 84
 - importing 85
- anode buffer jar, filling 28
- Any4Dye
 - creating a spectral calibration for 168
 - how to use 168
- Any5Dye
 - creating a spectral calibration for 168
 - how to use 168
- array port, illustration of 11
- Array View, viewing data in 153
- auto mode 142
 - features 144
 - scheduling runs with 144

B

- biohazard safety 202
- buffer fill-line 11
- buffer jar, illustration of 11
- buffer reservoir
 - assembly 25
 - filling 24
- buffer valve pin, illustration of 11

C

- capillary array
 - illustration of 11
 - installing 18
- capillary array knob, illustration of 11

- capillary array tip, illustration of 11
- Change Polymer Type Wizard, using 21
- check valve, diagram of 11
- chemical safety 201
- cleaning safety 199
- computer
 - configuration requirement 204
 - start up and log on 12
 - technical support for altered configuration 204
- customizing run modules 112

D

- Data Collection software
 - starting 102
- data, viewing in array view 153
- decontamination safety 199
- double-tapered ferrule, illustration of 11
- dye set G5-RCT
 - recommendations for use 166, 176, 189
- dyeset/primer files
 - list of 81

E

- electrical safety 198
- electrode, illustration of 11
- electromagnetic compatibility (EMC) standards 200

F

- file naming
 - acceptable characters 70, 102
 - invalid characters 70, 102
- fill down special 96, 124
- fragment analysis, creating required settings for 109

G

- GeneMapper 106
 - plate record 105
- GeneMapper-Generic, .fsa files 102

H

- heat-sealed plates 129

- I**
- instrument
 - illustrated parts of 9
 - operation, manual vs auto mode 142
 - startup 12
 - instrument protocol
 - creating for fragment analysis 109
 - creating for sequencing 74
 - importing 110
 - interconnect tube 11
- K**
- KB[®] Basecaller 80, 81, 175
- L**
- laser safety 199
 - license. *See* software license 189
 - load maps for 48-capillaries and 98-capillaries 137-139
 - lower polymer block, illustration of 11
 - Luer fitting, illustration of 11
- M**
- magnifying
 - spatial profiles 34
 - spectral profiles 60
 - manual mode
 - scheduling runs using 142
 - versus auto mode 142
 - mounting pin, illustration of 11
- O**
- O-ring, illustration of 11
 - overflow hole, illustration of 11
- P**
- password 12
 - pausing a run 147
 - PDP motor cover, illustration of 11
 - physical injury safety 197
 - piston, illustration of 11
 - plate record
 - creating for sequencing analysis 93
 - creating GeneMapper 121
 - for fragment analysis 105
 - for sequencing analysis 70
 - GeneMapper elements of 105
 - when to create 70, 105
 - plate run
 - spectral, using Any4Dye 168
 - stopped before autoextraction is complete 156
 - using Any4Dye with 172
 - using Any5Dye with 172
 - plates
 - assembling 128
 - components 128
 - heat-sealed 129
 - septa-sealed 130
 - polymer
 - adding 18
 - changing to new type 21
 - replenishing or changing 20
 - Polymer Delivery Pump (PDP), illustration of 11
 - polymer supply bottle cap with hole, illustration of 11
 - polymer supply bottle, illustration of 11
 - polymer supply tube illustration of 11
 - profile
 - passing spatial, examples of 37
 - spatial calibration, evaluating 34
 - pump block, illustration of 11
 - pump chamber, illustration of 11
- R**
- Replenish Polymer Wizard
 - using 21
 - reservoirs
 - filling 22, 24
 - placing into instrument 27
 - results group
 - creating for autoanalysis 114
 - creating for sequencing 85
 - exporting 119
 - importing 119
 - run
 - starting, stopping, skipping, pausing 147
 - run buffer
 - preparing 22
 - run history view, viewing data in 154
 - run modules
 - customizing 112
 - editable parameters 78
 - selecting for sequencing 75
- S**
- safety
 - alerts on instrument 197
 - biohazard 202
 - chemical 201
 - cleaning and decontamination 199
 - electrical 198
 - laser 199
 - physical injury 197

- standards 200
- sample file name, creating 117
- sample run
 - adding 99
 - adding for fragment analysis 126
- septa-sealed plates 129, 130
- sequencing
 - plate editor 71
 - run modules 75
- sequencing spectral calibrations
 - passing, examples of 61
- service console, using 17
- settings
 - required for automated fragment analysis 109
 - required for automated sequencing analysis 74
- software, Data Collection 102
- software, license
 - manage 189
 - obtain 189
 - renew 191
- spatial calibration
 - evaluating profile 34
 - performing 32
 - what it tells you 32
 - when to perform 32
- spatial profile
 - magnifying 34
 - passing, examples of 37
- spectral calibration
 - evaluating results 57
 - performing 43
 - spectral viewer 57
 - starting a run 54
 - troubleshooting 67
- spectral calibration, passing 63
- spectral profile
 - magnifying 60
- spectral run
 - using Any4Dye 168
- spectral viewer, selecting active spectrals 64
- starting
 - Data Collection software 102
 - instrument 12
 - run 147
- stopping a run 147

T

- technical support, for computers with altered configuration 204
- toolbar 147
- troubleshooting
 - solid red light 15
 - solid yellow light 15

U

- user name 12

W

- warranty
 - for computers with altered configuration 204
- waste reservoir assembly, illustration of 25
- water reservoir assembly
 - filling 24
 - illustration of 25
- water seal, illustration of 11
- water trap, illustration of 11
- wizard
 - Replenish Polymer Wizard, using 21



Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.lifetechnologies.com/support

www.lifetechnologies.com

