

SNaPshot® Multiplex System for SNP genotyping

One system, many applications

Features of the SNaPshot® Multiplex System:

- Useful for a multitude of applications—SNP and methylation analysis, fingerprinting, quantitative allele frequency, validating NGS data, and more
- Customized for your target
- Offers multiplexing capability (up to 10-plex)
- Sensitive allele-frequency detection (typically 5%)
- Compatible with all CE instruments
- Automated analysis using specific GeneMapper® software data analysis module

The SNaPshot® Multiplex System is a primer extension–based method developed for the analysis of single nucleotide polymorphisms (SNPs) (Figure 1). Through its multiplexing capability, up to 10 SNPs can be analyzed in a single reaction, regardless of their positions on the chromosome or the amount of separation from neighboring SNP loci. The ability to use unlabeled, user-defined primers allows researchers to incorporate SNPs of interest cost-effectively. The Multiplex Ready Reaction Mix (included in the system) helps ensure

SNaPshot® Kit Single-Base Extension Labeling Chemistry




Figure 1. SNaPshot® labeling chemistry relies on single-base extension and termination. The SNaPshot® Multiplex Kit uses a single-tube reaction to interrogate SNPs at known locations. The chemistry is based on the dideoxy single-base extension of an unlabeled oligonucleotide primer (or primers). Each primer binds to a complementary template in the presence of fluorescently labeled ddNTPs and DNA polymerase. The polymerase extends the primer by one nucleotide, adding a single ddNTP to its 3' end. The fluorescence color readout reports which base was added.

robust, reproducible analyses of multiplexed samples. Researchers can analyze more than 23,000 SNP genotypes per day on just one 3730xl Genetic Analyzer (Table 1).

SNP analysis

SNPs have been identified in all genomes and can be used for a multitude of analyses, including studying mutations implicated in various cancers, genetic disease research, mitochondrial DNA investigations, scrapie susceptibility in sheep, loss of heterozygosity, assessing performance

in food animal production, and even differentiating drug and non-drug forms of *Cannabis*. Table 2 lists some of the many ways our customers use the SNaPshot® Multiplex System in their research (Table 2). A more extensive list of publications demonstrating the utility of the SNaPshot® Multiplex System is available at www.lifetechnologies.com/snapshot.

Table 1. SNP analysis throughput on various Applied Biosystems® instruments. Numbers stated are based on 24-hour operation and 10 SNPs per capillary.

	310	3130	3130xl	3500	3500xL	3730	3730xl
SNP throughput per day	480	3,840	15,360	3,840	11,520	21,120	42,240

Table 2. Selected published articles that cite the SNaPshot® Multiplex System. Find more publications at www.lifetechnologies.com/snapshot.

Type of analysis	Application	Title	Journal
SNP analysis	Differentiating drug and non-drug forms of <i>Cannabis</i>	Differentiation of drug and non-drug <i>Cannabis</i> using a single nucleotide polymorphism (SNP) assay.	<i>Forensic Sci Int</i> 207:193 (2011)
	Antiretroviral resistance mutation	Pre-screening HIV-1 reverse transcriptase resistance mutations in subtype B patients using a novel multiplex primer extension assay.	<i>Curr HIV Res</i> 7:398 (2009)
	<i>Arabidopsis</i> markers	Establishment of a high-efficiency SNP-based framework marker set for <i>Arabidopsis</i> .	<i>Plant J</i> 36:122 (2003)
	Assess meat performance in pigs	SNaPshot minisequencing and a panel of candidate genes for complex routine testing of meat performance traits in pigs.	<i>Anim Biotechnol</i> 18:109 (2007)
	Blood typing	Detection of blood group genes using multiplex SNaPshot method.	<i>Transfusion</i> 49:2012 (2009)
	Method for degraded samples	Developing multiplexed SNP assays with special reference to degraded DNA templates.	<i>Nat Protoc</i> 1:1370 (2006)
	Heteroplasmy validation after NGS analysis	Detecting heteroplasmy from high-throughput sequencing of complete human mitochondrial DNA genomes.	<i>Am J Hum Genet</i> 87:237 (2010)
	Identifying species in <i>M. tuberculosis</i> complex	Identification and genotyping of <i>Mycobacterium tuberculosis</i> complex species by use of a SNaPshot Minisequencing-based assay.	<i>J Clin Microbiol</i> 48:1758 (2010)
	Loss of heterozygosity (LOH)	Multiplex SNaPshot genotyping for detecting loss of heterozygosity in the mismatch-repair genes MLH1 and MSH2 in microsatellite-unstable tumors.	<i>Clin Chem</i> 54:1844 (2008)
	Mitochondrial variants	Mitochondrial variants in schizophrenia, bipolar disorder, and major depressive disorder.	<i>PLoS One</i> 4:e4913 (2009)
	Pathogen resistance	Single nucleotide polymorphisms and haplotypes in the IL10 region associated with HCV clearance.	<i>Genes Immun</i> 6:347 (2005)
	Scrapie susceptibility screening	Primer extension assay for prion protein genotype determination in sheep.	<i>Mol Cell Probes</i> 18:33 (2004)
	Selective breeding in horses	Identification of horse chestnut coat color genotype using SNaPshot.	<i>BMC Res Notes</i> 2:255 (2009)
	Species identification within <i>Lactobacillus casei</i> isolates	Application of the SNaPshot minisequencing assay to species identification in the <i>Lactobacillus casei</i> group.	<i>Mol Cell Probes</i> 25:153 (2011)
	Subspecies identification in tigers	The development and validation of a single SNaPshot multiplex for tiger species and subspecies identification—Implications for forensic purposes.	<i>Forensic Sci Int Genet</i> [2011] [Epub ahead of print]
	Tracking cows and their products	SNPmplexViewer—toward a cost-effective traceability system.	<i>BMC Res Notes</i> 4:146 (2011)
Wheat breeding	Insertion site-based polymorphism markers open new perspectives for genome saturation and marker-assisted selection in wheat.	<i>Plant Biotechnol J</i> 8:196 (2010)	

Type of analysis	Application	Title	Journal
SNP quantification	Quantitative mutant measurement	Putative precursor cancer cells in human colorectal cancer tissue.	<i>Int J Clin Exp Pathol</i> 2:154 (2009)
	Quantitative SNP allele frequency measurement	Universal, robust, highly quantitative SNP allele frequency measurement in DNA pools.	<i>Hum Genet</i> 110:471 (2002)
	Allele-specific gene expression	Allelic variation in human gene expression. <i>cis</i> -acting variation in the expression of a high proportion of genes in human brain.	<i>Science</i> 297:1143 (2003) <i>Hum Genet</i> 113:149 (2003)
BAC fingerprinting	High-throughput fingerprinting of bacterial artificial chromosomes	High-throughput fingerprinting of bacterial artificial chromosomes using the snapshot labeling kit and sizing of restriction fragments by capillary electrophoresis.	<i>Genomics</i> 82:378 (2003)
Methylation	Methylation status	Single nucleotide extension technology for quantitative site-specific evaluation of metC/C in GC-rich regions.	<i>Nucleic Acids Res</i> 33:e95 (2005)
	Methylation status	Capillary electrophoretic analysis of methylation status in CpG-rich regions by single-base extension of primers modified with N6-methoxy-2,6-diaminopurine.	<i>Anal Biochem</i> 380:13 (2008)

Methylation

The study of epigenetic effects, including methylation, is emerging as an important component of genetic research. In a typical assay to detect methylation, bisulfite treatment of DNA deaminates unmethylated cytosines (converting them to uracils) and leaves methylated cytosines unchanged. The subsequent PCR amplification step converts uracil bases to thymines. Researchers can use the SNaPshot® Multiplex System to quantify the cytosine-to-thymine changes in treated and untreated samples to determine the methylation status.

BAC fingerprinting

Fingerprinting of large insert genomic fragment libraries, also known as BAC (bacterial artificial chromosomes) clones leads to the construction of a genome-wide physical map. These maps are critical to genome sequencing, positional cloning, and understanding the relative organization of genes and markers. When BAC libraries are arranged into maps that reflect the DNA sequence in a chromosome, they provide maximal information and utility. The SNaPshot® Multiplex Kit is a widely used, efficient method to label BAC fragments. The labeled fragments can then be

separated and detected on any Applied Biosystems® capillary electrophoresis instrument. Sizing information from GeneMapper® Software is imported into subsequent editing and contig-assembly programs. The high-quality results you get from the SNaPshot® Multiplex Kit provide an easy-to-use and cost-effective solution for high-throughput BAC fingerprinting.

A complete solution

Across the SNP analysis workflow, Life Technologies offers solutions to help you get the research data you need. Key components include:

SNaPshot® Multiplex System

The system supplies the SNaPshot® Multiplex Ready Reaction Mix, control primer mix, control template, and protocol.

SNaPshot® Primer Focus® Kit

Designed to determine the approximate fragment sizes generated by various primers prior to performing SNP genotyping multiplex reactions (critical if two oligos will produce overlapping signals when run simultaneously) and enables the setting of tight loci windows (bins) in GeneMapper® Software.

The GeneScan™ 120 LIZ® Size Standard

This standard is a LIZ® dye-labeled size standard that is designed for reproducible sizing of small-fragment analysis data generated with the SNaPshot® Multiplex System. It accurately sizes samples ranging from 20 to 120 nucleotides. When used with GeneMapper® Software, the GeneScan™ 120 LIZ® Size Standard eliminates the need for manual allele calls.

Matrix Standard Set DS-02

This standard set is used to generate the “multicomponent matrix” needed to detect the dyes (dR110, dR6G, dTAMRA™, dROX™, and LIZ® dyes), critical for successful five-dye SNP analysis on all Applied Biosystems® genetic analyzers.

GeneMapper® Software

GeneMapper® Software specializes in multi-application functionality, including SNP genotyping analysis for data generated with the SNaPshot® Multiplex System.

Data analysis simplified using GeneMapper® Software

GeneMapper® Software is a flexible genotyping software package that provides DNA sizing and quality allele calls for all Life Technologies® electrophoresis-based genotyping systems. This software specializes in multi-application functionality,

including SNP genotyping analysis (Figure 2). GeneMapper® Software can help users increase data processing efficiency with remote auto-analysis and command line operation, and allows for multi-user, client-server deployment. The software uses Process Quality Values (PQVs) for automated identification that reduces

data review time for high-throughput genotyping. In addition, the security and audit features help users meet 21 CFR 11 requirements. Peak Scanner™ Software v1.0 is an alternative free analysis software for low-complexity analysis. Download for free at www.lifetechnologies.com.




Figure 2. Example of SNP genotyping results obtained from GeneMapper® Analysis Software. The vertical colored boxes are bins created automatically by the software using a sample or an artificial extension product created using the SNaPshot® Primer Focus® Kit. Each bin defines the minimum and maximum allowable size for each allele. GeneMapper® Software identifies each peak and assigns the corresponding allele. In this example, the sample is heterozygous G/A for SNP 28. GeneMapper® Software generates a customized report with genotype result (including size, height, peak area) and provides information about the quality of the sample (process quality values (PQVs) in columns OS through GQ), facilitating review by the operator. For example, a yellow triangle in the peak height ratio column (PHR) indicates that the ratio between the 2 alleles is outside what is defined in the default analysis settings (50% for each allele). Unfilled peaks (outlined only) represent the GeneScan™ 120 LIZ® Size Standard positions.

Ordering information

Product	Quantity	Cat. No.
SNaPshot® Primer Focus® Kit	100 rxns	4329538
SNaPshot® Multiplex Kit	100 rxns	4323159
SNaPshot® Multiplex Kit	1,000 rxns	4323161
SNaPshot® Multiplex Kit	5,000 rxns	4323163
GeneScan™ 120 LIZ® Size Standard	800 analyses	4324287
DS-02 Matrix Standard Set for Dye Set E5	8 analyses	4323014
GeneMapper® Software v4.1.1	1	4366925
GeneMapper® Software v4.1.1, 30-day Demo	1	4366851
Peak Scanner™ Software v1.0 (free download)	1	4381867

For research use only. Not intended for any animal or human therapeutic or diagnostic use. ©2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. C0110582 0112