# SeqStudio™ Genetic Analyzer Cathode Buffer Container

SeqStudio™ Genetic Analyzer

Catalog Numbers A33401

Pub. No. MAN0016142 Rev. B.0



**WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **thermofisher.com/support**.

#### Product description

The SeqStudio™ Genetic Analyzer Cathode Buffer Container (CBC) contains running buffer for capillary electrophoresis. The container has two compartments. The rear compartment provides the cathode buffer for electrophoresis. The front compartment is for capillary wash and waste.

The CBC is a ready-to-use, disposable container with a radio frequency identification (RFID) tag incorporated into the label. The instrument uses the RFID tag to track buffer usage and expiration.

See the SeqStudio™ Genetic Analyzer Instrument and Software User Guide (Pub. No. MAN0016138) for usage instructions.

## **Contents and storage**

Description	Amount	Storage
Cathode Buffer Container (CBC)	1	<ul> <li>Store at 2-8°C.</li> <li>See the expiration date on the package. Do not use an expired Cathode Buffer Container.</li> </ul>

### Required materials not supplied

Unless otherwise indicated, all materials are available through **thermofisher.com**. MLS: Fisher Scientific (**fisherscientific.com**) or other major laboratory supplier.

Item	Source	
Reservoir septa	A35640	

## Usage limits

**Note:** The system software automatically tracks the time that is elapsed from the first installation of the SeqStudio<sup>™</sup> Genetic Analyzer Cathode Buffer Container and the number of injections that are performed using the Cathode Buffer Container. The Cathode Buffer Container expires whenever a usage limit is reached.

**IMPORTANT!** The software displays a warning message when a usage limit is reached. If you intend to continue the run without replacing the Cathode Buffer Container, see "Important notice regarding use of consumables that exceed supported limits" on page 2 before proceeding.


Criterion	On-instrument usage limit
Time	14 days from first installation
Injections	125 injections

# Assemble the SeqStudio™ Genetic Analyzer Cathode Buffer Container (CBC)

Equilibrate the CBC to room temperature (15 minutes to overnight) before assembling.

On a clean and level surface:


- 1. Ensure that the buffer is at or above the fill line.
- 2. Carefully peel off the seal.



3. Wipe off any buffer on top of the CBC with a lint-free tissue. Ensure that the top of the container is dry.



4. Place the reservoir septum on the CBC, then press firmly so that there are no gaps in the seal.



Note: The CBC is filled significantly above Fill Line to account for evaporation. Replace the CBC when the fluid level drops below the Fill Line.

# Important notice regarding use of consumables that exceed supported limits

BEFORE DISMISSING THE WARNING THAT THE CONSUMABLES HAVE REACHED SUPPORTED LIMITS AND CONTINUING WITH OPERATION OF THE INSTRUMENT, PLEASE READ AND UNDERSTAND THE FOLLOWING IMPORTANT NOTICE AND INFORMATION:

Life Technologies does not recommend the use of consumables that exceed supported limits. The recommended limits are designed to promote the production of high quality data and minimize instrument downtime. Reagent and consumable lifetime minimum performance are based on testing and studies that use reagents and consumables that have not exceeded supported limits.

The use of consumables beyond the supported limits may impact data quality or cause damage to the instrument or capillary array. The cost of repairing such damage is NOT covered by any Life Technologies product warranty or service plan. Customer use of expired consumables is at customer's own risk and without recourse to Life Technologies. For example, product warranties do not apply to defects resulting from or repairs required due to misuse, neglect, or accident including, without limitation, operation outside of the environmental or use specifications or not in conformance with Life Technologies instructions for the instrument system, software, or accessories.

Please see your specific service contract or limited product warranty for exact language regarding coverage and ask your Life Technologies representative if you have further questions.

#### Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Manufacturer: Life Technologies Corporation | 35 Wiggins Avenue | Bedford, MA 01730

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0016142

Revision	Date	Description
B.0	16 June 2017	New document

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

