

Zero Blunt® TOPO® PCR Cloning Kit

Catalog Number K2800-20SC

Doc. Part no. 25-0649

Pub. Part no. MAN0003551

Rev. A.0

Introduction

Follow these instructions to TOPO® Clone your blunt-end PCR product into pCR™-Blunt II-TOPO® and transform the reaction into chemically competent *E. coli* cells. For transformation of electrocompetent cells, a diagram of the multiple cloning site, and a manual, see www.lifetechnologies.com/support or contact Technical Support.

Produce Blunt PCR Products

Use the appropriate protocol to produce blunt-end PCR products using a thermostable proofreading polymerase. Make sure that the final extension step is sufficient to fully extend your PCR products (7–30 minutes).

TOPO® Cloning Reaction

1. Set up the following 6 μ L TOPO® Cloning reaction:

Reagent	Amount*
Fresh PCR Product	0.5–4 μ L
Salt Solution	1 μ L
Sterile Water	add to a total volume of 5 μ L
pCR™-Blunt II-TOPO®	1 μ L
Final Volume	6 μL

* For transformation of chemically competent *E. coli* only.

- Mix gently and incubate for 5 minutes at room temperature.
- Place tubes on ice. **Proceed to Transformation and Analysis.**

Transformation and Analysis

Follow the protocol in this section to transform chemically competent cells and to analyze positive clones. To transform electrocompetent cells, refer to the Zero Blunt® TOPO® PCR Cloning Kit manual.

One Shot® Chemical Transformation

1. Thaw 1 vial of One Shot® *E. coli* cells **on ice** for each transformation.
2. Add 2 µL of the TOPO® Cloning reaction to each vial of One Shot® cells to be transformed, and mix gently.
3. Incubate the vials on ice for 5–30 minutes.
4. Heat-shock the cells for 30 seconds at 42°C without shaking.
5. Add 250 µL of room temperature S.O.C. medium to the cells.
6. Cap the tubes and shake them at 37°C for 1 hour.
7. Spread 10–50 µL from each transformation on pre-warmed LB plates containing 50 µg/mL kanamycin or pre-warmed Low Salt LB plates containing 25 µg/mL Zeocin™ selective antibiotic. Refer to the Zero Blunt® TOPO® PCR Cloning Kit manual for a Low Salt LB medium recipe.
8. Incubate plates overnight at 37°C.


An efficient TOPO® Cloning reaction should produce several hundred colonies. Pick ~10 colonies for analysis. **Proceed to Analyze Positive Clones.**

Transformation and Analysis, Continued

Analyze Positive Clones

1. Culture the 10 colonies overnight in LB medium containing 50 µg/mL kanamycin or Low Salt LB medium containing 25 µg/mL Zeocin™ selective antibiotic.
2. Isolate plasmid DNA using your method of choice. For ultra-pure plasmid DNA, we recommend the PureLink® HQ Mini Plasmid Purification Kit (Cat. no. K2100-01).
3. Analyze the plasmid by restriction analysis or by sequencing to confirm the presence and orientation of the insert.

Map of pCR™-Blunt II-TOPO®



lac promoter/operator region: bases 95-216

M13 Reverse priming site: bases 205-221

lacZ-alpha ORF: bases 217-576

SP6 promoter priming site: bases 239-256

Multiple Cloning Site: bases 269-399

TOPO®-Cloning Site: bases 336-337

T7 promoter priming site: bases 406-425

M13 (-20) Forward priming site:
bases 433-448

Fusion joint: bases 577-585

ccdB lethal gene ORF: bases 586-888

kan promoter: bases 1099-1236

Kanamycin resistance ORF: bases 1237-2031

Zeocin™ resistance ORF: bases 2238-2612

pUC origin: bases 2724-3397

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

DISCLAIMER

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Life Technologies is a Thermo Fisher Scientific brand. © 2014 Thermo Fisher Scientific Inc. All rights reserved. Zeocin is a trademark of CAYLA, SA. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA

Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com

25 February 2014

