CERTEST

Salmonella

ONE STEP *Salmonella* CARD TEST

CERTEST BIOTEC S.L.

For information purposes only

CERTEST *Salmonella*

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ENGLISH

INTENDED USE

CerTest *Salmonella* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Salmonella* in stool samples and *Salmonella* suspected colonies in stool culture.

CerTest *Salmonella* card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of *Salmonella* infection (salmonellosis) and it could be used to identify of suspected isolates of *Salmonella* from selective media (stool culture).

INTRODUCTION

A total of 2501 different *Salmonella* serotypes have been identified up to 2004. *Salmonella enteritidis* and *Salmonella typhimurium*, the two most important serotypes for salmonellosis transmitted from animals to humans.

Salmonella infections are acquired by eating contaminated poultry, eggs or dairy products.

The *Salmonella spp.* also spreads easily from raw or undercooked poultry to innocent vegetables, fruit or other foods via contaminated hands, knives, countertops or cutting boards. Due to the ability of *Salmonella* to multiply in a wide variety of foods, it is important to be able to isolate the organisms even when present in very small numbers in the faeces.

The symptoms of *Salmonella* infection are abdominal pain, diarrhoea, mild fever, chills, headache, nausea and vomiting, developing 12–72 hours (but occasionally as long as 7 days) after infection. The discomfort generally lasts a few days. It can be dangerous for the elderly, infants and the immunocompromised, who may become extremely ill. *Salmonella* is also one of the leading predictors for reactive arthritis, a painful, chronic and potentially debilitating condition that causes joint inflammation.

TEST PRINCIPLE

CerTest *Salmonella* is based on the principle of a qualitative immunochromatographic assay for the determination of *Salmonella* in stool samples and *Salmonella* suspected colonies in stool culture.

The strip consists of a nitrocellulose membrane pre-coated with polyclonal antibodies on the test line (T), in the results window, against *Salmonella* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (polyclonal antibodies anti-*Salmonella*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Salmonella* polyclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Salmonella* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Salmonella* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Salmonella* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Salmonella* test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

Suspected *Salmonella* colonies in stool culture


The Hektoen Enteric Agar (HEA) is the method of choice for the isolation of *Salmonella* (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation in Hektoen the typical *Salmonella* colonies will be green with black centre due to H₂S production.

Salmonella colonies in HEA

1. Examine Hektoen plates after 1 day incubation. Select *Salmonella* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella* colonies (2) and add them to the collection tube.
2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

Note: Because it is possible for the *Salmonella* content in the stool to be too small, a negative result cannot be considered to be the same as a negative antigen finding. It is recommended to carry out enrichment in Rappaport-Vassiliadis medium (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation, spread in Hektoen Enteric Agar and incubate (aerobic atmosphere, 24h/37°C).


Salmonella Shigella Agar (SS) is another method for the isolation of *Salmonella* (aerobic atmosphere, 24-48 hours/35-37°C). After the incubation period in SS the typical *Salmonella* colonies will be transparent without any black center. The SS agar is a very inhibitor that sometimes prevents the growth of certain species of *Salmonella*. We recommend using this medium unless accompanied by another inhibitor (HEA).

Salmonella colonies in SS

1. Examine SS plates after incubation period. Select *Salmonella* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella* colonies (2) and add them to the collection tube.
2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

MATERIALS

MATERIALS PROVIDED

- CerTest *Salmonella* card tests
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer
- Hektoen Enteric Agar (HEA)/*Salmonella* Shigella Agar (SS)

TEST PROCEDURE


Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the collection tube in order to assure good sample dispersion.
2. Remove CerTest *Salmonella* card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>Salmonella</i>	Interpretation of the results
1.	- GREEN	There is no <i>Salmonella</i> presence. No infection caused by <i>Salmonella</i> .
2.	+ GREEN-RED	There is <i>Salmonella</i> presence. <i>Salmonella</i> infection, which might mean abdominal pain, diarrhoea, mild fever, chills, headache, nausea and vomiting, develops 12–72 hours (but occasionally as long as 7 days) after infection.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentrations to faint when the antigens concentration is close to the detection limit value of the test.
- CerTest *Salmonella* should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *Salmonella* in faecal samples, nevertheless, it can be due to several causes and/or species. A positive result should be followed up with additional laboratory techniques (biochemical and serological methods or by PCR) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a *Salmonella* determination should be carried out on a sample from an enrichment culture.
- Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Salmonellosis constitutes a major public health burden and represents a significant cost to society in many countries. Although outbreaks usually attract media attention, studies indicate that more than 80% of all salmonellosis cases occur individually rather than as outbreaks. In the United States of America, an estimated 1.4 million non-typhoidal *Salmonella* infections, resulting in 168000 visits to physicians, 15000 hospitalizations and 580 deaths annually.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit range for the different species are:

Salmonella enteritidis: $5 \cdot 10^7 - 3.12 \cdot 10^6$ CFU/mL of inactivated *Salmonella enteritidis* antigen cells culture.

Salmonella typhimurium: $10 \cdot 10^6 - 6.25 \cdot 10^5$ CFU/mL of inactivated *Salmonella typhimurium* antigen cells culture.

Salmonella typhi: $10 \cdot 10^7 - 6.25 \cdot 10^6$ CFU/mL of inactivated *Salmonella typhi* antigen cells culture.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Salmonella*, CerTest) and qPCR technique (VIASURE *Salmonella*, *Campylobacter* & *Y. enterocolitica* Real Time PCR Detection kit). The results were as follows:

		qPCR test: Viasure <i>Salmonella</i> , <i>Campylobacter</i> & <i>Y. enterocolitica</i> Real Time Detection kit		
		+	-	Total
IC test: CerTest <i>Salmonella</i>	+	47	2	49
	-	15	103	118
	Total	62	105	167

CerTest <i>Salmonella</i> vs Viasure <i>Salmonella</i> , <i>Campylobacter</i> & <i>Y. enterocolitica</i> Real Time PCR Detection Kit		
	Mean Value	95% CI (confidence interval)
Sensitivity	75.8%	63.3 – 85.8%
Specificity	98.1%	93.5 – 99.8%
PPV	95.9%	86.0 - 99.5%
NPV	87.3%	79.9 – 92.7%

The results showed a high sensitivity and specificity to detect *Salmonella* using CerTest *Salmonella*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Salmonella*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Adenovirus</i>	<i>Escherichia coli</i> O:026	<i>Norovirus</i> GII
<i>Astrovirus</i>	<i>Escherichia coli</i> O:157	<i>Rotavirus</i>
<i>Calprotectin</i>	<i>Entamoeba dispar</i>	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	<i>Entamoeba histolytica</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Giardia lamblia</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i> antigen GDH	<i>Helicobacter pylori</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i> Toxin A	<i>Hemoglobin</i>	<i>Streptococcus pneumoniae</i>
<i>Clostridium difficile</i> Toxin B	<i>Lactoferrin</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium perfringens</i>	<i>Legionella</i>	<i>Transferrin</i>
<i>Cryptosporidium</i>	<i>Listeria monocytogenes</i>	<i>Yersinia enterocolitica</i> O:3
<i>Escherichia coli</i> O:111	<i>Norovirus</i> GI	<i>Yersinia enterocolitica</i> O:9

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1. ABUBAKAR, I. "A systematic review of the clinical, public health and cost-effectiveness of rapid diagnostic tests for the detection and identification of bacterial intestinal pathogens in faeces and food". Health Technology Assessment 2007; Vol. 11: No. 36.
2. OCHOA REPÁRAZ, J. "Protective ability of subcellular extracts from *Salmonella* Enteritidis and from a rough isogenic mutant against salmonellosis in mice". Vaccine 23 (2005), 1495-1501.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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