---

# CERTEST

## *Clostridium difficile* antigen GDH

---

ONE STEP *Clostridium difficile* antigen  
GDH CARD TEST

---

CERTEST BIOTEC S.L.

---

# CERTEST *Clostridium difficile* antigen GDH

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## ENGLISH

### INTENDED USE

CerTest *Clostridium difficile* antigen GDH one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH) in stool samples.

CerTest *Clostridium difficile* antigen GDH card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection.

### INTRODUCTION

*Clostridium difficile* (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. *C. difficile* is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide.

Disease due to *C. difficile* develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. *C. difficile* can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death.

*Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the organism.

The toxigenic culture (TC) is used as the gold standard technique to determine *Clostridium difficile* infection. This method consists in culture and isolation of *C. difficile* from faeces, followed by toxin testing of the isolate, a labour-intensive assay to obtain a result.

### TEST PRINCIPLE

CerTest *Clostridium difficile* antigen GDH is based on the principle of a qualitative immunochromatographic assay for the determination of *Clostridium difficile* Glutamate Dehydrogenase (GDH) in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against GDH and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-GDH) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigen of the diluted sample reacts with the red-coloured conjugate complex (anti-GDH monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-GDH antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no GDH antigen presence and yet, the antigen may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-GDH antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Clostridium difficile* antigen GDH test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.

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
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

#### Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the stool collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

### MATERIALS

#### MATERIALS PROVIDED

- CerTest *Clostridium difficile* antigen GDH card tests
- Instructions for use
- Stool collection tubes with diluent

#### MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

### TEST PROCEDURE

Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Clostridium difficile* antigen GDH card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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### INTERPRETATION OF RESULTS (please refer to the illustration below)

	GDH	Interpretation of the results
1.	- GREEN	There is no <i>Clostridium difficile</i> presence. No infection caused by <i>Clostridium difficile</i> .
2.	+ GREEN-RED	There is <i>Clostridium difficile</i> presence. <i>Clostridium difficile</i> infection, which might mean milder or severe (colitis) diarrhoea caused by <i>C. difficile</i> or asymptomatic patient, carrier of disease.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

**NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

**INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigen concentration to faint when the antigen concentration is close to the detection limit value of the test.
4. CerTest *Clostridium difficile* antigen GDH should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of *Clostridium difficile* in faecal samples; nevertheless, it can be due to toxigenic or non-toxigenic strains of *Clostridium difficile*. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible the antigen concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a *Clostridium difficile* determination should be carried out, on a sample from an enrichment culture.
7. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
8. Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

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### EXPECTED VALUES

*Clostridium difficile* is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases.

### PERFORMANCE CHARACTERISTICS

#### Analytical sensitivity (detection limit)

Detection limit value of CerTest *Clostridium difficile* antigen GDH is 0.78- 0.10 ng/mL of *Clostridium difficile* GDH recombinant protein.

#### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Clostridium difficile* GDH, CerTest) and an Evaluation Criteria. The Evaluation Criteria consisted on the evaluation in parallel of two rapid test (CerTest *Clostridium difficile* GDH) vs other rapid test from the competitor; the discrepant results were confirmed by qPCR technique.

The results were as follows:

		Evaluation Criteria		
		+	-	Total
IC test: CerTest <i>Clostridium difficile</i> GDH	+	51	0	51
	-	0	75	75
	Total	51	75	126

CerTest GDH vs Evaluation Criteria		
	Mean Value	95% CI (confidence interval)
<b>Sensitivity</b>	100.0%	93.0 – 100.0%
<b>Specificity</b>	100.0%	95.2 – 100.0%
<b>PPV</b>	100.0%	93.0 - 100.0%
<b>NPV</b>	100.0%	95.2 –100.0%

The results showed a high sensitivity and specificity to detect *Clostridium difficile* glutamate dehydrogenasa (GDH) using CerTest *Clostridium difficile* antigen GDH.

#### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Clostridium difficile* antigen GDH; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Salmonella typhimurium</i>
<i>Astrovirus</i>	<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Bovine Haemoglobin</i>	<i>Escherichia coli</i> O26	<i>Norovirus GI</i>	<i>Shigella dysenteriae</i>
<i>Bovine Lactoferrin</i>	<i>Escherichia coli</i> O157	<i>Norovirus GII</i>	<i>Shigella flexneri</i>
<i>Bovine Transferrin</i>	<i>Giardia</i>	<i>Pig haemoglobin</i>	<i>Shigella sonnei</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Rotavirus</i>	<i>Streptococcus pneumococcal</i>
<i>Campylobacter jejuni</i>	<i>Human Calprotectin</i>	<i>Salmonella enteritidis</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile</i> Toxin A	<i>Human Haemoglobin</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia Enterocolitica</i> O:3
<i>Clostridium difficile</i> Toxin B	<i>Human Lactoferrin</i>	<i>Salmonella typhi</i>	<i>Yersinia Enterocolitica</i> O:9
<i>Clostridium perfringens</i>	<i>Human Transferrin</i>		

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1. WREN, M.W.D, et al. "Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of *Clostridium difficile* associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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