

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 431411

Product Description: Glass Microfiber Prefilters, 49.5 mm Square, for Use with 250 mL Filter Funnels

Component Materials:
0.32 mm Thick Binder-Free Microfiber Glass, Grade 161

Attributes:
Product Dimensions - 48.5 mm Square x 0.32 mm Thick

Lot Number Designation:
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1