

Axygen® Automation Tips

For Molecular Devices® Automation Systems

AXYGEN


A Corning Brand

Axygen® automation tips for Molecular Devices workstations are designed and manufactured to provide the same high performance as those sold by the original equipment manufacturer. The tips are “out of the box” compatible with the workstations’ liquid handling equipment, providing precise and accurate automated sample delivery without the need for program (definition) changes.

Axygen automation tips are designed and manufactured by Corning as high quality alternatives to tips sold by original equipment manufacturers. Axigen products are manufactured in facilities that are ISO 9001 and 13485 certified with sealed, positive pressure, clean areas to minimize the risk of contamination. In addition, Axigen tips are validated at customer locations to ensure their compatibility and functionality.

Features

- ▶ Tips are ready to use “out of the box” – no definition changes are required
- ▶ 96-well and 384-well formats
- ▶ Clear or conductive (black) polypropylene
- ▶ Automated manufacturing provides tips free from defects and contamination
- ▶ RNase-/DNase-free
- ▶ Nonpyrogenic





Axygen® Automation Tips for Molecular Devices® Automation Systems

Ordering Information

Cat. No.	Workstation Compatibility	Description	Qty/Pk	Qty/Cs
96-well Format				
PK-200-R	FLIPR® Tetra	200 µL automation tip, clear	960	4,800
PK-200-BK-R	FLIPR Tetra	200 µL automation tip, black	960	4,800
RMT-200-BK-R	FlexStation®	200 µL automaton tip, black	960	4,800
RMT-200-C-R	FlexStation	200 µL automation tip, clear	960	4,800
TT-50-C-HTR	PatchXpress® 7000A	200 µL automation tip, black	96	2,304
TT-50-C-HTR-S	PatchXpress 7000A	200 µL automation tip, black, sterile	96	2,304

384-well Format

Cat. No.	Workstation Compatibility	Description	Qty/Pk	Qty/Cs
 FT-384-C-R	FLIPR Tetra	30 µL automation tip, clear	3,840	19,200
 FT-384-BK-R	FLIPR Tetra	30 µL automation tip, black	3,840	19,200
PK-384-R	FLIPR Tetra Type A	30 µL automation tip, clear	3,840	19,200
LT-384-R	FLIPR Liberty B (Liberty)	30 µL automation tip, black	3,840	19,200

For a full listing of all Axigen automation tips, visit www.corning.com/lifesciences.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional Axigen product or distributor information, e-mail us at CLSCustServ@corning.com, visit www.corning.com/lifesciences, or call 1.800.492.1110. Outside the United States, call 978.442.2200.

For Axigen technical information, e-mail us at AxgSupport@corning.com or call 1.800.429.9436. Outside the United States, call 510.494.8900.

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

CORNING | **FALCON** | **AXYGEN** | **GOSSSELIN** | **PYREX**

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks.
All other trademarks are the property of their respective owners.