


Premium Tips and Expert Service for All Your Workstations.

All From a Single Source.

AXYGEN


A Corning Brand



Simplify your supply chain by creating a custom consumable/service bundle.

Automation Advantage is the new service/consumables plan that offers premium pipet tips and expert maintenance for your automated liquid handling equipment at substantial savings.

Until now, labs have been mostly limited to the service/consumable plans offered by instrument manufacturers. Automation Advantage gives you the freedom to design a plan that best fits your needs—and your budget.



Premium Automation Tips

Axygen® products are the leading brand of replacement automation consumables, providing hundreds of high quality tips for all major brands and models including:

- ▶ Sterile tips for sensitive biological applications
- ▶ Filtered tips that reduce carryover contamination
- ▶ Conductive racks for 384-well tips to reduce static
- ▶ Nested rack configurations that optimize precious deck space

Proven Performance

All tip designs are tested in an in-house automation lab using the same automation equipment for which the tips are designed. The test results are available from your Corning Account Manager and can be downloaded from www.corning.com/axygen.

Advanced Manufacturing

All Axygen product molds are designed and built in-house to ensure that we have complete control over product specifications. Our advanced manufacturing processes include automated assembly to minimize contamination risks, electronic equipment for preventing bent tips, and deionizing technology to reduce static.

Rigorous Quality Control

In addition to our extensive in-house testing, Axygen automation tips are validated at customer locations to ensure proper function. The tips are supplied with a Certificate of Compliance.

Comprehensive Instrument Service

MRA has more than 50 years' experience servicing automated equipment. MRA's service engineers deliver a one-stop solution across a variety of platforms, providing a total support package that maximizes up-time, streamlines methods development, and anticipates the equipment service needs of today's busy laboratories.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

For additional product or technical information, please visit www.corning.com/lifesciences or call 800.492.1110. Customers outside the United States, please call +1.978.442.2200 or contact your local Corning sales office.

CORNING | **FALCON** | **AXYGEN** | **GOSSSELIN** | **PYREX**

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.

Full Suite of Services

MRA's goal is to be the most complete and trusted liquid handling robotics service and support group on the market. Their services include:

- ▶ 24 hour response time for service
- ▶ Preventative maintenance
- ▶ Operator qualification
- ▶ Product qualification
- ▶ Application support
- ▶ Method development
- ▶ Assay optimization

MRA sources parts from an ISO 13845 certified organization for all service components.

A Commitment to Quality

MRA maintains the highest quality of products and services for their customers by uniquely combining service engineering techniques, regulatory standards, and customer focused-support programs.

Coverage for All Major Liquid Handling Platforms

- ▶ Hamilton
- ▶ Tecan
- ▶ Beckman
- ▶ PerkinElmer

Set Yourself Free Today. Contact Us for More Information.

Don't get locked-down by OEM manufacturers. Let us design a service/consumable plan for your lab that will provide exactly the services and consumables you need at a competitive price. To get started, contact your Corning Account Manager or visit www.corning.com/axygen.

Automation Advantage is available in North America.