



Corning® 1536 Well Black Clear Bottom Low Base Microplate

Focused on Performance

- ▶ Designed for optimal compatibility with the MDS FLIPR Tetra®
- ▶ Microplate bottom is 1.8 mm from the Tetra detector
- ▶ Thin bottom for optical clarity
- ▶ The Corning CellBIND® Surface has allowed customers to switch from more expensive poly-D-lysine (PDL) coated microplates



CORNING

Corning® 1536 Well Black Clear Bottom Low Base Microplates

Product Specifications Chart

The information below can be used for setting up the plate with low base automation equipment.

Microplate Parameters	1536 Well Clear Bottom Low Base Microplates
Plate dimensions (L x W x H)	127.8 x 85.5 x 10.4 mm
Well A1 location (row offset)	7.87 mm
Well A1 location (column offset)	11.0 mm
Well center to well center distance	2.25 mm
Well diameter (top/bottom)	1.78/1.51 mm
Well depth	4.8 mm
Distance from well bottom to plate bottom	1.8 mm
Total well volume	12.5 µL
Working volume	Up to 10 µL

Ordering Information

Cat. No.	Description	Treatment	Sterile	Lid	Qty/Pk	Qty/Cs
3831	1536 Well Microplate, Black Clear Bottom, Low Base	Not Treated	No	No	10	50
3832	1536 Well Microplate, Black Clear Bottom, Low Base	Corning CellBIND® Surface	Yes	Yes	10	50
3832BC	1536 Well Microplate, Black Clear Bottom, Low Base with bar code labels	Corning CellBIND Surface	Yes	Yes	10	50
3833	1536 Well Microplate, Black Clear Bottom, Low Base	Corning CellBIND Surface	Yes	Yes	20	100
3835	1536 Well Microplate, Black Clear Bottom, Low Base	Not Treated	No	No	20	100
3836	1536 Well Microplate, Black Clear Bottom, Low Base	Tissue Culture Treated	Yes	Yes	20	100
3838	1536 Well Microplate, Black Clear Bottom, Low Base	Tissue Culture Treated	Yes	Yes	10	50
3838BC	1536 Well Microplate, Black Clear Bottom, Low Base with bar code labels	Tissue Culture Treated	Yes	Yes	10	50

For additional product or technical information, please visit www.corning.com/lifesciences or call 800.492.1110. Customers outside the United States, call +1.978.442.2200 or contact your local Corning sales office listed below.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA / PACIFIC
Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA
Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700
Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

The Corning Family of Brands

CORNING **costar** **cellgro** **PYREX** **AXYGEN** **GOSSSELIN**

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks.
All other trademarks are the property of their respective owners.
Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001