

TaqMan[®] GTXpress[™] Master Mix

Catalog Numbers 4401857, 4401890, 4401892, 4403311

Pub. No. 4403089 Rev. D

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *TaqMan[®] GTXpress[™] Master Mix User Guide* (Pub. No. 4402746). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

TaqMan[®] GTXpress[™] Master Mix reduces by half the time typically required to run standard genotyping protocols. Use TaqMan[®] GTXpress[™] Master Mix to perform fast genotyping of single nucleotide polymorphisms (SNPs). TaqMan[®] GTXpress[™] Master Mix can be used with genomic DNA or sample lysates and any TaqMan[®] genotyping assay, including TaqMan[®] SNP Genotyping Assays, Custom TaqMan[®] SNP Genotyping Assays, TaqMan[®] Drug Metabolism Genotyping Assays, and TaqMan[®] Predesigned Assay Reagents for Allelic Discrimination.

Contents and storage

Cat. No.	Number of 5- μ L reactions	Amount	Storage ^[1]
4403311	400	1 mL	2–8°C
4401892	4,000	10 mL	
4401890	20,000	50 mL	
4401857	40,000	2 x 50 mL	

^[1] See packaging for expiration date.

Additional information

For detailed procedures, safety guidelines and component descriptions, see *TaqMan[®] GTXpress[™] Master Mix User Guide* (Pub. No. 4402746).

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4403089

Revision	Date	Description
D	15 October 2018	Updated for manufacturer, general style, formatting, and branding.
C	September 2011	Baseline for this revision.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license.