

CORNING	Life Sciences	Document Number: LSR00025	Rev.: 3
Description (Class): Reference, Customer Technical Data Sheet – Aspirating Pipets, 2 mL			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Guaranteed sterile claim (individual)
 - Contamination control
- Consistent performance
 - Reliability, contamination control
- Aseptic presentation
 - Contamination control
- Non-pyrogenic
 - Low endotoxins
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Non-wettable
 - Minimum meniscus
- Individually wrapped packaging: “Thermoform” polyethylene-coated paper on one side, clear plastic on the other
 - Reduces “static cling”, quick and easy to open. Two ways to open: Peel or Pop-through

STERILITY:

Product is gamma irradiated

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf

CORNING	Life Sciences	Document Number: LSR00025	Rev.: 3
		Description (Class): Reference, Customer Technical Data Sheet – Aspirating Pipets, 2 mL	

Falcon® 2 mL Aspirating Pipet

Material: Polystyrene

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	357558	Individually wrapped, 200/case

Dimensions: = XXX mm (XXX inches)

TIP:

Length: = N/A
 Interior Diameter: = 1.625 (.064)
 Outside Diameter: = 3.73 (.147)

MOUTHPIECE:

Length: = N/A
 Interior Diameter: = 4.165 (.164)
 Outside Diameter: = 6.27 (.247)
 Plug Length = N/A

EXTRUSION:

Interior Diameter: = N/A
 Outside Diameter: = N/A
 Length: = N/A

PIPET SIZE: 2 mL

INCREMENTS: N/A


GRADUATIONS: N/A

OVERFILL CAPACITY: N/A

COLOR CODE: N/A

GENERAL DIMENSIONS:

A = 282.58 (11.125)
 B = 3.73 (.147) Diameter
 C = 6.27 (.247) Diameter
 D = 6.27 (.247) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.