#### Corning Incorporated Life Sciences

Registered ISO 9001:2008

# **Product Description**

Catalog Number: 431733

Product Description: Corning® Storage Bottle, 500 mL, Octagonal, PET, with 31.7 mm HDPE Screw Cap,

Sterile, 12/Tray, 24/Case

## **Component Materials:**

Bottle - PET, clear, non-cytotoxic, non-hemolytic, meets USP, Class VI requirements for plastic containers and closures.

 Cap - HDPE, natural, non-cytotoxic, non-hemolytic, meets USP, Class VI requirements for plastic containers and closures.

#### Attributes:

Capacity - 500 mL; Tolerance: ±50 mL at any mark

Product Dimensions - 173 mm Height (with cap locked); 171 mm Height (w/o cap); 74 mm Width;

42 mm Cap OD; 31.7 mm Neck ID

## Autoclaving:

This product is not suitable for autoclaving.

### Heating / Cooling:

This product is designed for use at temperatures ranging from -40°C to 50°C.

#### Special Features:

This product has a tamper evident cap. The product is supplied with the cap in place but not locked. This product is manufactured in an environment maintained to the standards of a Class 100 clean room.

#### Sterilization:

This product is sterilized using an irradiation method compliant with the requirements of ANSI/AAMI/ISO 11137 to achieve a Sterility Assurance Level (SAL) of  $10^{-6}$ .

#### Pyrogens:

This product is evaluated for pyrogens according to USP 31 and meets the requirements for classification as non-pyrogenic.

Results less than 0.1 EU/mL.

## **Expiration Statement:**

This product is defined to have and expiration date of 3 years from the date of manufacture.

Rev No: 2