



The following features are present in the vector based on the nucleotide sequence.

Cytomegalovirus (CMV) promoter.....	8-635
Green <i>Renilla</i> luciferase gene.....	646-1581
BGH poly(A) signal.....	1590-1715
SV40 origin/promoter.....	1716-2280
Puromycin resistant gene	2281-2880
SV40 poly(A) signal	3042-3075
beta-lactamase (Amp ^R) gene	3184-4044
pUC replication origin (pUC Ori).....	4223-5027
Transcriptional terminator (Ter)	5028-5635
Lac operator 1 (Lac O1)	5636-5656
Transcriptional pause site (TPS).....	5789-5860