

CytoScan™ Consumables Kit for Applied Biosystems™ NIMBUS™

Catalog Number 902168

Pub. No. 703166 Rev. 2

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Product use

The Applied Biosystems™ CytoScan™ Consumables Kit for Applied Biosystems™ NIMBUS™ contains plastic consumables that are needed to run target preparation on the NIMBUS™ using the CytoScan™ Assay Kit (Cat. No. 901808). The kit contains consumables sufficient for 4 × 25 or 2 × 49 reaction runs. See the *CytoScan™ Assay Automated Target Preparation Solution on Applied Biosystems™ NIMBUS™ User Guide* (Pub. No. 703141) for instructions on using the consumables.

Storage

Store the CytoScan™ Consumables Kit for Applied Biosystems™ NIMBUS™ at room temperature.

Components

Component	Usage per run		Number per kit
	25 samples	49 samples	
96-well Semi-skirted Plate	5	7	20
Tube, Screw Cap, No Ribs	6	6	24
1.2 mL Square Well Storage Plate	1	1	4
Reservoir 4 ROW PART PYR PP 73 mL	6	6	25
96-well Full Skirt PCR Plate, Clear	7	7	28
2.2 mL Storage Plate	3	3	12

Images




Figure 1 96-well semi-skirted plate



Figure 2 Tube, screw cap, no ribs



Figure 3 1.2 mL square well storage plate



Figure 4 Reservoir 4 row part PYR PP 73 mL




Figure 5 96-well full skirt PCR plate, clear




Figure 6 2.2 mL storage plate

Customer and technical support

Visit thermofisher.com/support for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 703166

Revision	Date	Description
2	06 September 2017	Update document to current template.
1	01 April 2013	Initial release.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. All other trademarks are the property of their respective owners.

thermofisher.com/support | thermofisher.com/askaquestion

thermofisher.com

ThermoFisher
SCIENTIFIC