

Catalog Number: 430236

Product Description: Corning® 250 mL centrifuge tube cushion

Component Materials:
Cushion - Polyetherimide

Product Dimensions:

Outer diameter	-	2.35 in.	Height	-	1.80 in.
Inner diameter @ top	-	2.26 in.	Tolerances	-	+/- 0.05 in.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Dimensional: Fit of tube in cushion.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Run number for that year.

Revision Date:

06-27-18

Rev No: 2