

Pancreatic Elastase Turbilatex[®]

Novel
PETIA assay
for the determination
of faecal
**pancreatic
elastase.**

CerTest
BIOTEC

CerTest Biotec, S.L.
Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
Tel. (+34) 976 520 354
Fax (+34) 976 106 268
certest@certest.es | turbilatex@certest.es
www.certest.es



Pancreatic elastase

Pancreatic elastase is produced in the acinar cells of the pancreas. A reduced level of this protein in the gastrointestinal tract is related to pancreatic insufficiency. It typically occurs as a result of progressive pancreatic damage, because of several pathological conditions.

Human pancreatic elastase 1 (E1) remains undegraded during intestinal tract, which makes pancreatic elastase a perfect biomarker.

Pancreatic elastase Turbilatex® is a non-invasive in vitro diagnostic method used to monitorize pancreatic function, revealing pancreatic insufficiency in the patient.

Pancreatic elastase Turbilatex® is:

- Fecal testing
- Quantitative
- Based on monoclonal antibodies detection
- Liquid and ready to use
- Long shelf-life
- Adaptable to any analyser

References	Description
TL-022EL100ED	Pancreatic elastase Combo kit 100 det.
TL-022EL200ED	Pancreatic elastase Combo kit 200 det.
TL-022EL03E	Pancreatic elastase Sample diluent

Our provided kits:



FOB Turbilatex®




Calprotectin Turbilatex®



Transferrin Turbilatex®



Pancreatic elastase Turbilatex®



H. pylori Turbilatex®