

TECHNICAL SHEET

PRODUCT NAME: Entamoeba dispar recombinant outer membrane protein

Code: MT-25EDP, MT-25EDPs1

Physical State: Liquid

Source: Escherichia coli

Description: Entamoeba dispar recombinant outer membrane protein

expressed and purified from Escherichia coli.

Appearance: Clear

Purity: >95 %

Storage Conditions: Storage at 2-8°C.

Health & Safety Information: This product contains sodium azide as preservative. Although

the amount is 0.09%, this substance is poisonous and hazardous

and it should be handled by trained staff only

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in

accordance with the respective national regulations.

Presentation: Carbonate buffer with additives and salts, pH 8.3, 0.09%

Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 21/01/2019



www.certest.es