

## **TECHNICAL SHEET**

PRODUCT NAME: Human Haemoglobin protein (native extract)

Code: MT-29HHB, MT-29HHBs1

Physical State: Liquid

Source: Human

**Description:** Human native haemoglobin protein purified.

Appearance: Clear

**Purity:** >95 %

**Storage Conditions:** Storage at 2-8°C. Light protect.

Health & Safety Information: This product contains sodium azide as preservative. Although

the amount is 0.09%, this substance is poisonous and hazardous

and it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in

accordance with the respective national regulations.

**Presentation** Tris buffer with with additives and salts , pH: 8.1, 0.09% Sodium

Azide (NaN<sub>3</sub>)

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 22/01/2019



www.certest.es