

Laboratory Refrigerator



BPR-5V50(G)



BPR-5V100(G)

Technical Parameters

Model	BPR-5V50(G)	BPR-5V100(G)
External Size(W*D*H)mm	495*548*683	495*548*934
Internal Size(W*D*H)mm	415*405*455	415*405*706
Capacity	50L	100L
Temperature Range	2-8°C	
Temperature Accuracy	0.1°C	
Control System	Microprocessor Control	
Alarm	Audible and visual alarm for: High and low temperature,Door Ajar	
Refrigeration Type	Forced air refrigeration system	
Refrigerant	R600a	
Climate Type	N, SN	
Glass Door	Toughened with electric heating function,no condensation	
Construction	Constructure	Unibody design& High Pressure and High Density Cyclopentane Foaming
	Internal material	Cold-rolled steel coated with anti-bacteria powder
	External material	Cold-rolled steel coated with anti-bacteria powder
Consumption	85W	95W
Power Supply	AC220V±10%, 50/60Hz; AC110V±10%, 50/60Hz	
Shelves	2pcs, cold drawn steel wire impregnated shelf	3 pcs,cold drawn steel wire impregnated shelf
Door Lock	1 pc	
Accessory	Standard: USB Port Optional: Remote Alarm; temperature recorder, printer (USB port can not be selected at the same time)	
Package Size(W*D*H)mm	560*668*858	560*668*1109
Gross Weight(kg)	49	61

Laboratory Refrigerator



260L



360L

Technical Parameters

Model	BPR-5V260	BPR-5V360
External Size(W*D*H) mm	540*635*1900	660*655*1980
Capacity	260L	360L
Temp. Range	2°C~8°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor Control, Large Screen LED Display	
Alarm	Audible and visual alarm for: High and low temperature, System & sensor failure, Door ajar, Power failure alarm, Low backup battery alarm.	
Refrigeration Type	Forced air refrigeration system	
Refrigerant	R134a, CFC Free	
Condenser & Evaporator	Bundy tube condenser	
Defrost	Frostless design	
Construction	Structure	Unibody design & mono-assembly foaming (Rigid polyurethane insulation material)
	Internal material	HIPS engineering plastics, abrasion and acid resistant
	External material	Cold-rolled steel coated with anti-bacteria powder
Shelves	Adjustable shelves made of high quality steel wire	
	4 pcs	5 pcs
Consumption	200W	350W
Power Supply	AC220V±10%, 50/60Hz	
Standard Accessory	LED lamp & shelves, USB, probe access port, door keys	
Package Size(W*D*H) mm	600*670*2150	720*750*2150
Gross Weight(kg)	87	124