

**Corning Incorporated
Life Sciences**

Registered
ISO 9001:2008

Product Description

Catalog Number: 3479

Product Description: Costar ® Netwell™, 24mm, 74 µm mesh size, in a 6 well cluster plate with lid

Component Materials:

- Plate/Lid - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures.
- Netwell body - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures with grey concentrate.
- Mesh bottom - Polyester

Product Dimensions:

- Length of Plate - 5.030 in.
- Width of Plate - 3.365 in.
- Tolerances of Dimensions - +/- .050 in.
- Bottom of Netwell to plate - .050 in.
- Height with Lid - 1.100 in.

Sterilization:

The product is irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization*.

Sterility Assurance Level: SAL 10⁻³

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

- Visual Attributes: Visual examination of the product.
- Packaging: Inspection for seal and barrier integrity, accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 6