

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 26798-40

Product Description: PYREX® Closure, Polypropylene, Natural, 40 mm

Component Materials:
Polypropylene (PP)

Attributes:
Product Dimensions - 29 mm OD x 31 mm High

Autoclaving:
This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

Heating / Cooling:
Polypropylene may be heated to 140°C.

Polypropylene may be cooled to -196°C.

Vacuum / Pressure:
This product was not designed for and is not rated for vacuum applications.

This product was not designed for and is not rated for pressure applications.

Sterilization:
Shipped Not Sterile

Lot Number Designation:
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:
5 AUG 2013

Rev No: 1