

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 1730-2C

Product Description: Corning® Snap-Seal, Container, Sample, Opaque, Not Sterile, 1.5 OZ / 45 mL

Component Materials:
Polypropylene (PP)

Attributes:

Capacity - 45 mL \pm 2.2 mL (5%)
Product Dimensions - 30 mm OD x 84 mm Height

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes. NOTE: Autoclaving disposable products IS NOT recommended.

Heating / Cooling:

This product can be heated to 140°C.

This product can be cooled to -196°C.

Vacuum / Pressure:

This product is not designed for and is not rated for vacuum or pressure applications.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

17 JUN 2013

Rev No: 1