



The following features are present in the vector based on the nucleotide sequence.

Cytomegalovirus (CMV) promoter.....	8-635
<i>Cypridina</i> luciferase gene	646-2307
BGH poly(A) signal.....	2316-2441
SV40 origin/promoter.....	2442-3006
Puromycin resistant gene	3007-3606
SV40 poly(A) signal	3768-3801
beta-lactamase (Amp ^R) gene	3910-4770
pUC replication origin (pUC Ori).....	4949-5753
Transcriptional terminator (Ter)	5754-6361
Lac operator 1 (Lac O1)	6362-6382
Transcriptional pause site (TPS).....	6515-6586